





Thermo Scientific Laboratory Products Catalog



- Handheld Pipetting
- Automated Liquid Handling
- Microplate Instrumentation
- Nucleic Acid Purification and Electroporation
- PCR Instruments and Consumables
- Microcentrifuge Tubes and Racks
- Compliance & Calibration Services
- Leasing

The best products
and services
for demanding lab applications

Thermo
SCIENTIFIC






Thermo Scientific Laboratory Products Catalog

*Discover why more laboratories rely on
Thermo Scientific liquid handling solutions...*

This catalog is designed to help you select the best products to meet virtually any laboratory requirement. From handheld pipetting and automated liquid handling to microplate instrumentation, nucleic acid purification, PCR amplification and more – this resource makes it easy to meet your application challenges with a complete portfolio of industry-leading Thermo Scientific and Molecular BioProducts offerings.

For additional information, contact your sales representative or
visit www.thermoscientific.com



MBP[®]
Molecular BioProducts

Handheld Pipetting

Thermo Scientific Finnpiquette Pipetters

| | |
|--|----|
| Thermo Scientific Finnpiquette F1 Single Channel Pipetters | 8 |
| Thermo Scientific Finnpiquette F1 Multichannel Pipetters | 10 |
| Thermo Scientific Finnpiquette F1 GLP Kits | 11 |
| Thermo Scientific Finnpiquette F2 Single Channel Pipetters | 12 |
| Thermo Scientific Finnpiquette F2 Multichannel Pipetters | 14 |
| Thermo Scientific Finnpiquette F2 GLP Kits | 15 |
| Thermo Scientific Finnpiquette Stepper Pipetter | 16 |
| Thermo Scientific Finnpiquette Multistepper Pipetter | 17 |
| Thermo Scientific Finnpiquette Novus Single Channel Electronic Pipetters | 18 |
| Thermo Scientific Finnpiquette Novus Multichannel Electronic Pipetters | 19 |
| Thermo Scientific Finnpiquette Dispensers | 20 |
| Thermo Scientific Finnpiquette Accessories | 21 |
| Thermo Scientific Reagent Reservoirs | 22 |

Thermo Scientific Finntip® Pipette Tips

| | |
|--|----|
| Thermo Scientific Finntip Compatibility Chart | 24 |
| Thermo Scientific Finntip Pipette Tips | 25 |
| Thermo Scientific Finntip Extended Length Pipette Tips | 27 |
| Thermo Scientific Finntip Wide Orifice Pipette Tips | 28 |
| Thermo Scientific Finntip BioCon Pipette Tips | 29 |
| Thermo Scientific Finntip Filter Pipette Tips | 30 |
| Thermo Scientific Finntip Flex Pipette Tips | 31 |
| Thermo Scientific Finntip Flex Filter Pipette Tips | 33 |
| Thermo Scientific Finntip Stepper Tips | 34 |
| Thermo Scientific Finntip Multistepper Pipette Tips | 35 |

Thermo Scientific Matrix Electronic Pipetters and Accessories

| | |
|--|----|
| Thermo Scientific Matrix Single-Channel Electronic Pipetters | 36 |
| Thermo Scientific Matrix Multichannel Electronic Pipetters | 37 |
| Thermo Scientific Matrix EXP Electronic Pipetters | 39 |
| Thermo Scientific Matrix Equalizer Electronic Multichannel Pipetters | 40 |
| Thermo Scientific S1 Pipet Filler | 44 |
| Thermo Scientific Pipetter Accessories | 45 |
| Thermo Scientific Nunc Serological Pipets | 46 |
| Thermo Scientific Matrix Pipette Tip Compatibility Chart | 47 |
| Thermo Scientific Matrix Pipette Tips | 48 |
| Thermo Scientific Matrix TallTip Extended Length Pipette Tips | 50 |
| Thermo Scientific Matrix Filter Pipette Tips | 51 |
| Thermo Scientific Matrix TallTip Filter Pipette Tips | 52 |
| Thermo Scientific Matrix Reagent Reservoirs | 54 |

Molecular BioProducts Pipette Tips

| | |
|--|----|
| Molecular BioProducts and Pure Pipette Tips with MicroPoint Tip Design | 57 |
| Molecular BioProducts and Pure Pipette Tips with SoftFit Tip Design | 58 |
| Molecular BioProducts Pure Pipette Tips with SoftFit L Design | 59 |
| Molecular BioProducts ART® Barrier Tips with MicroPoint Design | 60 |
| Molecular BioProducts ART Barrier Ultra Micro Tips with MicroPoint Design | 61 |
| Molecular BioProducts Low Retention Pipette Tips with ART Self-Sealing Barrier | 62 |
| Molecular BioProducts Low Retention Pipette Tips | 64 |
| Molecular BioProducts Pipette Tips Compatibility Chart | 66 |

Automated Liquid Handling

Thermo Scientific Multidrop Reagent Dispensers and Accessories

| | |
|--|----|
| Thermo Scientific Multidrop Combi nL Reagent Dispenser | 76 |
| Thermo Scientific Multidrop Combi Reagent Dispenser | 77 |
| Thermo Scientific Multidrop 384 Reagent Dispenser | 78 |
| Thermo Scientific Multidrop DW Reagent Dispenser | 79 |
| Thermo Scientific Multidrop Combi, 384 and DW Dispensing Cassettes | 80 |
| Standard tube dispensing cassette | 80 |
| Small Tube Dispensing Cassette | 80 |
| SMART dispensing cassettes | 80 |
| Thermo Scientific FILLit Software for Multidrop Reagent Dispensers | 82 |
| Multidrop Combi nL FILLit Software | 82 |

Thermo Scientific Matrix Dispensers and Accessories

| | |
|---|----|
| Thermo Scientific Matrix WellMate Microplate Dispenser | 83 |
| Thermo Scientific Matrix WellMate Disposable Tubing Cartridges | 84 |
| Thermo Scientific Matrix WellMate Stacker Base Unit | 85 |
| Thermo Scientific Matrix PlateMate Plus/WellMate Stacker Chimneys | 86 |
| Thermo Scientific Matrix Hydra II Liquid Handling System | 87 |
| Thermo Scientific Matrix Hydra DT Pipetting Workstation | 88 |

Thermo Scientific Versette Pipetting Workstations and Accessories

| | |
|--|----|
| Thermo Scientific Versette Pipetting Workstation | 89 |
| Thermo Scientific Versette Pipetting Heads | 90 |
| Thermo Scientific Versette Accessories | 91 |

Automation Tips and Accessories

| | |
|--|-----|
| Thermo Scientific Versette ClipTip Automation Tips | 92 |
| Thermo Scientific Matrix D.A.R.T.s Tips | 93 |
| Thermo Scientific Matrix Filtered D.A.R.T.s Tips | 94 |
| Thermo Scientific PocketTip D.A.R.T.s Automation Tips | 95 |
| Thermo Scientific Matrix D.A.R.T.s Tip Transfer Tool | 96 |
| Thermo Scientific Matrix Disposable Automation Reservoirs | 97 |
| Molecular BioProducts BioRobotix ART Filter Tips for Automated Workstations | 98 |
| Molecular BioProducts BioRobotix Pipet Tips, Standard, Black | 100 |
| Molecular BioProducts BioRobotix Pipette Tips: Standard, Natural, Nonsterile | 101 |

Microplate Stackers and Accessories

| | |
|--|-----|
| Thermo Scientific RapidStak Automated Microplate Stacker with Polara RS Software | 103 |
| Thermo Scientific RapidStak Accessories: Microplate Stacks | 104 |
| Thermo Scientific RapidStak Instrument Drivers | 104 |
| Thermo Scientific Polara RS | 105 |
| Thermo Scientific RapidStak DLL Programming Kits | 105 |

Microplate Movers

| | |
|---|-----|
| Thermo Scientific Orbitor RS Microplate Mover | 107 |
|---|-----|

Microplate Handlers and Accessories

| | |
|---|-----|
| Thermo Scientific CataLyst Express Microplate Handler | 108 |
|---|-----|

Plastic Reservoirs • Storage Blocks

| | |
|---|-----|
| Thermo Scientific Nalgene Disposable Robotic Reservoirs | 109 |
| Thermo Scientific Nunc Disposable Plastic Reservoirs | 110 |
| Thermo Scientific Matrix Deepwell Storage Blocks | 111 |

Microplate Instrumentation

Incubator/Shakers

| | |
|--|-----|
| Thermo Scientific iEMS Incubator/Shaker | 114 |
| Thermo Scientific iEMS Incubator/Shaker HT | 115 |

Multimode Readers

| | |
|--|-----|
| Thermo Scientific Appliskan Multimode Reader | 116 |
| Thermo Scientific Varioskan Flash Multimode Reader | 118 |

Microplate Fluorometer and Luminometers

| | |
|---|-----|
| Thermo Scientific Fluoroskan Ascent FL Microplate Fluorometer and Luminometer | 120 |
| Thermo Scientific Fluoroskan Ascent Microplate Fluorometer | 122 |
| Thermo Scientific Luminoskan Ascent Microplate Luminometer | 124 |

Microplate Photometers

| | |
|--|-----|
| Thermo Scientific Multiskan EX Microplate Photometer | 126 |
| Thermo Scientific Multiskan FC Microplate Photometer | 127 |

Microplate Spectrophotometers

| | |
|---|-----|
| Thermo Scientific Multiskan GO Microplate Spectrophotometer | 128 |
| Thermo Scientific Multiskan Spectrum Microplate Spectrophotometer | 130 |

Microplate Washers

| | |
|--|-----|
| Thermo Scientific Wellwash 4 Mk2 Microplate Washer | 132 |
| Thermo Scientific Wellwash Microplate Washer | 133 |
| Thermo Scientific Wellwash Versa Microplate Washer | 134 |

Software

| | |
|---|-----|
| Thermo Scientific Ascent Software | 136 |
| Thermo Scientific SkanIt Software | 137 |

Nucleic Acid Purification and Electroporation

Thermo Scientific KingFisher Systems and Accessories

| | |
|--|-----|
| Thermo Scientific KingFisher Flex Magnetic Particle Processors | 140 |
| Thermo Scientific Consumables for KingFisher Flex Systems | 141 |
| Thermo Scientific KingFisher mL Magnetic Particle Processors | 142 |
| Thermo Scientific Consumables for KingFisher mL Systems | 143 |
| Thermo Scientific KingFisher Magnetic Particle Processors | 144 |
| Consumables for KingFisher Systems | 144 |
| Thermo Scientific BindIt Software for KingFisher Instruments | 145 |
| Thermo Scientific KingFisher Kits | 146 |

Electroporation Cuvettes

| | |
|--|-----|
| Molecular BioProducts Electroporation Cuvettes | 149 |
|--|-----|

PCR and Real-Time PCR Instruments and Consumables

Thermal Cyclers

| | |
|---|-----|
| Thermo Scientific Arktik Thermal Cycler | 152 |
| Thermo Scientific Piko Thermal Cycler | 153 |

Real-Time PCR System

| | |
|--|-----|
| Thermo Scientific Piko Real Real-Time PCR System | 154 |
|--|-----|

PCR Tubes, Plates and Caps

| | |
|--|-----|
| Thermo Scientific Piko PCR Plates | 156 |
| Thermo Scientific Piko Plate Illuminator | 157 |
| Sealing Films for Piko | 157 |
| Molecular BioProducts PCR Tubes | 159 |
| Molecular BioProducts PCR Plates and Caps | 160 |
| Molecular BioProducts PCR Strip Tubes and Caps | 160 |

Surface Decontamiant

| | |
|--|-----|
| Molecular BioProducts DNA AWAY Surface Decontaminant | 161 |
| Molecular BioProducts RNase AWAY Surface Decontaminant | 161 |
| Molecular BioProducts EasyStart PCR Mix-in-a-Tube | 162 |

Storage Reaction Tubes

| | |
|---|-----|
| Molecular BioProducts HotStart Storage Reaction Tubes | 163 |
|---|-----|

Microcentrifuge Tubes and Racks

Centrifuge and Microcentrifuge Tubes

| | |
|---|-----|
| Snap-Cap Centrifuge Tubes | 166 |
| Locking Lid Microcentrifuge Tubes | 166 |
| Capless Microcentrifuge Tubes | 167 |
| Screw Cap Microcentrifuge Tubes – Conical | 167 |
| Screw Cap Microcentrifuge Tubes – Free Standing | 168 |
| Screw Caps for Microcentrifuge Tubes | 168 |

Microcentrifuge and Microtiter Tubes

| | |
|---------------------------------------|-----|
| Microtiter Tubes | 169 |
| Specialty Microcentrifuge Tubes | 169 |

Flipper Racks • Cryogenic and Microtube Racks

| | |
|---|-----|
| Molecular BioProducts 4-Way Flipper Racks | 170 |
| Molecular BioProducts 81-Well Cryogenic Rack with Lid | 170 |
| Molecular BioProducts FlipStrip Microtube Racks with Lids | 171 |
| Molecular BioProducts Reversible Microtube Racks with Lids | 172 |
| Molecular BioProducts 96-Well Flipper Microtube Racks with Lids | 172 |

Compliance and Calibration Services

| | |
|---|-----|
| Compliance and Calibration Services | 174 |
|---|-----|

Leasing

| | |
|--|-----|
| Benefits of Leasing | 177 |
| Standard Lease and Finance Options | 178 |

Note: All trademarks and service marks noted with an * (e.g. Finnpiquette*) are the property of Thermo Fisher Scientific Inc. and its subsidiaries. All other trademarks and registered marks are the properties of their respective owners. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

For additional trademark information, please see page 184.



Handheld Pipetting

Proven leadership in the palm of your hand.

For 40 years, we have led the way in liquid handling and handheld pipetting products. Advanced ergonomics, dependable accuracy and precision, uncompromised safety – all are the hallmarks of our innovative product design. In 1971, we introduced the Thermo Scientific Finnpipette pipetter, the world's first continuously variable micropipetter. In 1976, we introduced the world's first multichannel pipetter. Since then, we have continuously enhanced our products to meet any pipetting challenge across the widest range of liquid handling applications.

Thermo Scientific Matrix pipetters are ergonomically designed to enhance user comfort and productivity while assuring accuracy of results. Thermo Scientific Matrix pipette tips and reservoirs are carefully designed to complement our pipetters, providing maximum performance and superior results. For the latest news and information regarding good laboratory pipetting, visit www.thermoscientific.com/GLP

With Molecular BioProducts (MBP) products, you can now choose from the largest selection of pipette tips in the world for virtually all of the leading pipetter brands in the industry.

For additional information about Thermo Scientific pipetting solutions, visit thermoscientific.com/pipette.





Handheld Pipetting

Thermo Scientific FinnpiPETTE Pipetters

| | |
|---|----|
| Thermo Scientific FinnpiPETTE F1 Single Channel Pipetters | 8 |
| Thermo Scientific FinnpiPETTE F1 Multichannel Pipetters | 10 |
| Thermo Scientific FinnpiPETTE F1 GLP Kits | 11 |
| Thermo Scientific FinnpiPETTE F2 Single Channel Pipetters | 12 |
| Thermo Scientific FinnpiPETTE F2 Multichannel Pipetters | 14 |
| Thermo Scientific FinnpiPETTE F2 GLP Kits | 15 |
| Thermo Scientific FinnpiPETTE Stepper Pipetter | 16 |
| Thermo Scientific FinnpiPETTE MultistepPetter | 17 |
| Thermo Scientific FinnpiPETTE Novus Single Channel Electronic Pipetters | 18 |
| Thermo Scientific FinnpiPETTE Novus Multichannel Electronic Pipetters | 19 |
| Thermo Scientific FinnpiPETTE Dispensers | 20 |
| Thermo Scientific FinnpiPETTE Accessories | 21 |
| Thermo Scientific Reagent Reservoirs | 22 |

Thermo Scientific Finntip Pipette Tips

| | |
|--|----|
| Thermo Scientific Finntip Compatibility Chart | 24 |
| Thermo Scientific Finntip Pipette Tips | 25 |
| Thermo Scientific Finntip Extended Length Pipette Tips | 27 |
| Thermo Scientific Finntip Wide Orifice Pipette Tips | 28 |
| Thermo Scientific Finntip BioCon Pipette Tips | 29 |
| Thermo Scientific Finntip Filter Pipette Tips | 30 |
| Thermo Scientific Finntip Flex Pipette Tips | 31 |
| Thermo Scientific Finntip Flex Filter Pipette Tips | 33 |
| Thermo Scientific Finntip Stepper Tips | 34 |
| Thermo Scientific Finntip MultistepPetter Pipette Tips | 35 |

Thermo Scientific Matrix Pipetters and Accessories

| | |
|--|----|
| Thermo Scientific Matrix Single-Channel Electronic Pipetters | 36 |
| Thermo Scientific Matrix Multichannel Electronic Pipetters | 37 |
| Thermo Scientific Matrix EXP Electronic Pipetters | 39 |
| Thermo Scientific Matrix Equalizer Electronic Multichannel Pipetters | 40 |
| Thermo Scientific S1 Pipet Filler | 44 |
| Thermo Scientific Pipetter Accessories | 45 |
| Thermo Scientific Nunc Serological Pipets | 46 |
| Thermo Scientific Matrix Pipette Tip Compatibility Chart | 47 |
| Thermo Scientific Matrix Pipette Tips | 48 |
| Thermo Scientific Matrix TallTip Extended Length Pipette Tips | 50 |
| Thermo Scientific Matrix Filter Pipette Tips | 51 |
| Thermo Scientific Matrix TallTip Filter Pipette Tips | 52 |
| Thermo Scientific Matrix Reagent Reservoirs | 54 |

Molecular BioProducts Pipette Tips

| | |
|--|----|
| Molecular BioProducts and Pure Pipette Tips with MicroPoint Tip Design | 57 |
| Molecular BioProducts and Pure Pipette Tips with SoftFit Tip Design | 58 |
| Molecular BioProducts Pure Pipette Tips with SoftFit L Design | 59 |
| Molecular BioProducts ART [®] Barrier Tips with MicroPoint Design | 60 |
| Molecular BioProducts ART Barrier Ultra Micro Tips with MicroPoint Design | 61 |
| Molecular BioProducts Low Retention Pipette Tips with ART Self-Sealing Barrier | 62 |
| Molecular BioProducts Low Retention Pipette Tips | 64 |
| Molecular BioProducts Pipette Tips Compatibility Chart | 66 |



Thermo Scientific Finnpiquette F1 Single Channel Pipetters



Finnpiquette* F1 Fixed- and Adjustable-Volume Pipetters set the standard in pipetting with superior comfort and performance.

Our F1 pipetter incorporates state-of-the-art innovations for ease of use and proven results. The Advanced Volume Gearing mechanism (AVG) is a self-supporting modular volume adjustment mechanism that enhances pipetting performance. It also delivers dependable accuracy and precision with real durability with its separate pipetter body.

The F1 Adjustable-Volume Pipetter offers a simple "set-and-forget" lock mechanism to adjust the volume, adjustable finger rest, and soft touch tip ejection while an innovative antimicrobial treatment protects the instrument against contamination.

Details

- **Antimicrobial surface:** Silver (Ag) ions create a naturally inhibiting surface to prevent the growth of bacteria, microbes and fungus
- **"Set-and-forget" pipetting button** enables an easy locking mechanism for volume adjustment, and rotates freely when not locked to prevent accidental volume drift during pipetting
- **Adjustable finger rest** can be adjusted 120° for optimal pipetting position
- **Super blow-out** piston with volumes of 50 μL and below ensures delivery of micro-size drops
- **ID tag space** for personal markings, blank tags supplied with the pipetter

Ergonomics

- **Extremely light** and smooth plunger action reduces the risk of Repetitive Strain Injury (RSI)
- **Lightweight construction** enables longer pipetting periods without fatigue
- **Soft-touch tip ejector** reduces the tip ejection forces by up to 50%
- **Adjustable finger rest** for superior comfort

Ordering Information: Ideal for use with Thermo Scientific Finntip pipette tips, as specified in the ordering table

Includes: Finnpiquette F1 pipetter, calibration tool, blank ID tags, shelf hanger, Instructions manual, and calibration certificate

Warranty: Five years with Web registration

F1 Fixed-Volume Pipettors

| Cat. No. | Description | Accuracy | Precision |
|----------|---------------|------------|-----------|
| 4651000 | 1 μ L | \pm 4.0% | 4.0% |
| 4651010 | 5 μ L | \pm 1.4% | 1.4% |
| 4651020 | 10 μ L | \pm 0.9% | 0.8% |
| 4651130 | 20 μ L | \pm 0.7% | 0.5% |
| 4651030 | 25 μ L | \pm 0.6% | 0.5% |
| 4651040 | 50 μ L | \pm 0.6% | 0.4% |
| 4651050 | 100 μ L | \pm 0.4% | 0.3% |
| 4651140 | 200 μ L | \pm 0.4% | 0.3% |
| 4651060 | 250 μ L | \pm 0.4% | 0.3% |
| 4651070 | 500 μ L | \pm 0.3% | 0.3% |
| 4651080 | 1000 μ L | \pm 0.3% | 0.3% |
| 4651090 | 2000 μ L | \pm 0.3% | 0.2% |
| 4651100 | 3000 μ L | \pm 0.3% | 0.2% |
| 4651110 | 5000 μ L | \pm 0.3% | 0.2% |
| 4651120 | 10000 μ L | \pm 0.3% | 0.2% |

F1 Adjustable-Volume Pipettors

| Cat. No. | Range | Increments | Accuracy | Precision | Finntip | Color Code |
|----------|----------------------|---------------|--------------------|--------------|------------------------------------|------------|
| 4641010 | 0.2-2 μ L, micro | 0.002 μ L | \pm 12.0 to 2.5% | 10.0 to 2.0% | Flex 10, 10, 50 | Pink |
| 4641020 | 0.5-5 μ L, micro | 0.01 μ L | \pm 6.0 to 1.5% | 5.0 to 1.0% | Flex 10, 10, 50 | Pink |
| 4641030 | 1-10 μ L, micro | 0.02 μ L | \pm 2.5 to 1.0% | 2.0 to 0.5% | Flex 10, 10, 50 | Pink |
| 4641040 | 1-10 μ L | 0.02 μ L | \pm 3.5 to 1.0% | 3.0 to 0.8% | Flex 200, 250 Univ. | Yellow |
| 4641050 | 2-20 μ L, micro | 0.02 μ L | \pm 3.0 to 1.0% | 2.5 to 0.4% | 50 micro | Turquoise |
| 4641060 | 2-20 μ L | 0.02 μ L | \pm 3.0 to 1.0% | 2.5 to 0.4% | Flex 200, 250 Univ. | Yellow |
| 4641130 | 5-50 μ L, micro | 0.02 μ L | \pm 3.0 to 0.6% | 2.5 to 0.3% | 50 micro | Turquoise |
| 4641140 | 5-50 μ L | 0.02 μ L | \pm 3.0 to 0.6% | 2.5 to 0.3% | Flex 200, 300, 250 Univ, 200 Ext | Yellow |
| 4641070 | 10-100 μ L | 0.2 μ L | \pm 3.0 to 0.8% | 1.0 to 0.2% | Flex 200, 250 Univ., 300, 200 Ext. | Yellow |
| 4641080 | 20-200 μ L | 0.2 μ L | \pm 1.8 to 0.6% | 0.7 to 0.2% | Flex 200, 250 Univ., 300, 200 Ext. | Yellow |
| 4641090 | 30-300 μ L | 1 μ L | \pm 1.5 to 0.6% | 0.6 to 0.2% | Flex 300, 300 | Orange |
| 4641100 | 100-1000 μ L | 1 μ L | \pm 1.0 to 0.6% | 0.6 to 0.2% | Flex 1000, 1000, 1000 Ext. | Blue |
| 4641110 | 0.5-5 mL | 0.01 mL | \pm 2.0 to 0.5% | 0.8 to 0.2% | 5 mL | Green |
| 4641120 | 1-10 mL | 0.02 mL | \pm 2.0 to 0.5% | 0.8 to 0.2% | 10 mL, Flex 10 mL Ext. | Red |

Thermo Scientific Finn timer F1 Multichannel Pipettors



Finn timer F1 Multichannel Pipettors offer the perfect combination of design innovation and user comfort for maximum productivity.

As with the single-channel Finn timer F1 pipettors, the AVG mechanism ensures a high level of accuracy and precision. The soft-touch tip ejection is especially useful on multichannel models, significantly reducing the effort needed to eject multiple tips. In addition, the super blow-out function in low volume models ensures accurate dispensing for even the smallest volumes.

Our multichannel pipettors are available in 8-, 12- and 16-channel models, with a selection of volume ranges, and work best with Thermo Scientific Finntips.

Details

- **Antimicrobial surface:** Silver (Ag) ions create a naturally-inhibiting surface to prevent the growth of bacteria, microbes and fungus
- **“Set-and-forget” pipetting button** enables an easy locking mechanism for volume adjustment, and rotates freely when not locked, preventing accidental volume drift during pipetting
- **Adjustable finger rest** can be adjusted 120° for optimal pipetting position
- **ID tag space** for personal markings, blank tags supplied with the pipetter
- **Super blow-out piston** with volumes of 50 µL and below ensures delivery of micro-size drops
- Compatible with 96-well and 384-well for both high and low volume microplates

Ergonomics

- **Soft-touch tip ejector** reduces the tip ejection forces by up to 50%
- **Extremely light** and smooth plunger action reduces the risk of Repetitive Strain Injury (RSI)
- **Lightweight construction** enables longer pipetting periods without fatigue

Ordering Information: Ideal for use with Thermo Scientific Finntip pipette tips, as specified in the ordering table

Includes: Finn timer F1 multichannel pipetter, calibration tool, module removal tool, blank ID tags, shelf hanger, Instruction manual and calibration certificate

Warranty: Five years with Web registration

| Cat. No. | Range | Increments | Accuracy | Precision | Finntip | Color Code |
|-------------------|-----------|------------|---------------|-------------|-------------------------------|------------|
| 8-Channel | | | | | | |
| 4661000 | 1-10 µL | 0.02 µL | ±12.0 to 2.4% | 8.0 to 1.6% | Flex 10, 10 micro, 50 micro | Pink |
| 4661010 | 5-50 µL | 0.1 µL | ±5.0 to 1.5% | 2.0 to 0.7% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4661020 | 10-100 µL | 0.2 µL | ±5.0 to 1.3% | 2.0 to 0.5% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4661030 | 30-300 µL | 1 µL | ±5.0 to 1.0% | 2.0 to 0.3% | Flex 300, 300 | Orange |
| 12-Channel | | | | | | |
| 4661040 | 1-10 µL | 0.02 µL | ±12.0 to 2.4% | 8.0 to 1.6% | Flex 10, 10 micro, 50 micro | Pink |
| 4661050 | 5-50 µL | 0.1 µL | ±5.0 to 1.5% | 2.0 to 0.7% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4661060 | 10-100 µL | 0.2 µL | ±5.0 to 1.3% | 2.0 to 0.5% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4661070 | 30-300 µL | 1 µL | ±5.0 to 1.0% | 2.0 to 0.3% | Flex 300, 300 | Orange |
| 16-Channel | | | | | | |
| 4661080 | 1-10 µL | 0.02 µL | ±12.0 to 2.4% | 8.0 to 1.6% | Flex 10 (384), 50 micro | Purple |
| 4661090 | 5-50 µL | 0.1 µL | ±5.0 to 1.5% | 2.0 to 0.7% | 50 micro | Turquoise |

Thermo Scientific Finnpiquette F1 GLP Kits

Finnpiquette F1 GLP Kits include state-of-the-art Finnpiquette F1 pipettors along with everything you need for superior pipetting.

All Finnpiquette F1 Good Laboratory Pipetting kits offer superior performance plus unique features, including adjustable finger rest for user comfort, "set and forget" pipetting button, and an antimicrobial surface.

Each complete kit provides everything you need for good laboratory pipetting.

Details

- **Adjustable finger rest** for improved ergonomics and comfortability
- **"Set-and-forget" pipetting button** offers easy and secured volume adjustment
- **Antimicrobial surface** provides enhanced protection
- **Large easy-to-read display** prevents eye strain
- **Soft-touch tip ejector** enables light end effortless tip ejection

Ordering Information: The Finnpiquette F1 GLP kit also includes samples of Finntip Flex tips for maximum performance and ergonomics, and a Good Laboratory Pipetting guide with extensive information about pipetting and factors that affect pipetting results.

Warranty: Five years with Web registration

| Cat. No. | Description | Volume Range |
|----------|----------------------------|---|
| 4700850 | F1 Kit 1 (1 to 1000 µL) | 1-10 µL, 10-100 µL, 100-1000 µL |
| 4700860 | F1 Kit 2 (0.2 to 1000 µL) | 0.2 µL, 2-20 µL, 20-200 µL, 100-1000 µL |
| 4700865 | F1 Kit 3 (10 to 10,000 µL) | 10-100 µL, 100-1000 µL, 1000-10000 µL |
| 4701060 | F1 Kit 4 (2 to 2000 µL) | 2-20 µL, 20-200 µL, 100-1000 µL |



F1 Kit 1 (1-1000 µL)

Finnpiquette F1, 1-10 µL
 Finnpiquette F1, 10-100 µL
 Finnpiquette F1, 100-1000 µL
 Finntip Flex tips: 10, 1 × 96; 200, 1 × 96; 1000, 1 × 96
 F-stand; F1 pen; F1 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir samples

F1 Kit 2 (0.2 to 1000 µL)

Finnpiquette F1, 0.2-2 µL
 Finnpiquette F1, 2-20 µL univ
 Finnpiquette F1, 20-200 µL
 Finnpiquette F1, 100-1000 µL
 Finntip Flex tips: 10, 1 × 96; 200, 2 × 96; 1000, 1 × 96
 F-stand; F1 pen; F1 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir samples

F1 Kit 3 (10 to 10,000 µL)

Finnpiquette F1, 10-100 µL
 Finnpiquette F1, 100-1000 µL
 Finnpiquette F1, 1000-10000 µL
 Finntip Flex tips: 200, 1 × 96; 1000, 1 × 96;
 FT 10 mL, 1 × 24
 F-stand; F1 pen; F1 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir samples

F1 Kit 4 (2 to 1000 µL)

Finnpiquette F1, 2-20 µL
 Finnpiquette F1, 20-200 µL
 Finnpiquette F1, 100-1000 µL
 Flex tips: (2) 200 µL, 10 × 96;
 100-1000 µL, 10 × 96
 F-stand; F1 pen; F1 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir samples

Thermo Scientific Finnpiquette F2 Single Channel Pipetters



Finnpiquette F2 Fixed- and Adjustable-Volume Pipetters are durable and fully autoclavable for everyday pipetting needs.

Experience superior comfort, performance, reliability, and repeatability in one of the lightest pipetters available. These pipetters include a Advanced Volume Gearing mechanism (AVG) – a self-supporting modular volume adjustment mechanism that enhances pipetting performance.

Details

- **Large Ergovisio display** with white numbers on black background for increased visibility
- **Double-action pipetting button** prevents accidental volume changes
- **Super blow-out piston** with volumes of 50 μL and below ensures delivery of micro-size drops
- **Fully autoclavable** for dependable contamination prevention
- **ID tag** space for personal markings; blank tags supplied with the pipetter

Ergonomics

- **Extremely light** and smooth plunger action reduces the risk of Repetitive Strain Injury (RSI)
- **Lightweight construction** enables longer pipetting periods without fatigue
- **Soft-touch tip ejector** reduces the tip ejection forces by up to 50%
- **Wide and supportive** finger rest for comfortable pipetting

Ordering Information: Ideal for use with Thermo Scientific Finntip pipette tips, as specified in the ordering table

Includes: Finnpiquette F2 pipetter, calibration tool, blank ID tags, shelf hanger, instruction manual and calibration certificate

Warranty: Five years with Web registration

F2 Fixed-Volume Pipetters

| Cat. No. | Description | Accuracy | Precision |
|----------|---------------------|-------------|-----------|
| 4652000 | 1 μL | $\pm 4.0\%$ | 4.0% |
| 4652010 | 5 μL | $\pm 1.4\%$ | 1.4% |
| 4652020 | 10 μL | $\pm 0.9\%$ | 0.8% |
| 4652030 | 25 μL | $\pm 0.6\%$ | 0.5% |
| 4652040 | 50 μL | $\pm 0.6\%$ | 0.4% |
| 4652050 | 100 μL | $\pm 0.4\%$ | 0.3% |
| 4652060 | 250 μL | $\pm 0.4\%$ | 0.3% |
| 4652070 | 500 μL | $\pm 0.3\%$ | 0.3% |
| 4652080 | 1000 μL | $\pm 0.3\%$ | 0.3% |
| 4652090 | 2000 μL | $\pm 0.3\%$ | 0.2% |
| 4652100 | 3000 μL | $\pm 0.3\%$ | 0.2% |
| 4652110 | 5000 μL | $\pm 0.3\%$ | 0.2% |
| 4652120 | 10000 μL | $\pm 0.3\%$ | 0.2% |
| 4652130 | 20 μL | $\pm 0.6\%$ | 0.5% |
| 4652140 | 200 μL | $\pm 0.4\%$ | 0.3% |



F2 Adjustable-Volume Pipetters

| Cat. No. | Range | Increments | Accuracy | Precision | Compatible Tips | Color Code |
|----------|----------------------|---------------|--------------------|--------------|---|------------|
| 4642010 | 0.2-2 μ L, micro | 0.002 μ L | \pm 12.0 to 2.5% | 10.0 to 2.0% | Flex 10, 10, 50 | Pink |
| 4642020 | 0.5-5 μ L, micro | 0.01 μ L | \pm 6.0 to 1.5% | 5.0 to 1.0% | Flex 10, 10, 50 | Pink |
| 4642030 | 1-10 μ L, micro | 0.02 μ L | \pm 2.5 to 1.0% | 2.0 to 0.5% | Flex 10, 10, 50 | Pink |
| 4642040 | 1-10 μ L | 0.02 μ L | \pm 3.5 to 1.0% | 3.0 to 0.8% | Flex 200, 250 Univ. | Yellow |
| 4642050 | 2-20 μ L, micro | 0.02 μ L | \pm 3.0 to 1.0% | 2.5 to 0.4% | 50 micro | Turquoise |
| 4642060 | 2-20 μ L | 0.02 μ L | \pm 3.0 to 1.0% | 2.5 to 0.4% | Flex 200, 250 Univ. | Yellow |
| 4642120 | 5-50 μ L, mico | 0.1 μ L | \pm 0.3 to 0.15% | 2.5 to 0.3% | 50 micro | Turquoise |
| 4642130 | 5- 50 μ L | 0.1 μ L | \pm 0.3 to 0.15% | 2.5 to 0.3% | Flex 200, Flex 300, 250 Univ., 200 Ext. | Yellow |
| 4642070 | 10-100 μ L | 0.2 μ L | \pm 3.0 to 0.8% | 1.0 to 0.2% | Flex 200, 250 Univ., 300, 200 Ext. | Yellow |
| 4642080 | 20-200 μ L | 0.2 μ L | \pm 1.8 to 0.6% | 0.7 to 0.2% | Flex 200, 250 Univ., 300, 200 Ext. | Yellow |
| 4642090 | 100-1000 μ L | 1 μ L | \pm 1.0 to 0.6% | 0.6 to 0.2% | Flex 1000, 1000, 1000 Ext. | Blue |
| 4642100 | 0.5-5 mL | 0.01 mL | \pm 2.0 to 0.5% | 0.8 to 0.2% | 5 mL | Green |
| 4642110 | 1-10 mL | 0.02 mL | \pm 2.0 to 0.5% | 0.8 to 0.2% | 10 mL, Flex 10 mL Ext. | Red |

Thermo Scientific Finn timer F2 Multichannel Pipettors



The Finn timer F2 Multichannel Pipettor includes an **Advanced Volume Gearing mechanism (AVG)** that provides a high level of accuracy and precision.

AVG is a self-supporting modular volume adjustment mechanism that enhances pipetting performance. Because it is separated from the pipettor body, its accuracy, precision and durability are significantly improved. The AVG mechanism is thermally isolated from the pipettor body, minimizing the effects of hand warmth on the accuracy.

These multichannel pipettors are available in 8-, 12- and 16-channel models with a selection of volume ranges. As with the single-channel Finn timer F2 pipettors, the AVG mechanism ensures a high level of accuracy and precision. In addition, the super blow-out function in the low volume models provides accurate dispensing for even the smallest volumes.

Details

- **Large Ergovisio** display with white numbers on black background for increased visibility
- **Double-action pipetting button:** prevents accidental volume changes
- **ID tag** space for personal markings, blank tags supplied with the pipettor
- **Super blow-out piston** with volumes of 50 μL and below ensures delivery of micro-size drops
- **Fully autoclavable:** entire pipettor can be autoclaved in one piece for convenient decontamination prevention
- **Wide, supportive finger rest** for user comfort

Ergonomics

- **Soft-touch tip ejector** reduces the tip ejection forces by up to 50%
- **Extremely light** and smooth plunger action reduces the risk of Repetitive Strain Injury (RSI)
- **Lightweight construction** enables longer pipetting periods without fatigue

Ordering Information: Ideal for use with Thermo Scientific Finntip pipette tips, as specified in the ordering table

Includes: Finn timer F2 multichannel pipettor, calibration tool, module removal tool, blank ID tags, shelf hanger, instruction manual and calibration certificate

Warranty: Five years with Web registration

| Cat. No. | Range | Increment | Accuracy | Precision | Finntip | Color Code |
|-------------------|----------------------|--------------------|--------------------|-------------|-------------------------------|------------|
| 8-Channel | | | | | | |
| 4662000 | 1-10 μL | 0.02 μL | ± 12.0 to 2.4% | 8.0 to 1.6% | Flex 10, 10, 50 | Pink |
| 4662010 | 5-50 μL | 0.1 μL | ± 5.0 to 1.5% | 2.0 to 0.7% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4662020 | 10-100 μL | 0.2 μL | ± 5.0 to 1.3% | 2.0 to 0.5% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4662030 | 30-300 μL | 1 μL | ± 5.0 to 1.0% | 2.0 to 0.3% | Flex 300, 300 | Orange |
| 12-Channel | | | | | | |
| 4662040 | 1-10 μL | 0.02 μL | ± 12.0 to 2.4% | 8.0 to 1.6% | Flex 10, 10, 50 | Pink |
| 4662050 | 5-50 μL | 0.1 μL | ± 5.0 to 1.5% | 2.0 to 0.7% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4662060 | 10-100 μL | 0.2 μL | ± 5.0 to 1.3% | 2.0 to 0.5% | Flex 200, 250 Univ., 200 Ext. | Yellow |
| 4662070 | 30-300 μL | 1 μL | ± 5.0 to 1.0% | 2.0 to 0.3% | Flex 300, 300 | Orange |
| 16-Channel | | | | | | |
| 4662080 | 1-10 μL | 0.02 μL | ± 12.0 to 2.4% | 8.0 to 1.6% | Flex 10 (384), 50 micro | Purple |
| 4662090 | 5-50 μL | 0.1 μL | ± 5.0 to 1.5% | 2.0 to 0.7% | 50 micro | Turquoise |

Thermo Scientific Finnpiquette F2 GLP Kits

Finnpiquette F2 GLP Kits provide durable, high-performance Finnpiquette F2 pipetters in an all-inclusive kit containing everything you need for convenient, superior pipetting.

Finnpiquette F2 pipetters feature an AVG (Advanced Volume Gearing mechanism) and a large display for improved usability and comfort.

Each complete kit offers everything you need for good laboratory pipetting.

Details

- **AVG: Volume gearing mechanism** for accuracy and precision
- **Fully autoclavable** for dependable protection
- **Soft-touch tip ejection** enables light tip ejection
- **Large display** prevents eye strain
- **Very light** pipetting forces for greater accuracy with less fatigue

Ordering Information: The Finnpiquette F2 GLP kit includes samples of Finntip Flex tips – the most advanced traditional tip – offering very low attachment and ejection forces. The Good Laboratory Pipetting Guide contains comprehensive pipetting information, from choosing the right type to detailed decontamination procedures.

Warranty: Five years with Web registration

| Cat. No. | Description | Volume Range |
|----------|----------------------------|---|
| 4700870 | F2 Kit 1 (1 to 1000 µL) | 1-10 µL, 10-100 µL, 100-1000 µL |
| 4700880 | F2 Kit 2 (0.2 to 1000 µL) | 0.2 µL, 2-20 µL, 20-200 µL, 100-1000 µL |
| 4700885 | F2 Kit 3 (10 to 10,000 µL) | 10-100 µL, 100-1000 µL, 1000-10000 µL |
| 4701070 | F2 Kit 4 (2 to 2000 µL) | 2-20 µL, 20-200 µL, 100-1000 µL |



F2 Kit 1 (1-1000 µL)

Finnpiquette F2, 1-10 µL
 Finnpiquette F2, 10-100 µL
 Finnpiquette F2, 100-1000 µL
 Flex tips: 10, 1 × 96; 200, 1 × 96;
 1000, 1 × 96
 F-stand; F2 pen; F2 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir demo pack

F2 Kit 2 (0.2 to 1000 µL)

Finnpiquette F2, 0.2-2 µL
 Finnpiquette F2, 2-20 µL univ
 Finnpiquette F2, 20-200 µL
 Finnpiquette F2, 100-1000 µL
 Flex tips: 10, 12 × 96; 200, 2 × 96;
 1000, 1 × 96
 F-stand; F2 pen; F2 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir demo pack

F2 Kit 3 (10 to 10,000 µL)

Finnpiquette F2, 10-100 µL
 Finnpiquette F2, 100-1000 µL
 Finnpiquette F2, 1000-10000 µL
 Flex tips: 200, 1 × 96; 1000, 1 × 96;
 FT 10 mL, 1 × 24
 F-stand; F2 pen; F2 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir demo pack

F2 Kit 4 (2 to 2000 µL)

Finnpiquette F2, 2-20 µL
 Finnpiquette F2, 20-200 µL
 Finnpiquette F2, 100-1000 µL
 Flex tips: (2) 200 µL, 10 × 96; 100-1000 µL,
 10 × 96
 F-stand; F2 pen; F2 brochure
 Good Laboratory Pipetting Guide
 Reagent reservoir demo pack

Thermo Scientific Finnpiquette Stepper Pipetter



The Finnpiquette Stepper Pipetter is a lightweight, easy-to-use repeater pipetter designed for one-handed dispensing.

Our lightweight, ergonomically-designed pipetter allows you to rapidly dispense up to 45 times in succession without refilling. Operating on the positive displacement principle, the Finnpiquette Stepper is ideally suited for work with aggressive and viscous liquids.

Details

- **Makes up to 45 deliveries** before refilling is necessary
- **Positive displacement design** ensures accurate results, even for viscous, volatile or other problematic liquids
- **Lightweight, ergonomic handle** and wide finger rest minimize operator hand strain
- **Operation is simple** – just attach the appropriate tip size, fill and prime the tip, and dial the selected volume
- **Volume is easy to set**, with the aid of the convenient volume-setting chart located on the handle

Ergonomics

- **Universal handle** also fits the multichannel Finntip Steppers
- **Lightweight construction** minimizes hand fatigue during repetitive dispensing
- **Separate tip ejector** for safe and effortless tip disposal
- **Wide finger rest** minimizes operator hand strain
- **Module rotates a full 360°**, optimizing efficiency for both right and left-handed users

Ordering Information: Universal Stepper handle is also available with an eight-channel module. For convenient storage, pipetter stand is available (Cat. No. 14-387-14)

Includes: Unit comes with one adapter (nonsterile) necessary for use with 25 or 50 mL tips

Required Accessories: Uses the Thermo Scientific Finntip Stepper tip, a plunger-style tip available in nonsterile or sterile (individually wrapped) form. This tip comes in seven sizes from 0.5 to 50 mL, each with five different delivery volumes, for an overall dispensing range of 10 to 5000 μ L. Only the disposable tip comes in contact with the sample, eliminating the risk of contaminating the pipetter body.

| Cat. No. | Description | Volume Range |
|----------|----------------------|-----------------|
| 4540 000 | Finnpiquette Stepper | 10-5000 μ L |

Thermo Scientific FinnpiPETTE MultistepPetter

The FinnpiPETTE MultistepPetter reduces the risk of error and increases productivity in immunoassay applications.

This lightweight, easy-to-use repeater pipetter provides simultaneous dispensing via eight channels. Designed especially for work with microplates, the FinnpiPETTE MultistepPetter accelerates routine laboratory procedures for greater productivity.

Details

- **8-channel air-displacement pipetter**, designed especially for work with microplates
- **Easy to use** – simply attach, fill, and prime the module and dial the selected volume
- **Delivers 50, 100, 150, 200, or 250 µL**
- **With a single loading** it can deliver 24 x 50 µL=192 wells, 12 x 100 µL=96 wells, 8 x 150 µL=64 wells, 6 x 200 µL=48 wells, or 5 x 250 µL=40 wells
- **All eight channels are calibrated** to simultaneously dispense exactly the same volume of liquid
- **Pipetter capacity:** 1250 µL
- **Accuracy:** ±5.0 to 2.0%; precision: 5.0 to 2.0%
- **Handy volume chart** on the handle lets you know the maximum number of deliveries possible at a given volume setting

Ergonomics

- **Universal handle** also fits the seven single-channel Finntip Steppers
- **Lightweight construction** minimizes hand fatigue during repetitive dispensing
- **Separate tip ejector** for safe and effortless tip disposal
- **Wide finger rest** minimizes operator hand strain
- **Module rotates a full 360°**, optimizing efficiency for both right and left-handed users

Ordering Information: Accepts Finntip MultistepPetter 1500 µL tips that are extra long, with more air space between the liquid and tip cone to prevent contamination

| Cat. No. | |
|----------|-----------------------|
| 4540-500 | Handle and module |
| 2206-590 | 8-Channel Module only |

| Cat. No. | Channels | Capacity | Disp.vol. | Increment | Volume | Inaccuracy | | Imprecision | | Finntip |
|-------------|----------|----------|-----------|-----------|--------|------------|------|-------------|-------|-----------------|
| | | | | | | µL | % | s.d.* µL | CV %* | |
| 4540 500 1) | 8 | 1250 µL | 50-250 | 50 µL | 250 µL | ±5.0 | ±2.0 | 5.0 | 2.0 | MultistepPetter |
| 2206 590 2) | 8 | 1250 µL | 50-250 | 50 µL | 250 µL | ±5.0 | ±2.0 | 5.0 | 2.0 | MultistepPetter |

* s.d. = Standard Deviation, CV = Coefficient of Variation ** Factory calibration limits. 1) Handle and module 2) Module only



Thermo Scientific Finnpiquette Novus Single Channel Electronic Pipetters



Finnpiquette Novus Single Channel Electronic Pipetters offer performance with simple operation and reduced stress for comfortable, productive pipetting.

Finnpiquette Novus delivers performance with greater ease of use and reduced risk of injury. Unique to Novus, the adjustable index-finger trigger action reduces common repetitive stress injuries caused by thumb-driven pipettes, and the easy-to-follow prompts guide Novus users through 10 pipetting functions.

Details

- **Intuitive user interface** features a simple button layout and no abbreviations
- **Seven languages:** Operate in the user's native language (English, French, German, Italian, Japanese, Spanish, Swedish)
- **Backlit technology** eliminates surrounding light reflections and improves contrast in low-light conditions
- **10 pipetting functions and nine aspirate/dispense speeds:** Ultimate in flexibility with options not available in manual pipettes
- **Personalize up to nine programs** for your most common protocols: Save time and ensure accurate and precise pipetting
- **Easy in-lab, 2- or 1-point calibration:** Provides accurate and precise results in different applications
- **Fully autoclavable tip cone:** Prevents cross-contamination
- **Long-life lithium-ion battery** allows approximately 4,000 pipetting operations; recharges in approximately one hour
- **A large volume range from 1 µL up 10 mL:** Meets virtually any application need

Ergonomic Design

- **Natural index finger pipetting operation:** Lets the thumb relax while pipetting
- **Adjustable finger rest:** For balance and comfort
- **Soft-touch tip ejection:** Minimizes thumb strain
- **Extremely lightweight construction:** Work longer without fatigue

Ordering Information: Ideal for use with Thermo Scientific Finntip pipette tips as specified in the ordering table. Also compatible with most universal-fit tips.

Includes: Finnpiquette Novus, Universal charger, service tool, instruction manual and calibration certificate

Warranty: Two years with Web registration

| Cat. No. | Volume Range (calibration) | Volume Range (functional**) | Increment | Accuracy %* | Precision %* | Color code | Compatible Finntip |
|----------|----------------------------|-----------------------------|-----------|--------------|--------------|------------|---------------------------------|
| 46200000 | 1-10 µL micro | 0.5-10 µL | 0.01 µL | ±3.5 to ±1.0 | 3.0 to 0.5 | Pink | Flex 10, 10, 20 micro, 50 micro |
| 46200100 | 1-10 µL | 0.5-10 µL | 0.01 µL | ±7.0 to ±1.0 | 6.0 to 0.5 | Yellow | Flex 200, 250 Univ., 200 Ext. |
| 46200200 | 5-50 µL micro | 2.5-50 µL | 0.1 µL | ±3.0 to ±0.8 | 2.5 to 0.3 | Turquoise | 50 |
| 46200300 | 5-50 µL | 2.5-50 µL | 0.1 µL | ±3.0 to ±0.8 | 2.5 to 0.3 | Yellow | Flex 200, 250 Univ., 200 Ext. |
| 46200400 | 10-100 µL | 5-100 µL | 0.1 µL | ±3.0 to ±0.8 | 1.0 to 0.2 | Yellow | Flex 200, 250 Univ., 200 Ext. |
| 46200500 | 30-300 µL | 15-300 µL | 1 µL | ±3.0 to ±0.6 | 0.7 to 0.2 | Orange | Flex 300, 300 |
| 46200600 | 100-1000 µL | 50-1000 µL | 1 µL | ±3.0 to ±0.6 | 0.6 to 0.2 | Blue | Flex 1000, 1000, 1000 Ext. |
| 46200700 | 0.5-5 mL | 0.25-5 mL | 0.01 mL | ±3.0 to ±0.6 | 0.8 to 0.2 | Green | 5 mL |
| 46200800 | 1-10 mL | 0.5-10 mL | 0.01 mL | ±3.0 to ±0.6 | 0.8 to 0.2 | Red | 10 mL, Flex 10 mL Ext. |

*Factory calibration limits achieved under strictly controlled conditions (ISO 8655).

**The Functional volume range indicates the volume range that the pipette can cover in the Stepper, Sequential Stepper and Sequential Aspirate functions.

Thermo Scientific Finnpiquette Novus Multichannel Electronic Pipettors

Finnpiquette Novus Multichannel Pipettors offer superior performance for microplate applications.

Finnpiquette Novus pipettors deliver outstanding performance with greater ease of use and less risk of injury. Unique to Novus pipettors, the adjustable index-finger trigger action prevents common repetitive stress injuries caused by thumb-driven pipettes, and the easy-to-follow prompts guide users through 10 pipetting functions.

Repetitive dispensing with stepper functionality excels at filling microplates by saving time and effort compared to manual pipettors. Dispensing of two selected volumes with an air cap between, followed by a mix step with DILUTE+MIX function, is excellent for the preparation of standard curves. The Novus 1200 µL 8-channel pipette is ideal for high volume stepping of an entire plate, (e.g., dispense 100 µL up to 12 columns). The pipettor is also a practical tool for filling high volume deep well plates. For 384-microplate applications, choose the Novus 16-channel model.

Details

- **Intuitive user interface** features a simple button layout and no abbreviations
- **Seven languages:** Operate in the user's native language (English, French, German, Italian, Japanese, Spanish, Swedish)
- **Backlit technology** eliminates surrounding light reflections and improves contrast in low-light conditions
- **10 pipetting functions and nine aspirate/dispense speeds:** Ultimate in flexibility and variety not available in manual pipettes
- **Personalize up to nine programs** for your most common protocols: Save time and ensure accurate and precise pipetting
- **Easy in-lab, 2- or 1-point calibration:** For accurate and precise results in different applications
- **Fully autoclavable tip cones***:** Prevents cross-contamination
- **Long-life lithium-ion battery** allows approximately 4,000 pipetting operations; recharges in approximately one hour
- **A large volume range from 1 µL up 1200 µL:** Accommodates many applications

Ergonomic Design

- **Natural index finger pipetting operation:** Lets the thumb relax while pipetting
- **Adjustable finger rest:** For balance and comfort
- **Soft-touch tip ejection:** Minimizes thumb strain
- **Extremely lightweight construction:** Work longer without fatigue



Ordering Information: Ideal for use with Thermo Scientific Finntip pipette tips, as specified in the ordering table. Also compatible with most universal-fit tips.

Includes: Finnpiquette Novus, universal charger, service tool, instruction manual and calibration certificate

Warranty: Two years with Web registration

| Cat. No. | Channels | Range (calibration) | Range (functional**) | Increment | Accuracy %* | Precision %* | Color code | Compatible Finntip |
|----------|----------|---------------------|----------------------|-----------|---------------|--------------|------------|---------------------------------|
| 46300000 | 8 | 1-10 µL | 0.5-10 µL | 0.1 µL | ±12.0 to ±2.4 | 8.0 to 1.6 | Pink | Flex 10, 10, 20 micro, 50 micro |
| 46300100 | 12 | 1-10 µL | 0.5-10 µL | 0.1 µL | ±12.0 to ±2.4 | 8.0 to 1.6 | Pink | Flex 10, 10, 20 micro, 50 micro |
| 46300200 | 8 | 5-50 µL | 2.5-50 µL | 0.1 µL | ±5.0 to ±1.5 | 2.0 to 0.7 | Yellow | Flex 200, 200 Ext, 250 |
| 46300300 | 12 | 5-50 µL | 2.5-50 µL | 0.1 µL | ±5.0 to ±1.5 | 2.0 to 0.7 | Yellow | Flex 200, 200 Ext, 250 |
| 46300700 | 16 | 5-50 µL | 2.5-50 µL | 0.1 µL | ±5.0 to ±1.5 | 2.0 to 0.7 | Turquoise | 50 |
| 46300400 | 8 | 30-300 µL | 15-300 µL | 1 µL | ±5.0 to ±1.0 | 2.0 to 0.3 | Orange | Flex 300, 300 |
| 46300500 | 12 | 30-300 µL | 15-300 µL | 1 µL | ±5.0 to ±1.0 | 2.0 to 0.3 | Orange | Flex 300, 300 |
| 46300800 | 8 | 100-1200 µL | 50-1200 µL | 1 µL | ±3.0 to ±1.0 | 0.9 to 0.2 | Turquoise | Flex 1200 |

*Factory calibration limits achieved under strictly controlled conditions (ISO 8655). **The Functional volume range indicates the volume range that the pipette can cover in the Stepper, Sequential Stepper and Sequential Aspirate functions. ***Except the 8-channel 1200 µL model

Thermo Scientific Finnpiquette Dispensers



Finnpiquette Dispensers – a superior tool for dosing liquids from reagent bottles.

Finnpiquette Dispensers protect users even when handling aggressive liquids. Six models are available, covering volumes between 0.2 and 60 mL.

Details

- Ergonomically designed, non-slip piston handle
- Excellent chemical and thermal resistance
- Fully autoclavable without disassembling
- Precise self-locking system in the volume setting
- Volumes are absolutely reproducible
- Dispensing head cap prevents dripping
- Universal fit for reagent bottles
- Ergonomic design prevents work fatigue
- Comprehensive volume range
- Threaded suction hose for bubble-free suction

Ergonomics:

- Ergonomic and non-slip handle knob
- Dispenser surface roughened for non-slip grip
- Easy to adjust volume adjustment wheel

Warranty: One year

Certifications: Individual calibration certificate with serial number

Notes: Limitations include HF, liquids that attack Halar, FEP or Hastelloy suspensions because solid particles might block the valves.

| Cat. No. | Description | Accuracy | Precision |
|----------|-------------------|----------|-----------|
| 4421120 | Range: 0.2-1.0 mL | ±6.0 µL | 2.0 µL |
| 4421130 | Range: 0.4-2.0 mL | ±12.0 µL | 4.0 µL |
| 4421140 | Range: 1-5 mL | ±30.0 µL | 10.0 µL |
| 4421150 | Range: 2-10 mL | ±60.0 µL | 20.0 µL |
| 4421160 | Range: 5-30 mL | ±180 µL | 60.0 µL |
| 4421170 | Range: 10-60 mL | ±360 µL | 120 µL |



Thermo Scientific Finnpiquette Accessories

Finnpiquette Stands provide safe and convenient manual and electronic pipetter storage.

The Finnpiquette F-stand accommodates six pipettes. The sturdy crossbar structure with a slot for each pipetter enables it to sit snugly in the stand. The Finnpiquette F-stand is suitable for both single channel and multichannel models.

The Finnpiquette multichannel stand is a multi-position stand for a multichannel pipetter. The pipetter can be stored in an upright or horizontal position. The stand has a small footprint and can hold any multichannel Finnpiquette.

The Finnpiquette Novus stand is designed for convenient storage of an electronic Finnpiquette Novus pipetter. Both single channel and multichannel Finnpiquette Novus models can be accommodated with ease.



Finnpiquette F-stand

| Cat. No. | Stands for Manual and Electronic Finnpiquettes | Qty./ Unit |
|----------|---|------------|
| 9420400 | F-Stand, Linear Stand, white, 6 positions | 1 |
| 9420390 | Multichannel Pipetter Stand, 1 position | 1 |
| 9420290 | Pipette Stand, gray, for 6 pipettors | 1 |
| 9420320 | Pipette Mini-stand, gray, for 3 pipettors | 1 |
| 9420340 | Carousel Stand, for 6 manual or 3 electronic pipettors, including cover | 1 |
| 2206040 | Finnpiquette Shelf Hanger | 1 |
| 9420360 | Novus Single Pipetter Stand | 1 |
| 2209480 | Novus Adapter for Carousel Stand (3 pcs) | 3 |
| 2209490 | Novus Adapter for Carousel Stand (1 pc) | 1 |

| Cat. No. | Finnpiquette Dispenser Accessories | Qty./ Unit |
|----------|---|------------|
| 4421260 | Bottle Adapters A25 & A38 (one of each) | 2 |
| 4421270 | Extended Suction Tube, 0,2-10 mL | 1 |
| 4421280 | Extended Suction Tube, 5-60 mL | 1 |
| 4421290 | Extended Discharge Tubing, 0,2-10 mL | 1 |
| 4421300 | Extended Discharge Tubing, 5-60 mL | 1 |



Finnpiquette multichannel stand

Finnpiquette Novus stand



Thermo Scientific Reagent Reservoirs



Our Reagent Reservoirs provide a wide variety of solutions to support all of your pipetting procedures.

Our product innovation extends beyond Thermo Scientific pipettors and tips. Even our reagent reservoirs have been designed to provide tailored solutions. Four variations of reservoirs are available: 25 mL, 25 mL with divider, 100 mL and a reusable polypropylene 75 mL reservoir.

Details

- **Single use reservoirs:** 25 mL, 25 mL divided and 100 mL
- **Trough within a trough** maximizes the amount of liquid accessible to pipette tips when using small amounts of reagent
- **Pour-off spouts on all four corners** reduce spillage when pouring reagents out of reservoir
- **Graduations on inside wall** enable quick measurement of remaining liquid
- **Extra-wide base** adds rigidity and stability to reservoir, helping to avoid spills
- **25 mL divided reservoir** allows the pipetting of two different reagents with up to eight channels on one side, and up to four on the other; maximizes recovery if using only single channel and a small amount of reagent
- **Sterile reservoirs** are packaged in snap-and-tear bags for sterility and easy opening

60 mL reagent basin

- Unique V-shape: ideally suited for multichannel pipettors
- Autoclavable polypropylene

8-Channel reagent vessel

- Capacity: 10 mL/channel
- Autoclavable polypropylene

100 mL reagent basin

- Capacity: 100 mL
- Autoclavable polypropylene



| Cat. No. | Description |
|----------|---|
| 9510027 | Reagent Basin, 60 mL |
| 9510037 | Reservoir, 8-well, 10 mL |
| 9510047 | Reservoir, 1-well, 100 mL |
| 95128085 | Reagent Reservoir, Disposable, 100 mL |
| 95128095 | Reagent Reservoir, Disposable, w/Divider, 25 mL |
| 95128080 | Reagent Reservoir Demo Package |
| 95128093 | Reagent Reservoir, Disposable, 25 mL |

▶ Centrifugation

*Depend on
the full range
of superior
Thermo Scientific
solutions...*

Thermo Scientific general purpose benchtop centrifuges combine exceptional capacity and maximum throughput with proven ergonomic design and outstanding energy efficiency. Accelerate your sample preparation with Thermo Scientific innovative technologies, including:

- **ClickSeal® biocontainment rotor lids for one-handed, certified sample protection**
- **Secure, push-button Auto-Lock® rotor exchange for application versatility**
- **Broad rotor selection, including our lightweight Thermo Scientific Fiberlite® carbon fiber rotors, providing high performance and increased capacity**

From microcentrifuges and benchtop instruments to advanced floor models and innovative carbon fiber rotors – our centrifuge systems deliver outstanding performance and reliability in the lab.



Thermo Scientific Finntip Pipette Tips

Finntip* pipette tips are manufactured using high quality materials and the latest molding techniques for proven reliability.

Thermo Scientific Finntip pipette tips are designed for full compatibility with their matching FinnpiPETTE pipettors, providing a perfect sealing of the tip along with optimum pipetting accuracy and precision. Manufactured with highest quality materials in automated production facilities, Finntip pipette tips deliver a superior pipetting solution for discerning users.

The Finntip product line covers a wide volume range from 0.2 µL to 10 mL. Available offerings include standard pipette tips as well as macro volume, extended, wide orifice and filtered designs in various packaging options. Also available are sterile racked tips that are certified free from DNA, DNase, RNase and endotoxin.

Details

- Manufactured from highest quality virgin polypropylene raw materials
- Versatile product offering including macro volume, extended, wide orifice and filtered design
- Wide volume range of 0.2 µL to 10 mL
- Sterile product versions certified free from DNA, DNase, RNase and endotoxin available
- Color-coded racks to match the compatible FinnpiPETTE pipette

Pipetter Compatibility: Finntip pipette tips are optimized for use with FinnpiPETTE pipettors. The universal design enables the compatibility also with pipettors from other manufacturers.

Certifications: Manufactured according to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD certified).



| Cat. No. | Description | Packaging Type | Compatible FinnpiPETTES |
|---|-------------------------------|--------------------------------|--|
| Finntip 10, 0.2-10 µL, 3.2 cm, Pink | | | |
| 9400310 | Finntip 10 | 1000/bag | F1 (4641010, 4641020, 4641030, 4661000, 4661040), F2 (4642010, 4642020, 4642030, 4662000, 4662040), Novus (46200000, 46300000, 46300100) |
| 9400300 | Finntip 10 | 10 × 96/rack | |
| 9400303 | Finntip 10 sterile | 10 × 96/rack | |
| 9400326 | Finntip 10 refill starter kit | 4 × 192/rack + 10 × 192/refill | |
| 9400327 | Finntip 10 refill | 20 × 192/refill | |
| Finntip 20, 0.2-20 µL, 3.3 cm, Purple | | | |
| 9400620 | Finntip 20 | 1000/bag | F1 (4641010, 4641020, 4641030, 4641050, 4661000, 4661040, 4661080), F2 (4642010, 4642020, 4642030, 4642050, 4662000, 4662040, 4662080), Novus (46200000, 46300000, 46300100) |
| 9400610 | Finntip 20 | 10 × 384/rack | |
| 9400613 | Finntip 20 sterile | 10 × 384/rack | |
| Finntip 50, 0.2-50 µL, 4.8 cm, Turquoise | | | |
| 9400360 | Finntip 50 | 1000/bag | F1 (4641010, 4641020, 4641030, 4641050, 4641130, 4661000, 4661040, 4661080, 4661090), F2 (4642010, 4642020, 4642030, 4642050, 4642120, 4662000, 4662040, 4662080, 4662090), Novus (46200000, 46200200, 46300000, 46300100, 46300700) |
| 9400370 | Finntip 50 | 10 × 384/rack | |
| 9400373 | Finntip 50 sterile | 10 × 384/rack | |

Thermo Scientific Finntip Pipette Tips, continued

| Cat. No. | Description | Packaging Type | Compatible Finnpipettes |
|--|--|------------------------------|---|
| Finntip 200 Ext, 5-200 µL, 7.8 cm, Orange | | | |
| 9400100 | Finntip 200 Ext | 400/box | F1 (4641040, 4641060, 4641140, 4641070, 4641080, 4661010, 4661020, 4661050, 4661060), F2 (4642040, 4642060, 4642130, 4642070, 4642080, 4662010, 4662050, 4662020, 4662060), Novus (46200100, 46200300, 46200400, 46300200, 46300300) |
| 9400130 | Finntip 200 Ext | 10 × 96/rack | |
| 9400133 | Finntip 200 Ext sterile | 10 × 96/rack | |
| Finntip 250 Universal, 0.5-250 µL, 5.2 cm, Yellow | | | |
| 9400250 | Finntip 250 Universal | 500/box | F1 (4641040, 4641060, 4641140, 4641070, 4641080, 4661010, 4661020, 4661050, 4661060), F2 (4642040, 4642060, 4642130, 4642070, 4642080, 4662010, 4662050, 4662020, 4662060), Novus (46200100, 46200300, 46200400, 46300200, 46300300) |
| 9400230 | Finntip 250 Universal | 1000/bag | |
| 9400220 | Finntip 250 Universal | 20,000/bulk | |
| 9400260 | Finntip 250 Universal | 10 × 96/rack | |
| 9400263 | Finntip 250 Universal sterile | 10 × 96/rack | |
| 9400266 | Finntip 250 Universal refill starter kit | 4 × 96/rack + 10 × 96/refill | |
| 9400267 | Finntip 250 Universal refill | 20 × 96/refill | |
| Finntip 300, 5-300 µL, 5.2 cm, Orange | | | |
| 9401240 | Finntip 300 | 20,000/bulk | F1 (4641040, 4641060, 4641140, 4641070, 4641080, 4641090, 4661010, 4661020, 4661050, 4661060, 4661030, 4661070), F2 (4642040, 4642060, 4642130, 4642070, 4642080, 4662010, 4662050, 4662020, 4662060, 4662030, 4662070), Novus (46200100, 46200300, 46200400, 46200500, 46300200, 46300300, 46300400, 46300500) |
| 9401250 | Finntip 300 | 10 × 96/rack | |
| 9401253 | Finntip 300 sterile | 10 × 96/rack | |
| 9401255 | Finntip 300 refill | 20 × 96/refill | |
| 9401260 | Finntip 300 | 1000/bag | |
| Finntip 1000, 100-1000 µL, 7.1 cm, Blue | | | |
| 9401030 | Finntip 1000 | 1000/bag | F1 (4641100), F2 (4642090), Novus (46200600) |
| 9401010 | Finntip 1000 | 10,000/bulk | |
| 9401070 | Finntip 1000 | 200/box | |
| 9401110 | Finntip 1000 | 10 × 96/rack | |
| 9401113 | Finntip 1000 sterile | 10 × 96/rack | |
| 9401115 | Finntip 1000 refill | 10 × 96/refill | |
| Finntip 1000 Ext, 100-1000 µL, 10.5 cm, Blue | | | |
| 9401410 | Finntip 1000 Ext | 1000/bag | F1 (4641100), F2 (4642090), Novus (46200600) |
| 9401420 | Finntip 1000 Ext | 5 × 96/rack | |
| 9401423 | Finntip 1000 Ext sterile | 5 × 96/rack | |
| Finntip 5 mL, 0.5-5 mL, 14.7 cm, Green | | | |
| 9402050 | Finntip 5 mL | 75/box | F1 (4641110), F2 (4642100), Novus (46200700) |
| 9402030 | Finntip 5 mL | 500/bag | |
| 9402010 | Finntip 5 mL | 3000/bulk | |
| 9402070 | Finntip 5 mL | 5 × 54/rack | |
| 9402073 | Finntip 5 mL sterile | 5 × 54/rack | |
| Finntip 10 mL, 1-10 mL, 15.0 cm, Red | | | |
| 9402150 | Finntip 10 mL | 40/box | F1 (4641120), F2 (4642110), Novus (46200800) |
| 9402151 | Finntip 10 mL | 100/bag | |
| 9402160 | Finntip 10 mL | 5 × 24/rack | |
| 9402163 | Finntip 10 mL sterile | 5 × 24/rack | |



Thermo Scientific Finntip Extended Length Pipette Tips

Finntip Extended Length Pipette Tips are designed for use with extremely narrow or deep vessels.

Extended length tips allow you to access the bottom of test tubes, reagent bottles, flasks and other vessels without touching the barrel of your pipetter against the side of the tube.

The longer tip length of these tips allows you to reach the bottom of long or narrow vessels that standard tips cannot reach. The selection also includes an extended 10 mL tip, which is among the longest reaching tips available. The tip is very narrow and has been designed to reach to the bottom of common 1 L laboratory bottles and most volumetric flasks.

Details

- Prevents cross-contamination from vessel walls
- Allows reach to the bottom of particularly long or narrow vessels that standard tips cannot reach
- Versatile packaging options, with sterile and filtered versions available
- Chemical-resistant tips made of high-quality polypropylene
- Available in three volume ranges: 1-200 μ L, 100-1000 μ L and 1-10 mL

Pipetter Compatibility: Finntip pipette tips are optimized for use with Finnpiptette pipettors. The universal design enables compatibility also with pipettors from other manufacturers.

Certifications: Manufactured according to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD certified).



| Cat. No. | Description | Packaging Type | Compatible Finnpiptettes |
|--|--------------------------------|---------------------|--|
| Finntip 200 Ext, 5-200 μL, 7.8 cm, Orange | | | |
| 9400100 | Finntip 200 Ext | 400/box | F1 (4641040, 4641060, 4641140, 4641070, 4641080, 4661010, 4661020, 4661050, 4661060), F2 (4642040, 4642060, 4642130, 4642070, 4642080, 4662010, 4662050, 4662020, 4662060), Novus (46200100, 46200300, 46200400, 46300200, 46300300) |
| 9400130 | Finntip 200 Ext | 10 \times 96/rack | |
| 9400133 | Finntip 200 Ext sterile | 10 \times 96/rack | |
| Finntip 1000 Ext, 100-1000 μL, 10.5 cm, Blue | | | |
| 9401410 | Finntip 1000 Ext | 1000/bag | F1 (4641100), F2 (4642090), Novus (46200600) |
| 9401420 | Finntip 1000 Ext | 5 \times 96/rack | |
| 9401423 | Finntip 1000 Ext sterile | 5 \times 96/rack | |
| Finntip 10 mL Ext, 1-10 mL, 26.7 cm | | | |
| 94060970 | Finntip Flex 10 mL Ext | 100/bag | F1 (4641120), F2 (4642110), Novus (46200800) |
| 94060973 | Finntip Flex 10 mL Ext sterile | 50/bag | |

Thermo Scientific Finntip Wide Orifice Pipette Tips



Finntip Wide Orifice Pipette Tips are ideal for pipetting fragile cell suspensions and macromolecules including genomic DNA.

The wide orifice tips eliminate mechanical shearing that causes cell fragmentation. Depending on the model, the inside diameter of the orifice is 1.1 mm or 1.8 mm.

Details

- Designed for pipetting macromolecule
- Eliminates cell fragmentation caused by mechanical shearing
- Two different models covering a volume range of 10-1000 μ L
- Wide orifice tips with orifice inside diameters of 1.1 mm or 1.8 mm

Pipetter Compatibility: Finntip pipette tips are optimized for use with Finnpiptettes. The universal design enables compatibility also with pipettors from other manufacturers.

Certifications: Manufactured according to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD certified).

| Cat. No. | Description | Packaging Type | Compatible Finnpiptettes |
|---|---------------------------|---------------------|---|
| Finntip 250 Wide, 10-250 μL, 5.3 cm, White | | | |
| 9405020 | Finntip 250 Wide | 1000/bag | F1 (4641040, 4641060, 4641140, 4641070, 4641080, 4661010, 4661020, 4661050, 4661060), |
| 9405120 | Finntip 250 Wide | 10 \times 96/rack | F2 (4642040, 4642060, 4642130, 4642070, 4642080, 4662010, 4662050, 4662020, 4662060), |
| 9405123 | Finntip 250 Wide sterile | 10 \times 96/rack | Novus (46200100, 46200300, 46200400, 46300200, 46300300) |
| Finntip 1000 Wide, 100-1000 μL, 7.2 cm, White | | | |
| 9405050 | Finntip 1000 Wide | 400/bag | F1 (4641100), |
| 9405060 | Finntip 1000 Wide | 6500/bulk | F2 (4642090), |
| 9405160 | Finntip 1000 Wide | 10 \times 96/rack | Novus (46200600) |
| 9405163 | Finntip 1000 Wide sterile | 10 \times 96/rack | |



Thermo Scientific Finntip BioCon Pipette Tips

Finntip Biocon Pipette Tips are designed for applications that demand the utmost in biological purity.


Designed for applications that demand the highest levels of biological purity, the individually packed Finntip BioCon tips are guaranteed free of any biological contamination, including DNA, DNase, RNase and endotoxin. They are sterilized with irradiation, making them safe for medical uses (non-invasive), pharmaceutical and food industries, and PCR applications, as well as molecular biology and cell technology. A certificate of quality is provided in each box.

Details

- For applications that demand the highest biological purity
- Individually packed for maximum protection
- Available in three volume ranges: 0.2-10 µL, 0.5-250 µL and 100-1000 µL
- Certified free from DNA, DNase, RNase and endotoxin
- Certificate of quality is provided in each box

Pipetter Compatibility: Finntip pipette tips are optimized for use with Finnpiptette pipettors. The universal design enables compatibility also with pipettors from other manufacturers.

Certifications: Manufactured according to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD certified).



| Cat. No. | Description | Volume | Packaging Type |
|----------|---------------------|----------------|----------------|
| 94053100 | Finntip BioCon 10 | 0.2 to 10 µL | 100/box |
| 94053200 | Finntip BioCon 250 | 0.5 to 250 µL | 100/box |
| 94053300 | Finntip BioCon 1000 | 100 to 1000 µL | 50/box |



Thermo Scientific Finntip Filter Pipette Tips

Finntip Filter Pipette Tips are designed to protect your samples from contamination.

Finntip Filter Pipette Tips provide superior protection against contamination. These sterile tips are produced in a controlled environment, sterilized with irradiation and are certified free from DNA, DNase, RNase and endotoxin. Our filter tips are ideal for a wide range of applications, especially sensitive procedures including PCR based reactions.

Finntip Filter Pipette Tips are available in racks with an assortment of volumes from 0.2 µL to 10 µL. Made from an inert non-sealing polyethylene matrix, the filter is effective at eliminating contamination from the pipetter to the sample, thereby ensuring the integrity of any pipetting function.

Details

- Wide volume range from 0.2 µL to 10 µL
- Filters composed of inert non-sealing polyethylene matrix effectively eliminate carryover contamination from the pipetter to the sample
- Ideal for use in sensitive applications, including PCR
- 10 µL to 1000 µL product versions, certified free from DNA, DNase, RNase and endotoxin
- Supplied as vacuum-sealed sterilized tip racks

Pipetter Compatibility: Finntip pipette tips are optimized for use with Finnpiptette pipettors. The universal design enables compatibility also with pipettors from other manufacturers.

Certifications: Manufactured according to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD-certified)



| Cat. No. | Description | Volume | Length | Packaging Type |
|----------|-------------------------------------|-------------|---------|----------------|
| 94052000 | Finntip Filter Micro 10 µL, sterile | 0.2-10 µL | 3.2 cm | 10 × 96/rack |
| 94052020 | Finntip Filter Micro 20 µL, sterile | 0.2-20 µL | 3.3 cm | 10 × 384/rack |
| 94052060 | Finntip Filter Micro 50 µL, sterile | 0.2-50 µL | 4.8 cm | 10 × 384/rack |
| 94052100 | Finntip Filter 10 µL, sterile | 0.5-10 µL | 5.2 cm | 10 × 96/rack |
| 94052150 | Finntip Filter 20 µL, sterile | 0.5-20 µL | 5.2 cm | 10 × 96/rack |
| 94052160 | Finntip Filter 30 µL, sterile | 0.5-30 µL | 5.2 cm | 10 × 96/rack |
| 94052200 | Finntip Filter 100 µL, sterile | 0.5-100 µL | 5.2 cm | 10 × 96/rack |
| 94052310 | Finntip Filter 100 µL Ext, sterile | 5-100 µL | 7.8 cm | 10 × 96/rack |
| 94052300 | Finntip Filter 200 µL, sterile | 0.5-200 µL | 5.2 cm | 10 × 96/rack |
| 94052320 | Finntip Filter 200 µL Ext, sterile | 5-200 µL | 7.8 cm | 10 × 96/rack |
| 94052350 | Finntip Filter 300 µL, sterile | 5-300 µL | 5.2 cm | 10 × 96/rack |
| 94052410 | Finntip Filter 1000 µL, sterile | 100-1000 µL | 7.1 cm | 10 × 96/rack |
| 94052430 | Finntip Filter 1000 µL Ext, sterile | 100-1000 µL | 10.5 cm | 5 × 96/rack |
| 94052550 | Finntip Filter 5 mL, sterile | 0.5-5 mL | 14.7 cm | 5 × 54/rack |
| 94052600 | Finntip Filter 10 mL, sterile | 1-10 mL | 15.0 cm | 5 × 24/rack |



Thermo Scientific Finntip Flex Pipette Tips

Finntip Flex Pipette Tips offer a consistent seal for higher accuracy and precision.

Finntip Flex Pipette Tips are flexible, soft and sophisticated high-end tips for discerning users. All elements of this pipetting system have been designed to meet the most demanding requirements for ergonomics and flexibility in today's fast-paced laboratories.

Finntip Flex Pipette Tips are available in easy-to-use 10, 200, 300, 1000, 1200 μL and 10 mL sizes. The soft tip reduces the amount of force needed to attach and eject tips, significantly reducing the risk of Repetitive Strain Injury (RSI). It also ensures a better seal with the pipetter for greater accuracy and precision.

Details

- Excellent for single and multichannel pipetting
- Available in six volume ranges: 0.2-10 μL , 1-200 μL , 5-300 μL , 5-1000 μL , 50-1200 μL and 1-10 mL
- Sterile tips are supplied in vacuum sealed, sterilized, 10x96 rack and they are certified to be free from DNA, DNase, RNase and endotoxin
- Save space with refills

Softer Tip

- The soft tip reduces the amount of force needed to attach and eject tips, significantly lowering the risk of Repetitive Strain Injury (RSI)
- Provides a better seal with the pipetter for greater accuracy and precision

Slip-free Rack

- The see-through cover is easy to open with one hand, leaving the other hand free to perform other tasks
- The sturdy base and slip-proof feet keep the rack stationary so you can easily, yet firmly, attach the tip

Pipetter Compatibility: Finntip pipette tips are optimized for use with Finnpiptette pipettors. The universal design enables compatibility also with pipettes from other manufacturers.

Certifications: Manufactured according to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD certified)



Thermo Scientific Finntip Flex Pipette Tips, continued

| Cat. No. | Description | Packaging Type | Compatible Finnpipettes |
|---|--------------------------------------|------------------------------|---|
| Finntip Flex 10, 0.1-10 µL, 3.3 cm, Pink | | | |
| 94060120 | Finntip Flex 10 | 1000/bag | F1 (4641010, 4641020, 4641030, 4661000, 4661040), F2 (4642010, 4642020, 4642030, 4662000, 4662040), Novus (46200000, 46300000, 46300100) |
| 94060100 | Finntip Flex 10 | 10 × 96/rack | |
| 94060103 | Finntip Flex 10 sterile | 10 × 96/rack | |
| 94060116 | Finntip Flex 10 refill starter kit | 1 × 96/rack + 20 × 96/refill | |
| 94060117 | Finntip Flex 10 refill kit | 20 × 96/refill | |
| Finntip Flex 200, 1-200 µL, 6.0 cm, Yellow | | | |
| 94060320 | Finntip Flex 200 | 1000/bag | F1 (4641040, 4641060, 4641140, 4641070, 4641080, 4661010, 4661020, 4661050, 4661060), F2 (4642040, 4642060, 4642130, 4642070, 4642080, 4662010, 4662050, 4662020, 4662060), Novus (46200100, 46200300, 46200400, 46300200, 46300300) |
| 94060310 | Finntip Flex 200 | 10 × 96/rack | |
| 94060313 | Finntip Flex 200 sterile | 10 × 96/rack | |
| 94060316 | Finntip Flex 200 refill starter kit | 1 × 96/rack + 20 × 96/refill | |
| 94060317 | Finntip Flex 200 refill | 20 × 96/refill | |
| Finntip Flex 300, 5-300 µL, 6.15 cm, Orange | | | |
| 94060520 | Finntip Flex 300 | 1000/bag | F1 (4641040, 4641060, 4641140, 4641070, 4641080, 4641090, 4661010, 4661020, 4661050, 4661060, 4661030, 4661070), F2 (4642040, 4642060, 4642130, 4642070, 4642080, 4662010, 4662050, 4662020, 4662060, 4662030, 4662070), Novus (46200100, 46200300, 46200400, 46200500, 46300200, 46300300, 46300400, 46300500) |
| 94060510 | Finntip Flex 300 | 10 × 96/rack | |
| 94060513 | Finntip Flex 300 sterile | 10 × 96/rack | |
| 94060516 | Finntip Flex 300 refill starter kit | 1 × 96/rack + 20 × 96/refill | |
| 94060517 | Finntip Flex 300 refill | 20 × 96/refill | |
| Finntip Flex 1000, 50-1000 µL, 8.7 cm, Blue | | | |
| 94060720 | Finntip Flex 1000 | 1000/bag | F1 (4641100), F2 (4642090), Novus (46200600) |
| 94060710 | Finntip Flex 1000 | 10 × 96/rack | |
| 94060713 | Finntip Flex 1000 sterile | 10 × 96/rack | |
| 94060716 | Finntip Flex 1000 refill starter kit | 1 × 96/rack + 16 × 96/refill | |
| 94060717 | Finntip Flex 1000 refill | 16 × 96/refill | |
| Finntip Flex 1200, 50-1200 µL, 8.7 cm, Turquoise | | | |
| 94060820 | Finntip Flex 1200 | 1000/bag | F1 (4641100), F2 (4642090), Novus (46200600, 46300800) |
| 94060810 | Finntip Flex 1200 | 10 × 96/rack | |
| 94060813 | Finntip Flex 1200 sterile | 10 × 96/rack | |
| 94060816 | Finntip Flex 1200 refill starter kit | 1 × 96/rack + 16 × 96/refill | |
| 94060817 | Finntip Flex 1200 refill | 16 × 96/refill | |
| Finntip Flex 10 mL Ext, 1-10 mL, 26.7 cm | | | |
| 94060970 | Finntip Flex 10 mL Ext | 100/bag | F1 (4641120), F2 (4642110), Novus (46200800) |
| 94060973 | Finntip Flex 10 mL Ext sterile | 50/bag | |





Thermo Scientific Finntip Flex Filter Pipette Tips

Finntip Flex Filter Pipette Tips provide effective protection against contamination.

The Finntip Flex Filter Pipette Tips combine the ergonomics of a standard Finntip Flex pipette tip and the protection provided by a filtered tip. With the soft Flex tips, the amount of force needed to attach and eject tips is significantly reduced for enhanced pipetting ergonomics.

Finntip Flex Filter Pipette Tips effectively minimize the risk of passing contamination from the pipetter to the sample. Including a range of products certified free from DNA, DNase, RNase and endotoxin, these sterile tips suit all sensitive procedures, including PCR-based reactions where contamination is a concern. The filters are composed of an inert non-sealing polyethylene matrix and are effective at eliminating carryover contamination.

Details

- The soft tip reduces the amount of force needed to attach and eject tips, significantly lowering the risk of Repetitive Strain Injury (RSI)
- Filters made of inert non-sealing polyethylene matrix, effectively eliminate carryover contamination from the pipetter to the sample
- Ideal for use in sensitive applications, including PCR
- 10 µL to 1200 µL product versions certified-free from DNA, DNase, RNase and endotoxin
- Supplied as vacuum-sealed sterilized tip racks

Pipetter Compatibility: Finntip pipette tips are optimized for use with Finnpipette pipetters. The universal design enables compatibility also with pipetters from other manufacturers.

Certifications: Manufactured according to ISO 9001, ISO 14001 and ISO 13485 (CE/IVD certified).



| Cat. No. | Description | Volume | Length | Packaging Type |
|----------|--|-------------|---------|----------------|
| 94056980 | Finntip Flex Filter 10 µL, sterile | 0.2-10 µL | 3.3 cm | 10 × 96/rack |
| 94056510 | Finntip Flex Filter 30 µL, sterile | 1-30 µL | 6.0 cm | 10 × 96/rack |
| 94056520 | Finntip Flex Filter 100 µL, sterile | 1-100 µL | 6.0 cm | 10 × 96/rack |
| 94056380 | Finntip Flex Filter 200 µL, sterile | 1-200 µL | 6.0 cm | 10 × 96/rack |
| 94056580 | Finntip Flex Filter 300 µL, sterile | 5-300 µL | 6.15 cm | 10 × 96/rack |
| 94056710 | Finntip Flex Filter 1000 µL, sterile | 100-1000 µL | 8.7 cm | 10 × 96/rack |
| 94056810 | Finntip Flex Filter 1200 µL, sterile | 50-1200 µL | 8.7 cm | 10 × 96/rack |
| 94056970 | Finntip Flex Filter 10 mL Ext, sterile | 1-10 mL | 26.7 cm | 50/bag |

Thermo Scientific Finntip Stepper Tips



Finntip Stepper Tips are the perfect fit with Thermo Scientific Finnpipette Stepper Pipeters.

The Finntip Stepper Tips are plunger-style tips for positive displacement pipetting. Seven sizes are available in sterile and non-sterile versions.

Details

- Plunger-style tips
- Overall dispensing range of 10 to 5000 µL
- Nonsterile tips are bulk packaged; sterile tips are individually wrapped

Certifications: ISO 9001 and ISO 14001

| Volume Setting Chart for Finntip Stepper Tips | | | | | |
|---|---------|---------|---------|---------|---------|
| Wheel Position | 1 | 2 | 3 | 4 | 5 |
| Maximum No. of Deliveries | 44 | 22 | 15 | 11 | 9 |
| Volume/Delivery with | | | | | |
| 0.5 mL Tip | 10 µL | 20 µL | 30 µL | 40 µL | 50 µL |
| 1.25 mL Tip | 25 µL | 50 µL | 75 µL | 100 µL | 125 µL |
| 2.5 mL Tip | 50 µL | 100 µL | 150 µL | 200 µL | 250 µL |
| 5.0 mL Tip | 100 µL | 200 µL | 300 µL | 400 µL | 500 µL |
| 12.5 mL Tip | 250 µL | 500 µL | 750 µL | 1000 µL | 1250 µL |
| 25 mL Tip | 500 µL | 1000 µL | 1500 µL | 2000 µL | 2500 µL |
| 50 mL Tip | 1000 µL | 2000 µL | 3000 µL | 4000 µL | 5000 µL |

| Cat. No. | Description | Packaging Type |
|--|-----------------------------------|----------------|
| Finntip Stepper 0.5 mL, dispensing volume: 10, 20, 30, 40, 50 µL | | |
| 9404170 | Finntip Stepper, 0.5 mL | 100/box |
| 9404173 | Finntip Stepper, 0.5 mL, sterile | 50/box |
| Finntip Stepper 1.25 mL, dispensing volume: 25, 50, 75, 100, 125 µL | | |
| 9404180 | Finntip Stepper, 1.25 mL | 100/box |
| 9404183 | Finntip Stepper, 1.25 mL, sterile | 50/box |
| Finntip Stepper 2.5 mL, dispensing volume: 50, 100, 150, 200, 250 µL | | |
| 9404190 | Finntip Stepper, 2.25 mL | 100/box |
| 9404193 | Finntip Stepper, 2.25 mL, sterile | 50/box |
| Finntip Stepper 5.0 mL, dispensing volume: 100, 200, 300, 400, 500 µL | | |
| 9404200 | Finntip Stepper, 5.0 mL | 50/box |
| 9404203 | Finntip Stepper, 5.0 mL, sterile | 25/box |
| Finntip Stepper 12.5 mL, dispensing volume: 250, 500, 750, 1000, 1250 µL | | |
| 9404210 | Finntip Stepper, 12.5 mL | 50/box |
| 9404213 | Finntip Stepper, 12.5 mL, sterile | 25/box |
| Finntip Stepper 25 mL, dispensing volume: 500, 1000, 1500, 2000, 2500 µL | | |
| 9404220 | Finntip Stepper, 25 mL | 20/box |
| 9404223 | Finntip Stepper, 25 mL, sterile | 10/box |
| Finntip Stepper 50 mL, dispensing volume: 1000, 2000, 3000, 4000, 5000 µL | | |
| 9404230 | Finntip Stepper, 50 mL | 10/box |
| 9404233 | Finntip Stepper, 50 mL, sterile | 10/box |

Thermo Scientific Finntip Multistepper Pipette Tips

Finntip Multistepper Pipette Tips optimize the performance of Thermo Scientific Finnpipette Multistepper Pipetters.


Specifically designed for use with the Thermo Scientific Finnpipette Multistepper Pipetter, this tip has a nominal volume of 1500 μL and accommodates standard 9 mm microplate spacing. Finntip Multistepper Tips are available in racks of 96 tips, in both sterile and non-sterile versions.

Details

- For use with Thermo Scientific Finnpipette Multistepper Pipetter
- Length: 8.8 cm
- Volume: 1500 μL

Certifications: ISO 9001 and ISO 14001

| Cat. No. | Description | Packaging Type |
|--|-------------------------------|---------------------|
| Finntip Multistepper, 1500 μL , 8.5 cm, Turquoise | | |
| 9401300 | Finntip Multistepper | 400/box |
| 9401330 | Finntip Multistepper | 10 \times 96/rack |
| 9401333 | Finntip Multistepper, sterile | 10 \times 96/rack |



Thermo Scientific Matrix Single-Channel Electronic Pipettors



Matrix* Single-Channel Electronic Pipettors combine intuitive step-based programming with superior ergonomic design to provide power and flexibility for enhanced pipetting.

Step-based programming allows you to perform pipetting routines that would either be cumbersome or virtually impossible with manual pipettors. What used to take many motions of the thumb and manual adjustment of a volume knob can now be accomplished with the light touch of a button.

This ergonomically-designed pipettor performs all the tasks you might do with a manual pipettor – only faster and with less effort. Step-based programming makes it logical, intuitive and easy to create programs, while onboard memory saves five programs with up to 40 steps each.

Details

- **Step-based programming** enables logical, intuitive and easy to create programs
- **Onboard memory** saves five programs (up to 40 steps each), allowing you to spend more time pipetting and less time programming
- **Individual speed control:** five variable speeds are available for each step of your procedure (aspirate, dispense or mix)
- **Trigger-based operation** and tip ejection eliminate use of the thumb and reduces strain when ejecting tips
- **Five volume ranges** are available, increasing the flexibility to perform everything from micro to large volume pipetting
- **Operates while plugged in**, so there is no downtime during recharging
- **Light tip application and ejection:** Designed in conjunction with Matrix pipette tips for superior fit with minimal application/ejection forces

Includes: AC power supply, warranty card, instruction manual and calibration report

Warranty: One year

| Cat. No. | Volume | Accuracy* | Precision** | Color Code |
|--|-------------|------------------|------------------|------------|
| 1029 | 0.5-12.5 µL | ±1.0% or 0.05 µL | 0.4% or 0.05 µL | Red |
| 1020 | 1.0-30 µL | ±1.0% or 0.15 µL | 0.3% or 0.05 µL | Purple |
| 1021 | 2-125 µL | ±0.6% or 0.3 µL | 0.2% or 0.1 µL | Yellow |
| 1022 | 5-250 µL | ±0.6 or 0.5 µL | 0.15% or 0.15 µL | Blue |
| 1024 | 15-1250 µL | ±0.5% or 3.0 µL | 0.13% or 0.6 µL | Green |
| Thermo Scientific Matrix Single Channel Short Barrel Pipettor | | | | |
| 1122 | 5-250 µL | ±0.6% or 0.5 µL | 0.15% or 0.15 µL | Blue |
| 1124 | 15-1250 µL | ±0.5% or 3.0 µL | 0.13% or 0.6 µL | Green |

Note: For Accuracy and Precision, values shown are expressed as a percent (%) deviation or microliter (µL) value. When applied to desired volume, the greater of the two values will always apply.

**Both values represent the deviation from the mean.*

***Percentage values are expressed as the coefficient of variation.*

Microliter values are expressed as the standard deviation.

Thermo Scientific Matrix Multichannel Electronic Pipetters

Matrix Multichannel Electronic Pipetters deliver the power of electronic pipetting in microplate formats.

In addition to a well-balanced, ergonomic design, these pipetters feature step-based programming that allows you to perform pipetting routines that would be cumbersome or virtually impossible with manual pipetters. What once took many thumb motions and manual volume adjustments can now be accomplished with the light touch of a button. Available in 8-, 12-, or 16-channel models for performing a range of operations – from microvolume pipetting in a PCR plate to filling microplates.

Details

- **Step-based programming** for logical, intuitive, and easy to create programs
- **Onboard memory** saves five programs (up to 40 steps each), allowing you to spend more time pipetting and less time programming
- **Individual speed control:** Five variable speeds are available for each step of your procedure (aspirate, dispense or mix)
- **Large volume capacity** of 1250 μL per tip allows you to fill an entire microplate with up to 100 μL per well from a single aspiration
- **Lithium-ion battery technology** allows more time working, less time charging; no memory effect

Design

- **Trigger-based operation and tip ejection** eliminates use of the thumb
- **Unique design** keeps your hand close to the pipetting surface; tip ejector utilizes two or three fingers, not the thumb
- **Light tip application and ejection:** Designed in conjunction with Matrix pipette tips for superior fit requiring minimal application/ejection forces

Step-Based Programming

- Easy-to-use feature allows you to automate anything you would do with a manual multichannel pipette
- All functions (aspirate, dispense, mix) can be linked together (up to 40 steps) and saved for future use; up to five programs can be stored in memory
- Example: Aspirate 1200 μL , dispense 50 μL to the first two rows, then 100 μL to the next ten rows while setting specific dispense speeds for each step
- Alternatively, aspirate 20 μL and then 200 μL of diluent, then repeat per each row of a plate for an in-plate serial dilution
- “Scratch Pad” programming mode allows fast and simple program changes that are not stored in memory
- Paced-dispense function dispenses the programmed volume at the specified pace (fast or slow) automatically

Includes: Universal power supply, warranty card, instruction manual and calibration report

Warranty: One year



Thermo Scientific Matrix Multichannel Electronic Pipettors, continued



| Cat. No. | Volume | Channels | Accuracy* | Precision** | Color Code |
|----------|------------------------|----------|-----------------------------------|----------------------------|------------|
| 2009 | 0.5-12.5 μL | 8 | $\pm 2.5\%$ or 0.15 μL | 2.0% or 0.15 μL | Red |
| 2001 | 2-125 μL | 8 | $\pm 2.0\%$ or 1.0 μL | 1.0% or 0.60 μL | Yellow |
| 2002 | 5-250 μL | 8 | $\pm 2.0\%$ or 1.5 μL | 0.7% or 1.0 μL | Blue |
| 2004 | 15-1250 μL | 8 | $\pm 1.5\%$ or 6.0 μL | 0.6% or 3.0 μL | Green |
| 2019 | 0.5-12.5 μL | 12 | $\pm 2.5\%$ or 0.15 μL | 2.0% or 0.15 μL | Red |
| 2011 | 2-125 μL | 12 | $\pm 2.0\%$ or 1.0 μL | 1.0% or 0.60 μL | Yellow |
| 2012 | 5-250 μL | 12 | $\pm 2.0\%$ or 1.5 μL | 0.7% or 1.00 μL | Blue |
| 2014 | 15-850 μL | 12 | $\pm 1.5\%$ or 4.0 μL | 0.6% or 2.50 μL | Orange |
| 2069 | 0.5-12.5 μL | 16 | $\pm 2.5\%$ or 0.15 μL | 2.0% or 0.15 μL | Red |
| 2060 | 1.0-30 μL | 16 | $\pm 2.0\%$ or 0.30 μL | 1.7% or 0.25 μL | Purple |
| 2061 | 2-125 μL | 16 | $\pm 2.0\%$ or 1.0 μL | 1.0% or 0.60 μL | Yellow |

Note: For accuracy and precision, values shown are expressed as a percent (%) deviation or microliter (μL) value. When applied to desired volume, the greater of the two values will always apply.

*Both values represent the deviation from the mean.

**Percentage values are expressed as the coefficient of variation. Microliter values are expressed as the standard deviation.

Thermo Scientific Matrix EXP Electronic Pipetters

Matrix EXP Electronic Pipetters rapidly increase the throughput of sample transfers from tube racks to microplates.

Our innovative expandable tip spacing system increases productivity with ease. Simply pull the rod to expand the tips for accessing test tube or microcentrifuge tube racks, 24 or 48 well plates – and then push the rod to compress the spacing for pipetting into 96 well microplates. This allows you to transfer up to eight samples at once, and saves many arm/thumb motions when compared to performing this routine with single-channel pipetters.

Easy to use, step-based programming allows you to automate anything you would do with a manual multichannel pipetter.

Details

- **Step-based programming** for logical, intuitive, and easy to create programs
- **Onboard memory** saves five programs (up to 40 steps each), allowing you to spend more time pipetting and less time programming
- **Individual speed control:** Five variable speeds are available for each step of your procedure (aspirate, dispense, or mix)
- **Trigger-based operation and tip ejection** eliminates use of the thumb and reduces strain when ejecting tips
- **Unique design** keeps your hand close to the pipetting surface; tip ejector utilizes two or three fingers, not the thumb
- **Large volume capacity** of 1250 µL per tip allows you to fill an entire microplate with up to 100 µL per well from a single aspiration
- **Lithium-ion battery** technology allows more time working, less time charging; no memory effect
- **Operates while plugged in**, so there is no downtime during recharging

Expandable Tip Spacing

- Expandable tip spacing allows you to transfer multiple samples between different tube rack and microplate configurations
- The pipette tips can be expanded to access test tube racks, 24-well plates, or 48-well plates – or closed to work with 96-well plates
- Also operates as a standard pipetter for 24- or 48-well plates

Includes: Universal power supply, warranty card, instruction manual and calibration report

Warranty: One year



| Cat. No. | Volume | Channels | Accuracy* | Precision** | Color Code | Spacing (mm) | Microplate format | Tube format |
|----------|-------------|----------|------------------|-----------------|------------|--------------|-------------------|----------------------|
| 2621 | 2-125 µL | 6 | ±2.0% or 1.0 µL | 1.0% or 0.6 µL | Yellow | 9.00-19.81 | 96/24-well | 16 mm |
| 2622 | 5-250 µL | 6 | ±2.0% or 1.5 µL | 0.7% or 1.0 µL | Blue | 9.00-19.81 | 96/24 | 16 mm |
| 2624 | 15-1250 µL | 6 | ±1.5% or 6.0 µL | 0.6% or 3.0 µL | Green | 9.00-19.81 | 96/24-well | 16 mm |
| 2029 | 0.5-12.5 µL | 8 | ±2.5% or 0.15 µL | 2.0% or 0.15 µL | Red | 9.00-13.06 | 96-well | microcentrifuge |
| 2021 | 2-125 µL | 8 | ±2.0% or 1.0 µL | 1.0% or 0.6 µL | Yellow | 9.00-14.15 | 96/48-well | 12 or 13 mm |
| 2022 | 5-250 µL | 8 | ±2.0% or 1.5 µL | 0.7% or 1.0 µL | Blue | 9.00-14.15 | 96/48-well | 12 or 13 mm |
| 2024 | 15-1250 µL | 8 | ±1.5% or 6.0 µL | 0.6% or 3.0 µL | Green | 9.00-14.15 | 96/48-well | 12 or 13 mm |
| 2229 | 0.5-12.5 µL | 12 | ±2.5% or 0.15 µL | 2.0% or 0.15 µL | Red | 6.00-9.00 | 96-well | Terasaki, gels, 864s |

Note: For Accuracy and Precision, values shown are expressed as a percent (%) deviation or microliter (µL) value. When applied to desired volume, the greater of the two values will always apply.

**Both values represent the deviation from the mean. **Percentage values are expressed as the coefficient of variation. Microliter values are expressed as the standard deviation.*

Thermo Scientific Matrix Equalizer Electronic Multichannel Pipetters



Matrix Equalizer Electronic Multichannel Pipetters speed sample transfer between labware with different configurations.

The Matrix Electronic Pipette with equalizer tip spacing increases productivity while reducing the risk of Repetitive Strain Injury (RSI). The tip spacing mechanism expands and contracts with a simple slide adjustment to perform sample transfers between virtually any tube, rack, microplate or horizontal gel box.

Details

- Flexibility to **pipette between many types of labware** with five volume and three tip spacing ranges
- **Step-based programming** makes it logical, intuitive and easy to create programs
- Onboard memory saves five programs (up to 40 steps each), allowing you to spend more time pipetting and less time programming
- Choice of five speeds independently controls aspiration, dispense and mix steps
- Ergonomic design with **trigger-based operation** and tip ejection increases comfort and reduces the risk of repetitive stress injuries
- **Lithium-ion battery** technology allows more time working, less time charging; no memory effect
- Five volume ranges and three spacing ranges allow the flexibility to pipette across many types of labware
- Designed in conjunction with Matrix pipette tips, for superior tip fit with minimal application and ejection forces

Adjustable Equal Tip Spacing

- Tips can be spaced anywhere from 9 mm (96-well microplates) to 14.15 mm (12-13 mm tube rack spacing or 48-well plates) in increments of 0.1 mm
- Distance between tips is always equal – simply set the desired spacing, then pull or push the rod to expand or contract the tips
- Spacing can be set while tips are empty or full, enabling quick transfer of multiple samples at once between labware with different configurations
- Matrix Equalizer 384 models (0.5-12.5, 1-30 and 2-125 μ L) feature the same unique tip spacing mechanism but allow you to set the tip spacing down to 4.5 mm – perfect for multichannel dispensing from 96-well into 384-well plates or for multi-lane sample loading of agarose gels. The 8-channel pipetter expands to 14.15 mm; the 12-channel to 9 mm

Applications:

- Plate-to-plate transfers between 96- and 384-well microplates
- Multichannel sample transfers between tube racks and microplates
- Multichannel gel loading directly from microplates or tube racks
- Sample additions
- Serial dilutions

Includes: Universal power supply, warranty card, instruction manual and calibration report

Warranty: One year

| Cat. No. | Model | Range | Channels | Accuracy* | Precision** | Color Code | Spacing (mm) | Microplate format |
|----------|---------------|------------------|----------|----------------------------|----------------------|------------|--------------|-------------------|
| 2032 | Equalizer | 5-250 μ L | 8 | \pm 2.0% or 1.5 μ L | 0.7% or 1.0 μ L | Blue | 9.0-14.2 | 96/48-well |
| 2034 | Equalizer | 15-1250 μ L | 8 | \pm 1.5% or 6.0 μ L | 0.6% or 3.0 μ L | Green | 9.0-14.2 | 96/48-well |
| 2139 | Equalizer 384 | 0.5-12.5 μ L | 8 | \pm 2.5% or 0.15 μ L | 2.0% or .15 μ L | Red | 4.5-14.2 | 384/96/48-well |
| 2130 | Equalizer 384 | 1.0-30 μ L | 8 | \pm 2.0% or 0.30 μ L | 1.7% or 0.25 μ L | Purple | 4.5-14.2 | 384/96/48-well |
| 2131 | Equalizer 384 | 2-125 μ L | 8 | \pm 2.0% or 1.0 μ L | 1.0% or 0.6 μ L | Yellow | 4.5-14.2 | 384/96/48-well |
| 2239 | Equalizer 384 | 0.5-12.5 μ L | 12 | \pm 2.5% or 0.15 μ L | 2.0% or .15 μ L | Red | 4.5-9.0 | 384/96-well |
| 2230 | Equalizer 384 | 1.0-30 μ L | 12 | \pm 2.0% or 0.30 μ L | 1.7% or 0.25 μ L | Purple | 4.5-9.0 | 384/96-well |
| 2231 | Equalizer 384 | 2-125 μ L | 12 | \pm 2.0% or 1.0 μ L | 1.0% or 0.60 μ L | Yellow | 4.5-9.0 | 384/96-well |

Note: For accuracy and precision, values shown are expressed as a percent (%) deviation or microliter (μ L) value. When applied to desired volume, the greater of the two values will always apply.

*Both values represent the deviation from the mean.

**Percentage values are expressed as the coefficient of variation.

Microliter values are expressed as the standard deviation.



NEW




F1-ClipTip Pipetting System

You'll feel the difference the first time you hold the **NEW** Thermo Scientific F1-ClipTip Pipetting System.

Achieve newfound confidence knowing that once attached, your tips are locked firmly in place, enabling consistent, reproducible pipetting for higher quality results. Your tips will not loosen or fall off regardless of application pressure. Transform your daily pipetting with an airtight seal on every channel for security you can feel.


breakthrough ClipTip technology

FRICION SEALING SYSTEMS



- Increased tip attachment and ejection forces
- Loose tips which can drop off or leak air
- Decreased confidence in reproducibility of results
- Increased Tip Cone Wear

CLIPTIP SYSTEM



- Extremely low tip attachment & ejection forces
- Complete seal on every channel
- No loose tips
- Reproducibility you can count on

See it in action! www.thermoscientific.com/cliptip



*Depend on
the full range
of superior
Thermo Scientific
solutions...*

**Worry-free pipetting is finally here.
With the only tip that clips.**

Thermo Scientific Matrix Serological Pipet Filler



The Matrix Serological Pipet Filler is a lightweight, cordless pipet filler that allows longer, fatigue-free pipetting.

The Matrix Serological Pipet Filler is an intuitive, lightweight, cordless pipet filler for use with glass or plastic serological pipets. Easy and comfortable to operate, it features a heavy-duty rechargeable lithium-ion battery that is guaranteed for life and allows the pipet filler to operate approximately three to four times longer between charges than typical units using nickel cadmium batteries. Available in a choice of five colors, the Matrix Serological Pipet Filler can be wall mounted or placed in its table stand. It includes a 1 mL pipet support to eliminate “wobble” when using 1 mL pipets.

Details

- Large backlit LCD display provides visual confirmation of remaining battery charge and speed settings
- Separate aspirate and dispense speed controls provide precise control over pipetting speeds and reduce risk of overpipetting with smaller pipets
- Rechargeable lithium-ion battery comes with lifetime guarantee, and offers long runtime per charge while eliminating battery memory problems
- Battery has rapid charge rate of 80% in one hour and delivers up to 20 hours of continuous operation before recharge
- Nosepiece can accommodate pipets from 1 to 100 mL; included 1 mL pipet stabilizer prevents 1 mL pipets from wobbling or coming loose from nose piece
- Comes with weighted table stand and a convenient wall-mountable holder

Separate Aspirate and Dispense Speed Controls

- Buttons on the rear of the instrument allow the individual setting of aspiration (upper button) and dispense (lower button) speeds
- Eight speeds are available for each function: Simply push for faster or slower
- Selection displayed on the LCD screen
- A zero speed selection on dispense enables gravity dispensing, and extra slow aspiration speeds prevent overpipetting when using 1 mL pipets

Includes: Matrix Serological Pipet Filler, universal charger, table stand, wall-mounted holder, 1 mL pipet support, instruction manual and warranty card

Warranty: One year

Specifications

| | |
|---------------------------------|--|
| Aspirate Speeds | 8 |
| Dispense Speeds | 8 (+ Gravity Dispense) |
| Battery | Lithium-ion |
| Approx. Runtime Between Charges | 15 Continuous Hours |
| Recharge Time | 3 hrs. |
| Weight | 220 g |
| Nosepiece Pipet Holder | Autoclavable Silicone Rubber |
| Pipette Compatibility | All major brands of glass or plastic serological pipettes, 1 to 100 mL |
| Filter Assembly | 0.45 μ Standard |

| Cat. No. | Color |
|----------|--|
| 9531 | Matrix Serological Pipet Filler, Red |
| 9541 | Matrix Serological Pipet Filler, Green |
| 9501 | Matrix Serological Pipet Filler, White |
| 9511 | Matrix Serological Pipet Filler, Clear |
| 9521 | Matrix Serological Pipet Filler, Blue |

Thermo Scientific Matrix Pipetter Accessories

Matrix Pipetter Accessories are optimized for use with Matrix Electronic Pipettors and Matrix Serological Pipet Fillers.

Pipetter Stands for Matrix Electronic Pipettors

- **Single stand:** Free-standing metal unit holds one Matrix electronic single or multichannel pipetter
- **Three-position wall rack:** For Matrix electronic single or multichannel pipettors; includes mounting hardware
- **Six-position revolving stand:** Ideal for labs using multiple pipettors; rotates freely to keep your bench top organized

Tube Racks

- Helps you take full advantage of the expandable tip system on a Matrix EXP or Equalizer electronic pipettors
- Two different tube racks are available for simple reformatting from tubes to microplates
- Test tube racks work with either 12 × 75 mm or 13 × 100 mm test tubes in an 8 × 12 or 5 × 12 format
- Microcentrifuge tube rack works with either 0.5 mL or 1.5 mL tubes

Matrix Serological Pipet Filler Accessories

- **The wall-mounted holder and table stand** offer compact storage between applications. It is possible to recharge the pipet filler in both stands by plugging in the charger. The wall-mounted holder is excellent for storing pipet filler when working in laminar flow cabinets to save space. The wing stand** allows setting the Matrix Serological Pipet Filler down without removing or emptying the serological pipet.
- **The easy-to-change filter** prevents aerosol or liquid cross contamination of the pipet together with an autoclavable silicone adapter. The 1 mL Pipet Support ensures that 1 mL pipets will hold as securely as larger sizes.



Single Pipet Filler Wall-mounted Holder



1 mL Pipet Support



Matrix Single Pipette Stand

| Cat. No. | Description | Qty. |
|--|---|---------|
| Pipetter Stands for Matrix Electronic Pipettors | | |
| 1063 | Wall Mountable Matrix Pipette Stand, 3-position | 1 |
| 1066 | Revolving Carousel Matrix Pipette Stand, 6-position | 1 |
| 8066 | Matrix Single Pipette Stand | 1 |
| Tube Racks | | |
| 8813 | Sample Transfer Tube Rack, Autoclavable, 12 × 75 mm and 13 × 100 mm tubes, 8 × 12 | 1 |
| 8814 | Sample Transfer Tube Rack, Autoclavable, 0.5 mL or 1.5 mL Microcentrifuge Tubes, 8 × 12 | 1 |
| Matrix Serological Pipet Filler Accessories | | |
| 9066 | Single Pipet Filler Wall-mounted Holder | 1 |
| 9067 | Single Pipet Table Stand | 1 |
| 9069 | Wing Stand for Serological Pipet Filler** | 1 |
| 9057 | Hydrophobic Filters, 0.45 µm | 25/case |
| 9064 | Nosepiece (holder, gripper, filter) | 1 |
| 9065 | Silicone Pipet Gripper | 4/case |
| 9070 | 1 mL Pipet Support | 1 |
| ** Not available for use or sale within the United States or for importation into the United States. | | |

Thermo Scientific Nunc Serological Pipets


Nunc* Serological Pipets are accurate, disposable plastic pipets calibrated to deliver to the tip.

Sterilized and plugged for convenience, these pipets are an excellent choice for cell culture applications and have a certified Sterility Assurance Level (SAL) of 10^{-6} .

Details

- Volumes are color-coded for quick size identification, and packaging is color-coded for ease in sorting and selecting the correct size
- Includes "reverse" and "negative" graduations in easy-to-read black printed scale
- Convenient extra graduations to full pipet volume
- Top ends are beveled with cotton plugging
- Nontoxic, certified nonpyrogenic – suitable for cell culture work
- Bold clear graduations with easy-to-read scale markings
- Individually wrapped in convenient peel-open paper bags
- Printed with black scale
- Polystyrene
- Sterile

Ordering Information: Supplied sterile, individually packaged in easy-open paper/plastic wrapper. Overpacked in convenient dispensing shelf pack.



| Cat. No. | Capacity | Graduation | Neg. Grad. | Tolerance | Color Code | Qty./Bag |
|----------|----------|------------|------------|-----------|------------|----------|
| 159609 | 1 | 0.01 | 0.3 | ±0.02 | Black | 200 |
| 159617 | 2 | 0.01 | 0.3 | ±0.04 | Black | 125 |
| 159625 | 5 | 0.1 | 2 | ±0.10 | Blue | 50 |
| 159633 | 10 | 0.1 | 3 | ±0.20 | Orange | 50 |
| 159641 | 25 | 0.2 | 10 | ±0.50 | Green | 50 |
| 159668 | 50 | 0.5 | 10 | ±0.75 | Purple | 25 |



Thermo Scientific Matrix Pipette Tip Compatibility Chart

| | Vol. (µL) | Electronic Single Channel | | | | | Electronic Multichannel (6-, 8- and 12-Channel) Fixed, Equalizer, EXP | | | | | Electronic Multichannel 16-Channel and 384-Equalizer | | | Hybrid Single Channel | | | | | Hybrid Multichannel 8-Channel and 12-Channel | | | | | |
|--|------------------------|---------------------------|------|------|------|------|---|-----|-----|-----|------|--|--|-----|-----------------------|----|-----|-----|------|--|----|-----|---|---|----------------------------------|
| | | 12.5 | 30 | 125 | 250 | 1250 | 12.5 | 125 | 250 | 850 | 1250 | 12.5 | 30 | 125 | 12.5 | 30 | 125 | 300 | 1250 | 12.5 | 30 | 300 | | | |
| Matrix Pipette Tips | Cat. No. | 1029 | 1020 | 1021 | 1022 | 1024 | 2002 2019 2001 2012 2004 2009 2011 2022 2024 2029 2021 2032 2034 2229 2621 2622 2014 2624 | | | | | | 2069 2130 2131 2139 2230 2231 2239 2060 2061 | | | | | | | | | | | | 2310 2312 2315 2320 2322 2325 |
| 30 µL TallTip | 7600, 7631, 7632 | • | • | | | | • | | | | | | | | | | | | | | | | | | |
| 12.5 µL Filter TallTip | 7635 | • | • | | | | • | | | | | | | | | | | | | | | | | | |
| 250 µL Tip | 7250, 7151, 7152, 7156 | | | • | • | | • | • | | | | | | | | | | | | | | | | | |
| 200 µL Filter Tip | 7275 | | | • | • | | • | • | | | | | | | | | | | | | | | | | |
| 30 µL Filter Tip | 7155 | | A | | | | | | | | | | | | | | | | | | | | | | |
| 250 µL TallTip | 7280, 7281, 7282 | | | • | • | | • | • | | | | | | | | | | | | | | | | | |
| 200 µL Filter TallTip | 7285 | | | • | • | | • | • | | | | | | | | | | | | | | | | | |
| 300 µL Tip | 7320, 7321, 7322 | | | • | • | | • | • | | | | | | | | | | | | | | | | | |
| 250 µL Filter Tip | 7325 | | | • | • | | • | • | | | | | | | | | | | | | | | | | |
| 300 µL TallTip | 7080, 7081, 7082 | | | • | • | | • | • | | | | | | | | | | | | | | | | | |
| 250 µL Filter TallTip | 7085 | | | • | • | | | | | | | | | | | | | | | | | | | | |
| 1250 µL Tip | 8050, 8041, 8042, 8046 | | | | | • | | | • | • | | | | | | | | | | | | | | | |
| 1250 µL Filter Tip | 8045 | | | | | • | | | • | • | | | | | | | | | | | | | | | |
| 1250 µL TallTip | 8250, 8241, 8242, 8246 | | | | | • | | | • | • | | | | | | | | | | | | | | | |
| 1250 µL Filter TallTip | 8245 | | | | | • | | | • | • | | | | | | | | | | | | | | | |
| 12.5 µL 384 Tip | 7421, 7422 | | | | | | | | | | | • | | | | | | | | | | | | | |
| 12.5 µL 384 Filter Tip | 7425 | | | | | | | | | | | • | | | | | | | | | | | | | |
| 30 µL 384 Tip | 7431, 7432 | | | | | | | | | | | | • | | | | | | | | | | | | |
| 30 µL 384 Filter Tip | 7435 | | | | | | | | | | | | • | | | | | | | | | | | | |
| 125 µL 384 Tip | 7441, 7442 | | | | | | | | | | | | | • | | | | | | | | | | | |
| 125 µL 384 Filter Tip | 7445 | | | | | | | | | | | | | • | | | | | | | | | | | |
| 12.5 µL/30 µL Matrix Hybrid ClipTip | 7120, 7121, 7122 | | | | | | | | | | | | | • | • | | | | | | • | • | | | |
| 12.5 µL Filter Matrix Hybrid ClipTip | 7165 | | | | | | | | | | | | | • | | | | | | | • | | | | |
| 30 µL Filter Matrix Hybrid ClipTip | 7175 | | | | | | | | | | | | | | • | | | | | | | • | | | |
| 125 µL/300 µL Matrix Hybrid ClipTip | 7130, 7131, 7132 | | | | | | | | | | | | | | | • | • | | | | | | | • | |
| 125 µL/300 µL Filter Matrix Hybrid ClipTip | 7135 | | | | | | | | | | | | | | | • | • | | | | | | | • | |
| 1250 µL Matrix Hybrid ClipTip | 7140, 7141, 7142 | | | | | | | | | | | | | | | | | | | | | | • | | |
| 1250 µL Filter Matrix Hybrid ClipTip | 7145 | | | | | | | | | | | | | | | | | | | | | | • | | |

A = Requires special tip ejector (included with pipetter)

Thermo Scientific Matrix Pipette Tips

Matrix Pipette Tips offer a comprehensive selection of tips supplied in innovative Flo-Thru racks.

Matrix pipettors and pipette tips have been designed together to ensure that a perfect balance exists between the security of the seal and ease of application/ejection. Manufactured under tight QC specifications, these standard tips rival special low-retention tips and are capable of pipetting fluids, such as dH₂O, 50% glycerol, hybridized buffer and 0.1% SDS.

Details

- Sealing rings on tips ensure confident, leak-free seal without banging tips on the pipetter
- These tips also offer a secure seal on most other pipettors brands available, due to the compression qualities of the sealing rings
- Highly polished molds and rigorous maintenance procedures enable low liquid retention tips without special treatments or additives

Innovative Racks and Packaging

- Rigid rack design ensures solid, even tip seating; rack will not flex or bend when applying tips
- Unique Flo-Thru rack individually supports tips to prevent them from wobbling when trying to seat them on the pipetter
- Hinge prevents sterile lid from coming into contact with the lab bench
- Venting along the edge or around each tip improves heat flow during autoclaving, and reduces condensation that can form on the tips
- Dual-position lid can be propped open to prevent vacuum deformation during autoclaving
- Racks are constructed of clean recycled materials and are made as small as possible to reduce waste
- Available in various packaging formats, including the environmentally friendly ECOTIPS refill systems

Pipetter compatibility: Fits Thermo Scientific Matrix Electronic Pipettors



| Cat. No. | Description | Packaging Type |
|---|--------------------------|----------------|
| 10 µL Matrix Pipette Tips, 0.2-10 µL, 3.1 cm | | |
| 7610 | 10 µL Tip | 1000/bag |
| 7611 | 10 µL Tip | 10 × 96/rack |
| 7612 | 10 µL Tip, sterile | 10 × 96/rack |
| 12.5 µL Matrix Pipette Tips, 0.5-12.5 µL, 4.1 cm | | |
| 7421 | 12.5 µL 384 Tip | 10 × 384/rack |
| 7422 | 12.5 µL 384 Tip, sterile | 10 × 384/rack |
| 20 µL Matrix Pipette Tips, 0.5-20 µL, 4.6 cm | | |
| 7620 | 20 µL Tip | 1000 tips/bag |
| 7621 | 20 µL Tip | 10 × 96/rack |
| 30 µL Matrix Pipette Tips, 1-30 µL, 4.8 cm | | |
| 7431 | 30 µL 384 Tip | 10 × 384/rack |
| 7432 | 30 µL 384 Tip, sterile | 10 × 384/rack |
| 125 µL Matrix Pipette Tips, 2-125 µL, 5.4 cm | | |
| 7441 | 125 µL 384 Tip | 10 × 384/rack |
| 7442 | 125 µL 384 Tip, sterile | 10 × 384/rack |
| 250 µL Matrix Pipette Tips, 2-250 µL, 4.9 cm | | |
| 7250 | 250 µL Tip | 1000 tips/bag |
| 7151 | 250 µL Tip | 10 × 96/rack |
| 7152 | 250 µL Tip, sterile | 10 × 96/rack |
| 7156 | 250 µL ECOTIPS | 10 × 96/refill |
| 300 µL Matrix Pipette Tips, 2-300 µL, 5.8 cm | | |
| 7320 | 300 µL Tip | 1000 tips/bag |
| 7321 | 300 µL Tip | 10 × 96/rack |
| 7322 | 300 µL Tip, sterile | 10 × 96/rack |
| 1250 µL Matrix Pipette Tips, 15-1250 µL, 7.9 cm | | |
| 8050 | 1250 µL Tip | 500 tips/bag |
| 8041 | 1250 µL Tip | 10 × 96/rack |
| 8042 | 1250 µL Tip, sterile | 10 × 96/rack |
| 8046 | 1250 µL ECOTIPS | 10 × 96/refill |



Thermo Scientific Matrix TallTip Extended Length Pipette Tips



Matrix TallTip Pipette Tips provide extended length and a narrow shaft for maximum reach.

Access the bottom of test tubes, reagent bottles, flasks and other vessels without touching the barrel of your pipetter against the side of the tube.

The longer length of Matrix TallTip Extended Length tips allows you to reach the bottom of long or narrow vessels that standard tips cannot reach. Non-beveled, narrow tip ends provide cleaner touch-off of sidewalls for greater precision.

Details

- Prevents cross-contamination from vessel walls
- Enables the use of large volume tips in small spaces, resulting in less liquid displacement
- Non-beveled, narrow tip ends provide cleaner touch-off of sidewalls for greater precision
- Allows you to reach the bottom of particularly long or narrow vessels that standard tips cannot reach
- 102 mm maximum length offers safe access to 100 mm test tubes
- Flo-Thru rack reduces tip wobble and improves airflow while autoclaving
- Universal fit works with most major brands of pipetters

Pipetter compatibility: Fits Thermo Scientific Matrix Electronic Pipetters



| Cat. No. | Description | Packaging Type |
|---|-------------------------------|-----------------------|
| 30 μL Matrix TallTip, 0.5-30 μL, 4.1 cm | | |
| 7600 | 30 μ L TallTip | 1000/bag |
| 7631 | 30 μ L TallTip | 10 \times 96/rack |
| 7632 | 30 μ L TallTip, sterile | 10 \times 96/rack |
| 250 μL Matrix TallTip, 2-250 μL, 7.8 cm | | |
| 7280 | 250 μ L TallTip | 500/bag |
| 7281 | 250 μ L TallTip | 10 \times 96/rack |
| 7282 | 250 μ L TallTip, sterile | 10 \times 96/rack |
| 300 μL Matrix TallTip, 2-300 μL, 10.2 cm | | |
| 7080 | 300 μ L TallTip | 500/bag |
| 7081 | 300 μ L TallTip | 10 \times 96/rack |
| 7082 | 300 μ L TallTip, sterile | 10 \times 96/rack |
| 1250 μL Matrix TallTip, 15-250 μL, 10.2 cm | | |
| 8250 | 1250 μ L TallTip | 500/bag |
| 8241 | 1250 μ L TallTip | 10 \times 96/rack |
| 8242 | 1250 μ L TallTip, sterile | 10 \times 96/rack |
| 8246 | 1250 μ L TallTip ECOTIPS | 10 \times 96/refill |



Thermo Scientific Matrix Filter Pipette Tips

Matrix Filter Pipette Tips are rigorously tested to deliver optimum performance.

Just because a tip has a filter doesn't necessarily mean that it is effective in blocking aerosols. Matrix filter tips are tested for many types of potential filter failure, and every lot must pass a rigorous battery of tests. If accidental contact with the filter should occur due to over-pipetting, your sample can be recovered by simply dispensing back out. Because the filter is hydrophobic, the sample will remain uncontaminated and will not be trapped in the filter.

Whether you're pipetting genomic material, bacterial samples or other precious samples, the risk of cross-contamination from aerosols created during pipetting is always a concern. To make certain no contamination of the pipetter barrel occurs, the use of filter tips is recommended.

Details

- Stringent testing guarantees filter tips that effectively block aerosols
- Hydrophobic polyethylene filter material ensures no particulates contaminate samples
- Inadvertent contact with filter will not contaminate sample or result in sample loss

Stringent Filter Testing

- Every filter must meet certain porosity requirements and must form a complete barrier across face of tip

Recommended for:

- Bacteriology
- Gel loading
- PCR setup
- Preventing cross-contamination during routine pipetting

Pipetter Compatibility: Thermo Scientific Matrix Electronic Pipetters

| Cat. No. | Description | Volume | Length | Packaging Type |
|----------|---------------------------------|-------------|--------|----------------|
| 7615 | 10 µL Filter Tip, sterile | 0.2-10 µL | 3.1 cm | 10 × 96/rack |
| 7425 | 12.5 µL 384 Filter Tip, sterile | 0.5-12.5 µL | 4.1 cm | 10 × 384/rack |
| 7155 | 30 µL Filter Tip, sterile | 1-30 µL | 4.9 cm | 10 × 96/rack |
| 7435 | 30 µL 384 Filter Tip, sterile | 1-30 µL | 4.8 cm | 10 × 384/rack |
| 7445 | 125 µL 384 Filter Tip, sterile | 2-125 µL | 5.4 cm | 10 × 384/rack |
| 7275 | 200 µL Filter Tip, sterile | 2-200 µL | 5.4 cm | 10 × 96/rack |
| 7325 | 250 µL Filter Tip, sterile | 2-250 µL | 5.8 cm | 10 × 96/rack |
| 8045 | 1250 µL Filter Tip, sterile | 15-1250 µL | 7.9 cm | 10 × 96/rack |



Thermo Scientific Matrix TallTip Filter Pipette Tips

Matrix TallTip Filter Pipette Tips provide extended length for maximum reach and contamination protection for your samples and pipetters.

Extended-length tips allow you to access the bottom of test tubes, reagent bottles, flasks and other vessels without touching the barrel of your pipette against the side of the tube.

The longer tip length of Matrix TallTip Extended Length tips allows you to reach the bottom of long or narrow vessels that standard tips cannot reach. Non-beveled, narrow tip ends provide cleaner touch-off of sidewalls for greater precision. The hydrophobic filters effectively block aerosols and ensure that no particulates contaminate samples.



Details

- Prevents cross-contamination from vessel walls
- Allows the use of large volume tips in small spaces, resulting in less liquid displacement
- Non-beveled, narrow tip ends provide cleaner touch-off of sidewalls for greater precision
- Reach the bottom of particularly long or narrow vessels that standard tips cannot reach
- 102 mm maximum length offers safe access to 100 mm test tubes
- Flo-Thru rack reduces tip wobble and improves airflow while autoclaving
- Universal fit works with most major brands of pipetters

Hydrophobic Filter

- Stringent testing guarantees filter tips that effectively block aerosols
- Hydrophobic filter material ensures no particulates contaminate samples
- Easy-to-open packaging helps product to maintain sterility while being easy to access
- Hinged rack lids can be opened without laying lid on bench top, maintaining sterility

Pipetter Compatibility: Fits Thermo Scientific Matrix Electronic Pipetters

| Cat. No. | Description | Volume | Length | Packaging Type |
|----------|-------------------------------------|-------------|---------|----------------|
| 7635 | 12.5 µL TallTip Filter Tip, sterile | 0.5-12.5 µL | 4.1 cm | 10 × 96/rack |
| 7285 | 200 µL TallTip Filter Tip, sterile | 2-200 µL | 7.8 cm | 10 × 96/rack |
| 7085 | 250 µL TallTip Filter Tip, sterile | 2-250 µL | 10.2 cm | 10 × 96/rack |
| 8245 | 1250 µL TallTip Filter Tip, sterile | 15-1250 µL | 10.2 cm | 10 × 96/rack |



The Matrix Electronic Pipette with equalizer tip spacing increases productivity while reducing the risk of Repetitive Strain Injury (RSI).



Thermo Scientific Matrix Reagent Reservoirs

Matrix Reagent Reservoirs offer a unique, innovative design to accommodate a variety of applications.

These sterile reservoirs offer an extra-wide base to add rigidity and stability, helping to avoid spills. The 25 mL divided reservoir allows the pipetting of two different reagents with up to eight channels on one side, and up to four on the other, maximizing recovery if using only a single channel and a small amount of reagent.

Details

- Trough-within-a-trough maximizes amount of liquid accessible to pipette tips when using small amounts of reagent
- Pour-off spouts on all four corners reduce spillage when pouring reagents out of reservoir
- Graduations on inside wall enable quick measurement of remaining liquid
- Extra-wide base adds rigidity and stability to reservoir, helping to avoid spills
- Sterile reservoirs are packaged in snap-and-tear bags for sterility and easy opening
- 25 mL divided reservoir allows the pipetting of two different reagents with up to eight channels on one side, and up to four on the other; maximizes recovery if using only single channel and small amount of reagent



| Cat. No. | Capacity | Packaging | Mfr. No. |
|-----------------------------------|------------------------|----------------|----------|
| Sterile Polystyrene | | | |
| 8093 | 25 mL | 10 Bags of 10 | 8093 |
| 8093EA | 25 mL | | 8093EA |
| 8094 | 25 mL | Indiv. wrapped | 8094 |
| 8095 | 25 mL, with divider | 10 Bags of 10 | 8095 |
| 8096 | 25 mL, with divider | Indiv. wrapped | 8096 |
| 8085 | 100 mL | 10 Bags of 10 | 8085 |
| 8086 | 100 mL | Indiv. wrapped | 8086 |
| Autoclavable Polypropylene | | | |
| 8075 | 75 mL, without divider | 100 per case | 8075 |
| Reservoir Lids | | | |
| 8076 | 75 mL | 10 per case | 8076 |

Thermo Scientific Matrix Memowell Pipetting Aid

The Matrix Memowell Pipetting Aid is a 96-well LED light box that helps you to keep track of pipetting in a microplate, reducing errors.

This pipetting aid is ideally suited for tracking the pipetting of colorless reagents or minute liquid volumes, and marking progress when interrupted or when moving between sample tubes and a 96-well plate. An LED light is illuminated under specific wells, rows or columns of a microplate, which can be used to keep track of specific areas of a plate when filling, or indicate wells that need to be aspirated. To advance the light(s) to the next well or row, simply press a button on the unit, or use the optional foot switch for hands-free operation.

Details


- Microplate tracking reduces errors when filling microplates or performing plate-to-plate transfers; also facilitates hit-picking
- Contains six pre-programmed sequences; six additional programs can be accommodated
- Works manually or with optional MemoControl software for even greater tracking capabilities
- Plate reader results can be used to control lighting of specific columns, rows, or individual wells via software
- Multiple Memowell units can be linked together to keep track of plate-to-plate transfers
- Onboard batteries allow cordless operation without a power supply
- Optional foot switch control enables hands-free operation

MemoControl Software

- Microsoft® Excel plug-in; used in conjunction with a Windows® PC and plate reader
- Allows you to program Memowell to illuminate rows, columns or individual wells based on user-defined plate readings
- Manipulate reader data in Excel before using Memowell
- LED display can be programmed to indicate all wells containing hits with a light that is on, off or blinking
- The software supports the illumination of a four-well set within a 384-well microplate, which can be targeted and displayed on a PC monitor

Recommended for: Hit-picking, tracking plate filling, tracking plate-to-plate transfers

| Cat. No. | Description |
|----------|-------------------------------|
| 5000 | Memowell Pipetting Aid |
| 5008 | MemoControl Software |
| 6201 | Foot Switch Control |
| 5002 | Rechargeable Batteries |
| 5009 | PC Connection Cable, 9/25 Pin |



**Depend on
the full range
of superior
Thermo Scientific
solutions...**

▶ Blood Banking

We offer a complete range of quality products to support the key steps in the blood banking process – from blood separation, sample preparation and testing to storage, cryopreservation and blood related research. Our portfolio includes a wide choice of centrifuges, cold storage equipment, biological safety cabinets, media and reagents, and much more.

Rely on proven Thermo Scientific solutions and application expertise to meet your critical blood banking requirements for:

- **Capacity and reliability**
- **Security and safety**
- **Product tracking and compliance**

Learn more at www.thermoscientific.com/bloodbanking



Processing

A broad range of low speed centrifuges for the separation and preparation of whole blood and blood products



Sample Preparation

A broad range of versatile centrifuges and accessories for the efficient preparation of blood samples



Testing

An array of instruments, equipment and consumables to facilitate blood type testing, serological screening and nucleic acid testing



Storage

A comprehensive range of cold storage equipment and consumables for storing and safeguarding valuable blood products



Research

High-quality laboratory equipment and cell culture consumables for researchers working with blood products

Molecular BioProducts and Pure Pipette Tips with MicroPoint Tip Design


Molecular BioProducts* Tips with MicroPoint* design offer a tip profile that is significantly smaller and narrower than standard tips.

Molecular BioProducts and Pure Pipette Tips with MicroPoint Tip Design provide precise delivery, especially when working with low-volume samples. This is the result of a tip profile that is nearly 50 percent smaller than standard blunt-end tips and significantly more narrow than beveled-end tips. Crystal clear tips allow you to easily see transparent samples.

Details

- Same superior quality as ART* tips, except without the barrier
- Crystal clear tips allow you to easily see transparent samples
- Reference marks on many Molecular BioProducts Tips provide a quick and convenient guide to visually establish whether you are pipetting the desired volume of liquid
- All Pure tips are racked, pre-sterilized and certified free of RNase, DNase, DNA and pyrogen; other Molecular BioProducts racked and bulk tips are nonsterile and autoclavable. Pure 10 tips are ideal for low volume samples
- The extended length Pure 10 REACH* and Pure 1000 REACH prevent contamination of samples, tubes, and pipetters
- Ultra Micro and Pure Ultra Micro are also well suited for low volume samples
- 200 and Pure 200 are equivalent to the "Universal Yellow" tip and have a reference mark at 10 μL

| Cat. No. | Description | Volume | Packaging |
|----------|------------------|------------------------|-----------------|
| 3501 | Pure 10 | 0.1-10 μL | Racked, Sterile |
| 3502 | 10 | 0.1-10 μL | Racked |
| 3500 | 10 | 0.1-10 μL | Bulk |
| 3511 | Pure 10 Reach | 0.1-10 μL | Racked, Sterile |
| 3512 | 10 Reach | 0.1-10 μL | Racked |
| 3510 | 10 Reach | 0.1-10 μL | Bulk |
| 3521 | Pure Ultra Micro | 0.5-20 μL | Racked, Sterile |
| 3522 | Ultra Micro | 0.5-20 μL | Racked |
| 3520 | Ultra Micro | 0.5-20 μL | Bulk |
| 3551 | Pure 200 | 0.5-200 μL | Racked, Sterile |
| 3552 | 200 | 0.5-200 μL | Racked |
| 3550 | 200 | 0.5-200 μL | Bulk |
| 3541 | Pure XLP | 5-200 μL | Racked, Sterile |
| 3542 | XLP | 5-200 μL | Racked |
| 3540 | XLP | 5-200 μL | Bulk |
| 3531 | Pure 200G | 40-200 μL | Racked, Sterile |
| 3532 | 200G | 40-200 μL | Racked |
| 3530 | 200G | 40-200 μL | Bulk |
| 3571 | Pure 300 | 1-300 μL | Racked, Sterile |
| 3572 | 300 | 1-300 μL | Racked |
| 3570 | 300 | 1-300 μL | Bulk |
| 3561 | Pure 300E | 1-300 μL | Racked, Sterile |
| 3562 | 300 | 1-300 μL | Racked |
| 3591 | Pure 1000G | 100-1000 μL | Racked, Sterile |
| 3592 | 1000G | 100-1000 μL | Racked |
| 3590 | 1000G | 100-1000 μL | Bulk |
| 3791 | Pure 1000 Reach | 100-1000 μL | Racked, Sterile |
| 3792 | 1000 Reach | 100-1000 μL | Racked |



Notes: HLT tips are available in yellow and blue colors, graduated, beveled tips and large volumes.

Pipetter Compatibility: Universal fit for research-grade pipetters



Molecular BioProducts and Pure Pipette Tips with SoftFit Tip Design



Molecular BioProducts and Pure Pipette Tips feature SoftFit* design, which allows the tips to conform securely to the barrel of the pipetter.

Designed with thinner walls than ordinary pipette tips, SoftFit Tips conform more easily to the pipetter barrel, enabling a better seal and allowing the tip to eject with less force. Available as sterile or nonsterile offerings.

Details

- Same superior quality as ART tips, except without the barrier
- Crystal clear tips allow you to easily see transparent samples
- Reference marks on many Molecular BioProducts tips provide a quick and convenient guide to visually establish whether you are pipetting the desired volume of liquid
- All Pure tips are racked, pre-sterilized and certified free of RNase, DNase, DNA and pyrogen; other Molecular BioProducts racked and bulk tips are nonsterile and autoclavable
- SoftFit design provides universal fit and easy loading and ejecting; ideal for multichannel pipetting

Pipetter Compatibility: Universal fit for research-grade pipettes

| Cat. No. | Description | Volume | Packaging |
|----------|-------------|-------------|-----------------|
| 3771 | Pure 300U | 1.0-300 µL | Racked, sterile |
| 3772 | 300U | 1.0-300 µL | Racked |
| 3581 | Pure 1000 | 100-1000 µL | Racked, sterile |
| 3582 | 1000 | 100-1000 µL | Racked |
| 3580 | 1000 | 100-1000 µL | Bulk |



Molecular BioProducts Pure Pipette Tips with SoftFit L Design

Molecular BioProducts SoftFit L Pipette Tips are designed for use with Rainin® LTS® LiteTouch® Pipettors.

SoftFit L Tips have a positive stop feature that fits exclusively to Rainin LTS pipettors.

Designed with thinner walls than ordinary pipette tips, SoftFit L tips conform more easily to the pipettor barrel, providing a better seal and allowing the tip to eject with less force. Available as sterile or nonsterile offerings.

Details

- Molecular BioProducts tips are nonsterile; Pure are sterile (see below)
- Available in lift-off or hinged racks.

Hinged Racks

- Covers can be used as hinged or lift-off style
- Easy, one-hand operation of locking clasp
- Pipette tip inserts clip securely into racks but are easy to remove
- Empty rack part numbers allow for easy conversion to reload system
- All components are completely recyclable and contain recycling symbols
- Economical and environmental design reduces waste and saves bench space

Pipettor Compatibility: Rainin LTS LiteTouch Pipettor



| Cat. No. | Volume | Packaging |
|---------------------------|---------|-----------------------|
| MBP Tips | | |
| 3722 | 20 µL | Racks, 10 × 96 |
| 3722-HR | 20 µL | Hinged Racks, 10 × 96 |
| 3752 | 200 µL | Racks, 10 × 96 |
| 3752-HR | 200 µL | Hinged Racks, 10 × 96 |
| 3732 | 300 µL | Racks, 10 × 96 |
| 3732-HR | 300 µL | Hinged Racks, 10 × 96 |
| 3782 | 1000 µL | Racks, 8 × 96 |
| 3782-HR | 1000 µL | Hinged Racks, 8 × 96 |
| 3742 | 1200 µL | Racked, 8 × 96 |
| 3742-HR | 1200 µL | Hinged Racks, 8 × 96 |
| Pure Tips, Sterile | | |
| 3721 | 20 µL | Racks, 10 × 96 |
| 3721-HR | 20 µL | Hinged Racks, 10 × 96 |
| 3731 | 200 µL | Racks, 10 × 96 |
| 3751 | 200 µL | Racks, 10 × 96 |
| 3751-HR | 200 µL | Hinged Racks, 10 × 96 |
| 3731 | 300 µL | Racks, 10 × 96 |
| 3731-HR | 300 µL | Hinged Racks, 10 × 96 |
| 3781 | 1000 µL | Racks, 8 × 96 |
| 3781-HR | 1000 µL | Hinged Racks, 8 × 96 |
| 3741 | 1200 µL | Racks, 8 × 100 |
| 3741-HR | 1200 µL | Hinged Racks, 8 × 96 |

Molecular BioProducts ART* Barrier Tips with MicroPoint Design

Molecular BioProducts ART Tips are pre-sterilized and precision-molded, with the ART self-sealing barrier.

Ideal for use in genetic studies, forensics, PCR and radioisotope sampling, ART barrier tips with MicroPoint design are pre-sterilized and precision-molded. The ART self-sealing barrier prevents carryover contamination during aspiration and delivery of samples.

Details

- Pre-sterilized and precision-molded tips with the ART self-sealing barrier
- Certified RNase-, DNase-, and pyrogen-free
- Self-sealing barrier prevents carryover contamination during aspiration and delivery of samples
- MicroPoint design minimizes surface area around the tip orifice for reduced surface tension and greater accuracy
- Ideal for use in genetic studies, forensics, PCR and radioisotope sampling
- Available in numerous sizes and styles to meet all of your liquid handling needs

Pipetter Compatibility: Universal fit for research-grade pipettors

Note: Standard ART tips are also available in bulk and individually-wrapped packaging.



| Cat. No. | ART Tip | Volume | Packaging |
|---|-------------------------|-------------|---|
| Standard ART Tips with MicroPoint Design | | | |
| 2155P | Gel 20P, Gel Loading | 0.5-20 µL | Racked, sterile |
| 2155 | Gel 100, Gel Loading | 0.5-100 µL | Racked, sterile |
| 2065 | 100 | 0.5-100 µL | Racks, 10 × 96, Sterile |
| 2065-HR | 100 | 0.5-100 µL | Hinged racks, 10 × 96, Sterile |
| 2065-RI | 100 | 0.5-100 µL | Reload insert |
| 2065E | 100E | 1-100 µL | Racks, 10 × 96, Sterile |
| 2065E-HR | 100E | 1-100 µL | Hinged racks, 10 × 96, Sterile |
| 2065E-RI | 100E | 1-100 µL | Reload insert |
| 2069 | 200 | 1-200 µL | Racks, 10 × 96 Sterile |
| 2069-HR | 200 | 1-200 µL | Hinged Racks, 10 × 96, Sterile and 2069-RI, 200 1-200 µL, Reload Insert, 10 × 96 |
| 2159P | XLP, Ext. Length | 5-180 µL | Racks, 8 × 96, Sterile |
| 2160P | XLP 200, Ext. Length | 5-200 µL | Racks, 8 × 96, Sterile |
| 2069G | 200G, Genomic | 5-200 µL | Racks, 10 × 96, Sterile |
| 2070 | 300 | 1-300 µL | Racks, 10 × 96, Sterile |
| 2079E | 1000E | 100-1000 µL | Racks, 8 × 100 Sterile |
| 2079G | 1000G, Genomic | 100-1000 µL | Racks, 8 × 100, Sterile |
| 2079 | 1000 Reach, Ext. Length | 100-1000 µL | Racks, 8 × 100, Sterile |
| 2179-HR | 1000XL | 100-1000 µL | Hinged Racks, 8 × 96, Sterile and 2179-RI, 1000XL, 100-1000 µL, Reload Insert, 8 × 96 |
| 2180B | 5000 | 1-5 mL | Bag, 250 |
| 2149 | 20 | 0.1-20 µL | Racks, 10 × 96, sterile |
| 2149P | 20P | 0.1-20 µL | Racks, 10 × 96, Sterile |
| 2149P-HR | 20P | 1-20 µL | Hinged Racks, 10 × 96, Sterile and 2149P-RI, 20P, 1-20 µL, Reload Insert, 10 × 96 |

Molecular BioProducts ART Barrier Ultra Micro Tips with MicroPoint Design

Molecular BioProducts ART Tips are pre-sterilized and precision-molded, with the ART self-sealing barrier.

Certified RNase-, DNase-, and pyrogen-free, ART Ultra Micro Barrier Tips prevent carryover contamination during aspiration and delivery of samples. They are perfect for low-volume applications in genetic studies, forensics, PCR and radioisotope sampling.

Details

- Guaranteed to prevent carryover contamination during aspiration and delivery of samples
- Certified RNase-, DNase-, and pyrogen-free
- Ideal for low-volume applications in genetic studies, forensics, PCR and radioisotope sampling

Tip Types:

- The ART 10 and 10F have a 2 μL reference mark and are perfect for pipetting low volumes
- The extended tip length of the ART 10 REACH* prevents contamination of samples, tubes, and pipettors
- The ART 20E is an ultra-micro tip design for Eppendorf* and Nichiryo/BenchMate* pipettors
- The ART 20P has a 10 μL reference mark and is well suited for low volumes


Note: Standard ART tips are also available in bulk and individually-wrapped packaging.



| Cat. No. | ART Tip | Volume | Packaging |
|----------|----------|------------------|---------------------------------------|
| 2139 | 10 | 10 μL | Racked, 10 \times 96 |
| 2139-HR | 10 | 10 μL | Hinged Racks, 10 \times 96, Sterile |
| 2139-RI | 10 | 10 μL | Reload Insert, 10 \times 96 |
| 2140 | 10 Reach | 10 μL | Racked, 10 \times 96 |
| 2140-HR | 10 Reach | 10 μL | Hinged racks, 10 \times 96 |
| 2140-RI | 10 Reach | 10 μL | Reload insert |
| 2139F | 10F | 10 μL | Racked, 10 \times 96 |
| 2149E | 20E | 10 μL | Racked, 10 \times 96 |
| 2149E-HR | 20E | 10 μL | Hinged racks, 10 \times 96 |
| 2149E-RI | 20E | 10 μL | Reload insert |
| 2149P | 20P | 20 μL | Racked, 10 \times 96 |
| 2149P-HR | 20P | 20 μL | Hinged racks, 10 \times 96 |
| 2149P-RI | 20P | 20 μL | Reload insert |



Molecular BioProducts Low Retention Pipette Tips with ART Self-Sealing Barrier



Molecular BioProducts Low Retention Pipette Tips increase sample reproducibility and accuracy.

The extremely hydrophobic inner surface of Molecular BioProducts Low Retention Pipette Tips increases sample reproducibility and sample accuracy. Certified RNase-, DNase- and pyrogen-free, they are ideal for use in genetic studies, forensics, PCR and radioisotope sampling.

Details

- Specifically designed for eliminating sample retention within the pipette tip
- Extremely hydrophobic inner surface increases sample reproducibility and sample accuracy
- Pre-sterilized and precision-molded tips with the ART self-sealing barrier
- Prevent carryover contamination during aspiration and delivery of samples
- Certified RNase-, DNase- and pyrogen-free
- Ideal for use in genetic studies, forensics, PCR and radioisotope sampling
- SoftFit-L tips are compatible with LTS pipetters made by Rainin Instrument, LLC
- Racked, pre-sterilized

Note: ART Low Retention Tips are available in hinged racks and reload inserts. Please see www.thermoscientific.com for a full list of available products.


Pipette compatability: SoftFit L tips are compatible with LTS pipetters made by Rainin Instrument, LLC. All other tips have universal fit for research grade pipetters.



| Cat. No. | ART Tip | Volume | Packaging |
|-------------|-----------------------|-------------|---------------------------------|
| 2139-05 | 10 | 0.1-10 µL | Rack, 10 × 96 |
| 2139-05-HR | 10 | 10 µL | Hinged Racks, 10 × 96 |
| 2139-05-RI | 10 | 10 µL | Reload Insert, 10 × 96 |
| 2140-05 | 10 Reach*, MicroPoint | 0.1-10 µL | Rack, 10 × 96 |
| 2140-05-HR | 10 Reach, MicroPoint | 0.1-10 µL | Hinged Rack, 10 × 96 |
| 2140-05-RI | 10 Reach, MicroPoint | 0.1-10 µL | Reload insert |
| 2139F-05 | 10F, MicroPoint | 0.5-10 µL | Rack, 10 × 96 |
| 2149P-05 | 20P, MicroPoint | 0.5-20 µL | Rack, 10 × 96 |
| 2149P-05-HR | 20P, MicroPoint | 0.5-20 µL | Hinged Racks, 10 × 96 |
| 2149P-05-RI | 20P, MicroPoint | 0.5-20 µL | Reload insert |
| 2149E-05 | 20E, Ultra Micro | 0.5-10 µL | Rack, 10 × 96 |
| 2149E-05-HR | 20E, Ultra Micro | 0.5-10 µL | Hinged Racks, 10 × 96 |
| 2149E-05-RI | 20E, Ultra Micro | 0.5-10 µL | Hinged Racks, Reload insert |
| 2749-05 | 20L, SoftFit-L | 0.5-10 µL | Rack, 10 × 96 |
| 2749-05-HR | 20L, SoftFit-L | 0.5-10 µL | Hinged Rack, 10 × 96 |
| 2749-05-RI | 20L, SoftFit-L | 0.5-10 µL | Reload insert |
| 2065E-05 | 100E, MicroPoint | 0.1-100 µL | Rack, 10 × 96 |
| 2065E-05-HR | 100E MicroPoint | 1-100 µL | Hinged Rack, 10 × 96 |
| 2065E-05-RI | 100E MicroPoint | 1-100 µL | Reload insert, 10 × 96, Sterile |
| 2069-05 | 200, MicroPoint | 1.0-200 µL | Rack, 10 × 96 |
| 2069-05-HR | 200, MicroPoint | 200 µL | Hinged Rack, 10 × 96 |
| 2069-05-RI | 200, MicroPoint | 200 µL | Reload insert |
| 2769-05 | 200L, SoftFit L | 200 µL | Rack, 10 × 96 |
| 2769-05-HR | 200L, SoftFit L | 200 µL | Hinged Rack, 10 × 96 |
| 2769-05-RI | 200L, SoftFit L | 200 µL | Reload insert |
| 2739-05 | 300L, SoftFit L | 300 µL | Rack, 10 × 96 |
| 2739-05-HR | 300L, SoftFit L | 300 µL | Hinged Rack, 10 × 96 |
| 2739-05-RI | 300L, SoftFit L | 300 µL | Reload insert |
| 2279-05 | 1000, SoftFit | 100-1000 µL | Rack, 10 × 96 |
| 2179-05-HR | 1000XL | 100-1000 µL | Hinged Rack, 8 × 96 |
| 2179-05-RI | 1000XL | 100-1000 µL | Reload insert |
| 2779-05-HR | 1000L, SoftFit L | 100-1000 µL | Hinged Rack, 8 × 96 |
| 2779-05-RI | 1000L, SoftFit L | 100-1000 µL | Reload insert |
| 2779-05 | 1000L, SoftFit L | 100-1000 µL | Hinged Rack, 8 × 96 |
| 2789-05 | 1200L, SoftFit L | 100-1200 µL | Rack, 8 × 96 |
| 2789-05-HR | 1200L, SoftFit L | 1200 µL | Hinged Rack, 8 × 96 |
| 2789-05-RI | 1200L, SoftFit L | 1200 µL | Reload insert |



Molecular BioProducts Low Retention Pipette Tips



Molecular BioProducts Low Retention Pipette Tips are specifically designed for eliminating sample retention within the pipette tip.

The inner surface of the Molecular BioProducts Low Retention Pipette tip is extremely hydrophobic, reducing loss of sample from adhesion to tip interior.

Details

- Inner surface is extremely hydrophobic, increasing sample reproducibility and sample accuracy
- Crystal clear tips allow you to easily see transparent samples
- Reference marks on many Molecular BioProducts tips provide a quick and convenient guide to visually establish whether you are pipetting desired volume of liquid
- All Pure tips are racked, pre-sterilized and certified free of RNase, DNase, DNA and pyrogen
- Other Molecular BioProducts racked and bulk tips are nonsterile and autoclavable

Pipetter compatibility: SoftFit L tips are compatible with LTS pipetters made by Rainin Instrument, LLC. All other tips have universal fit for research grade pipetters.



| Cat. No. | ART Tip | Volume | Packaging |
|--------------------------|---------------------------|---------------|-------------------------------|
| Pure Sterile Tips | | | |
| 3501-05 | Pure 10 | 0.1-10 µL | Racked, 10 × 96 |
| 3511-05 | Pure 10, Reach | 0.1-10 µL | Racked, 10 × 96 |
| 3721-05 | Pure 20L, SoftFit L | 0.1-20 µL | Racked, 10 × 96 |
| 3721-05-HR | Pure 20L, SoftFit L | 0.1-20 µL | Hinged Rack, 10 × 96 |
| 3931-05 | Pure 200 | 0.5-200 µL | Racked, 10 × 96 |
| 3551-05 | Pure 250 | 0.5-250 µL | Racked, 10 × 96 |
| 3751-05 | Pure 200L, SoftFit L | 0.5-250 µL | Racked, 10 × 96 |
| 3751-05-HR | Pure 200L, SoftFit L | 0.5-250 µL | Hinged Rack, 10 × 96 |
| 3731-05 | Pure 300L, SoftFit L | 10-300 µL | Racked, 10 × 96 |
| 3731-05-HR | Pure 300L, SoftFit L | 10-300 µL | Hinged Rack, 10 × 96 |
| 3771-05 | Pure 300U | 1.0-300 µL | Racked, 10 × 96, Multichannel |
| 3781-05 | Pure 1000L, SoftFit L | 100-1000 µL | Racked, 8 × 96 |
| 3781-05-HR | Pure 1000L, SoftFit L | 100-1000 µL | Hinged Rack, 8 × 96 |
| 3581-05 | Pure 1000 | 100-1000 µL | Racked 8 × 96 |
| 3741-05 | Pure 1200L, SoftFit | L 100-1200 µL | Racked, 8 × 96 |
| 3741-05-HR | Pure 1200L, SoftFit | L 100-1200 µL | Hinged Rack, 8 × 96 |
| Nonsterile Tips | | | |
| 3502-05 | 10 Micro Tip, MicroPoint | 0.1-10 µL | Racked, 10 × 96 |
| 3500-05 | 10 Micro Tip, MicroPoint | 0.1-10 µL | Bulk |
| 3512-05 | 10 Reach, MicroPoint | 0.1-10 µL | Racked, 10 × 96 |
| 3510-05 | 10 Reach, MicroPoint | 0.1-10 µL | Bulk |
| 3722-05 | 20L, SoftFit L | 0.1-20 µL | Racked, 10 × 96 |
| 3722-05-HR | 20L, SoftFit L | 0.1-20 µL | Hinged Rack, 10 × 96 |
| 3932-05 | 200 | 0.5-200 µL | Racked, 10 × 96 |
| 3930-05 | 200 | 0.5-200 µL | Bulk |
| 3550-05 | 200 | 0.5-200 µL | Bulk |
| 3552-05 | 200L, SoftFit, MicroPoint | 0.5-200 µL | Racked, 10 × 96 |
| 3752-05 | 200L, SoftFit L | 10-250 µL | Racked, 10 × 96 |
| 3752-05-HR | 200L, SoftFit L | 10-250 µL | Hinged Rack, 10 × 96 |
| 3732-05 | 300L, SoftFit L | 10-300 µL | Racked, 10 × 96 |
| 3732-05-HR | 300L, SoftFit L | 10-300 µL | Hinged Rack, 10 × 96 |
| 3772-05 | 300U | 1.0-300 µL | Racked, 10 × 96, Multichannel |
| 3782-05 | 1000L, SoftFit L | 100-1000 µL | Racked, 8 × 96 |
| 3782-05-HR | 1000L, SoftFit L | 100-1000 µL | Hinged Rack, 8 × 96 |
| 3582-05 | 1000 | 100-1000 µL | Racked 8 × 96 |
| 3742-05 | 1200L, SoftFit L | 100-1200 µL | Racked, 8 × 100 |
| 3742-05-HR | 1200L, SoftFit L | 100-1200 µL | Hinged Rack, 8 × 96 |

Molecular BioProducts Pipette Tips Compatibility Chart

| ART Cross Reference | Gilson/Pipetman | | | | | | | | RAININ | | | | | | | | | | | | | | | |
|-------------------------|-----------------|---------------------|---------------------|-----------------------|-----------------------|-------------------------|----------------------|------------------------|-----------------------|-------------------|----------------------|----------------------|-------------------------|-------------------------|---------------------------|---------------------------|-----------------------------|-----------------------------|-------------------------------|---------------------------------------|---------------------------------------|---|---|--|
| | Pipetman P-2ul | Pipetman P-10 10 µL | Pipetman P-20 20 µL | Pipetman P-100 100 µL | Pipetman P-200 200 µL | Pipetman P-1000 1000 µL | Pipetman P-5000 5 mL | Pipetman P-10 mL 10 mL | EP-10/EZ-MIC-10 10 µL | EP-25/EZ-25 25 µL | EP-100/EZ-100 100 µL | EP-250/EZ-100 250 µL | EP-1000/EZ-1000 1000 µL | Pipet Lite LTS L-2 2 µL | Pipet Lite LTS L-10 10 µL | Pipet Lite LTS L-20 20 µL | Pipet Lite LTS L-100 100 µL | Pipet Lite LTS L-200 200 µL | Pipet Lite LTS L-1000 1000 µL | LTS 8 & 12 Multi-Channel L8-10 L12-10 | LTS 8 & 12 Multi-Channel L8-20 L12-20 | LTS 8 & 12 Multi-Channel L8-200 L12-200 | LTS 8 & 12 Multi-Channel L8-300 L12-300 | |
| 10 µL | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 10 | • | • | | | | | | | • | | | | | | | | | | | | | | | |
| ART 10 REACH | • | • | | | | | | | • | | | | | | | | | | | | | | | |
| ART 10F | | | | • | • | | | | | | | | | | | | | | | | | | | |
| ART 20E | | | | | | | | | | | | | | | | | | | | | | | | |
| 20 µL to 100 µL | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 20 | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 20P | | | | | | | | | • | | | | | | | | | | | | | | | |
| ART GEL 20P | | | • | • | • | | | | | | | | | | | | | | | | | | | |
| ART 20L | | | | | | | | | | | | | | | • | | | | | | | | | |
| ART 20 ERGO-F | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 50U | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 100 | | | | • | • | | | | | | | | | | | | | | | | | | | |
| ART 100E | | | | • | • | | | | | | | | | | | | | | | | | | | |
| ART 100 ERGO-G | | | | • | | | | | | | | | | | | | | | | | | | | |
| ART 100 ERGO-F | | | | | | | | | | | | | | | | | | | | | | | | |
| ART GEL 100 | | | | | | | | | | | | | | | | | | | | | | | | |
| ART XLP | | | | | | | | | | | | | | | | | | | | | | | | |
| 200 µL to 500 µL | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 200 | | | | • | • | | | | | | | • | | | | | | | | | | | | |
| ART 200U | | | | | • | | | | | | | • | | | | | | | | | | | | |
| ART 200G | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 200L | | | | | | | | | | | | | | | | | • | • | | | | | | |
| ART 200 ERGO-G | | | | | • | | | | | | | | | | | | | | | | | | | |
| ART 200 ERGO-F | | | | | | | | | | | | | | | | | | | | | | | | |
| ART XLP 200 | | | | | | | | | | | | | | | | | | | | | | | | |
| ART XLG | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 300 | | | | | | | | | | | | | | | | | | | | | | | | |
| 1 mL to 5 mL | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 1000 | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 1000E | | | | | | • | | | | | | • | | | | | | | | | | | | |
| ART 1000 REACH | | | | | | • | | | | | | | | | | | | | | | | | | |
| ART 1000G | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 1000L | | | | | | | | | | | | | | | | | | | | | • | | | |
| ART 1000 ERGO-G | | | | | | • | | | | | | | | | | | | | | | | | | |
| ART 1000 ERGO-F | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 1200 | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 1250 | | | | | | | | | | | | | | | | | | | | | | | | |
| ART 5000 | | | | | | | | | | | | | | | | | | | | | | | | |

Molecular BioProducts Pipette Tips Compatibility Chart

| ART Cross Reference | mLA | | | | BRANDTECH | | | | | | | | | | | | | | | | |
|-------------------------|-----------|------------|-------------|----------------------|-----------------------------------|-----------------------------------|-----------------------------------|-----------------------------------|------------------------------------|------------------------------------|------------------------------------|---------------------|---------------------------------|---------------------|---------------------|----------------------|----------------------|----------------------|-----------------------|-----------------------|---|
| | mLA 50 µL | mLA 200 µL | mLA 5000 µL | Digital Pipet 200 µL | Transferpette Multi-Channel 10 µL | Transferpette Multi-Channel 20 µL | Transferpette Multi-Channel 25 µL | Transferpette Multi-Channel 50 µL | Transferpette Multi-Channel 100 µL | Transferpette Multi-Channel 200 µL | Transferpette Multi-Channel 300 µL | Transferpette 10 µL | Transferpette 20 µL Ultra Micro | Transferpette 20 µL | Transferpette 50 µL | Transferpette 100 µL | Transferpette 200 µL | Transferpette 250 µL | Transferpette 1000 µL | Transferpette 5000 µL | |
| 10 µL | | | | | | | | | | | | | | | | | | | | | |
| ART 10 | | | | | • | • | | | | | | | • | • | | | | | | | |
| ART 10 REACH | | | | | • | • | | | | | | | • | • | | | | | | | |
| ART 10F | • | • | | | | | • | • | • | | | | | | | | | | | | |
| ART 20E | | | | | • | • | | | | | | | • | • | | | | | | | |
| 20 µL to 100 µL | | | | | | | | | | | | | | | | | | | | | |
| ART 20 | | | | | • | • | | | | | | | | | | | | | | | |
| ART 20P | • | | | | | | • | • | • | • | • | | | • | • | • | • | | | | |
| ART GEL 20P | | | | | | | | • | • | • | • | | | • | • | • | • | | | | |
| ART 20L | | | | | | | | | | | | | | | | | | | | | |
| ART 20 ERGO-F | | | | | | | | | | | | | | | | | | | | | |
| ART 50U | | | | | | | | | | | | | | | | | | | | | |
| ART 100 | • | • | | | | | • | • | • | | | | | | | | • | • | | | |
| ART 100E | • | | | | | | • | | • | • | • | | | • | • | | • | | | | |
| ART 100 ERGO-G | | | | | | | | | | | | | | | | | | | | | |
| ART 100 ERGO-F | | | | | | | | | | | | | | | | | | | | | |
| ART GEL 100 | | | | | | | | | • | • | • | | | | | | | | • | | |
| ART XLP | | | | | | | • | | • | • | • | | | | • | • | • | | | | |
| 200 µL to 500 µL | | | | | | | | | | | | | | | | | | | | | |
| ART 200 | • | • | | • | | | • | • | • | | | | | | | | | | • | | |
| ART 200U | | | | | | | | | | • | | | | | | | | | • | | |
| ART 200G | | | | | | | | | | | | | | | | | | | | | |
| ART 200L | | | | | | | | | | | | | | | | | | | | | |
| ART 200 ERGO-G | | | | | | | | | | | | | | | | | | | | | |
| ART 200 ERGO-F | | | | | | | | | | | | | | | | | | | | | |
| ART XLP 200 | | | | | | | • | | | • | • | | | | | | | | • | | |
| ART XLG | • | • | | | | | • | | | • | • | | | | | | | | • | | |
| ART 300 | | | | | | | • | | • | • | • | | | | | | | | • | | |
| 1 mL to 5 mL | | | | | | | | | | | | | | | | | | | | | |
| ART 1000 | | | | | | | | | | | | | | | | | | | | | • |
| ART 1000E | | | | | | | | | | | | | | | | | | | | | • |
| ART 1000 REACH | | | | | | | | | | | | | | | | | | | | | • |
| ART 1000G | | | | | | | | | | | | | | | | | | | | | • |
| ART 1000L | | | | | | | | | | | | | | | | | | | | | |
| ART 1000 ERGO-G | | | | | | | | | | | | | | | | | | | | | |
| ART 1000 ERGO-F | | | | | | | | | | | | | | | | | | | | | |
| ART 1200 | | | | | | | | | | | | | | | | | | | | | |
| ART 1250 | | | | | | | | | | | | | | | | | | | | | |
| ART 5000 | | | • | | | | | | | | | | | | | | | | | | |

Molecular BioProducts Pipette Tips Compatibility Chart

| CORNING | | | | | | HAMILTON | | | | | | CAPP | | | | Thermo Scientific | | | | | | | | | | | |
|--------------|------|------|-------|------|--------|----------|-------|------|--------|------|---------|----------------|-----------------|-----------------|------------------|-------------------|-------------------|-----------------|-----------------|------------------|-------------------|-----------------------|------------------------|--------------------------|--------------------------|--|--|
| 4959 | 2 µL | 4961 | 20 µL | 4963 | 200 µL | 4960 | 10 µL | 4962 | 100 µL | 4964 | 1000 µL | Precision 2 µL | Precision 10 µL | Precision 25 µL | Precision 100 µL | Precision 300 µL | Precision 1000 µL | Capp Aero 10 µL | Capp Aero 50 µL | Capp Aero 200 µL | Capp Aero 1000 µL | Matrix Impact 1250 µL | Matrix Impact2 1250 µL | Matrix Electronic 125 µL | Matrix Electronic 250 µL | | |
| 10 µL | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | | | | | | | | | | | | | | | | | | | | | | | | |
| | | | | </ | | | | | | | | | | | | | | | | | | | | | | | |

Molecular BioProducts Pipette Tips Compatibility Chart

| ART Cross Reference | | ART Filter (sterile) | | | | | | | | | | | | | | | | | | | | | |
|--|--------------------|----------------------|-------|----------|-------|-------|-----------|---------|------------------|------|-------|-------|------------|---------|------|------|------|------|-------|-------------|-------------|------------|---|
| | | 10 | 10F | 10 Reach | 20E | 20P | Ergo-F 20 | GEL 20P | Solvent Safe 20P | 100 | 100E | 100E | Ergo-F 100 | GEL 100 | 101 | 200 | 300 | 1000 | 1000E | Ergo-F 1000 | Ergo-G 1000 | 1000 Reach | |
| Thermo Scientific Pipetter Description | Cat. No. | 2139 | 2139F | 2140 | 2149E | 2149P | 2350 | 2155P | 5449P | 2065 | 2065E | 5479E | 2360 | 2155 | 5469 | 2069 | 2070 | 2279 | 2079E | 2380 | 2340 | 2079 | |
| F1 0.2-2 µL | 4641010 | • | | • | | | | | | | | | | | | | | | | | | | |
| F1 0.5-5 µL | 4641020 | • | | • | | | | | | | | | | | | | | | | | | | |
| F1 1-10 µL micro | 4641030 | • | | • | | | | | | | | | | | | | | | | | | | |
| F1 1-10 µL | 4641040 | | • | | | | | | | | | | | | | | | | | | | | |
| F1 2-20 µL micro | 4641050 | | | | | | | | | | | | | | | | | | | | | | |
| F1 2-20 µL | 4641060 | | | | | • | • | • | • | | | | | | | | | | | | | | |
| F1 10-100 µL | 4641070 | | | | | | | | | • | | | • | • | • | • | • | | | | | | |
| F1 20-200 µL | 4641080 | | | | | | | | | | | | | | | • | • | | | | | | |
| F1 30-300 µL | 4641090 | | | | | | | | | | | | | | | | • | | | | | | |
| F1 100-1000 µL | 4641100 | | | | | | | | | | | | | | | | | • | | • | • | • | |
| F1 0.5-5 µL | 4641110 | | | | | | | | | | | | | | | | | | | | | | |
| F1 1-10 µL | 4641120 | | | | | | | | | | | | | | | | | | | | | | |
| F1 8-ch 1-10 µL | 4661000 | • | | | | | | | | | | | | | | | | | | | | | |
| F1 12-ch 1-10 µL | 4661040 | • | | | | | | | | | | | | | | | | | | | | | |
| F1 16-ch 1-10 µL | 4661080 | | | | | | | | | | | | | | | | | | | | | | |
| F1 8-ch 5-50 µL | 4661010 | | | | | | | | | | | | • | | | | | | | | | | |
| F1 12-ch 5-50 µL | 4661050 | | | | | | | | | | | | • | | | | | | | | | | |
| F1 16-ch 5-50 µL | 4661090 | | | | | | | | | | | | | | | | | | | | | | |
| F1 8-ch 10-100 µL | 4661020 | | | | | | | | | | | | • | | | | | | | | | | |
| F1 12-ch 10-100 µL | 4661060 | | | | | | | | | | | | • | | | | | | | | | | |
| F1 8-ch 30-300 µL | 4661030 | | | | | | | | | | | | | | | | • | | | | | | |
| F1 12-ch 30-300 µL | 4661070 | | | | | | | | | | | | | | | | • | | | | | | |
| F2 0.2-2 µL | 4642010 | • | | • | | | | | | | | | | | | | | | | | | | |
| F2 0.5-5 µL | 4642020 | • | | • | | | | | | | | | | | | | | | | | | | |
| F2/F3 1-10 µL micro | 4642030 4640000 | • | | • | | | | | | | | | | | | | | | | | | | |
| F2/F3 1-10 µL | 4642040 4640010 | | • | | | • | • | | • | | | | | • | | | | | | | | | |
| F2/F3 2-20 µL micro | 4642050 4640020 | | | | | | | | | | | | | | | | | | | | | | |
| F2/F3 2-20 µL | 4642060 4640030 | | | | | • | • | | • | | | | | | | | | | | | | | |
| F2/F3 10-100 µL | 4642070 4640040 | | | | | | | | | • | | | • | • | • | • | • | | | | | | |
| F2/F3 20-200 µL | 4642080 4640050 | | | | | | | | | | | | | | | | • | | | | | | |
| F2/F3 100-1000 µL | 4642090 4640060 | | | | | | | | | | | | | | | | | | | | | | • |
| F2 0.5-5 mL | 4642100 | | | | | | | | | | | | | | | | | | | | | | |

Molecular BioProducts Pipette Tips Compatibility Chart

| ART Cross Reference | | ART Filter (sterile) | | | | | | | | | | | | | | | | | | | | | |
|--|--------------------|----------------------|-------|----------|-------|-------|-----------|---------|------------------|------|-------|-------|------------|---------|------|------|------|------|-------|-------------|-------------|------------|---|
| | | 10 | 10F | 10 Reach | 20E | 20P | Ergo-F 20 | GEL 20P | Solvent Safe 20P | 100 | 100E | 100E | Ergo-F 100 | GEL 100 | 101 | 200 | 300 | 1000 | 1000E | Ergo-F 1000 | Ergo-G 1000 | 1000 Reach | |
| Thermo Scientific Pipetter Description | Cat. No. | 2139 | 2139F | 2140 | 2149E | 2149P | 2350 | 2155P | 5449P | 2065 | 2065E | 5479E | 2360 | 2155 | 5469 | 2069 | 2070 | 2279 | 2079E | 2380 | 2340 | 2079 | |
| F2/F3 1-10 µL | 4642110 4640070 | | | | | | | | | | | | | | | | | | | | | | |
| F2 8-ch 1-10 µL | 4662000 | • | | | | | | | | | | | | | | | | | | | | | |
| F2 12-ch 1-10 µL | 4662040 | • | | | | | | | | | | | | | | | | | | | | | |
| F2 16-ch 1-10 µL | 4662080 | | | | | | | | | | | | | | | | | | | | | | |
| F2 8-ch 5-50 µL | 4662010 | | | | | | | | | | | | • | | | | | | | | | | |
| F2 12-ch 5-50 µL | 4662050 | | | | | | | | | | | | • | | | | | | | | | | |
| F2 16-ch 5-50 µL | 4662090 | | | | | | | | | | | | | | | | | | | | | | |
| F2 8-ch 10-100 µL | 4662020 | | | | | | | | | | | | • | | | | | | | | | | |
| F2 12-ch 10-100 µL | 4662060 | | | | | | | | | | | | • | | | | | | | | | | |
| F2 8-ch 30-300 µL | 4662030 | | | | | | | | | | | | | | | | • | | | | | | |
| F2 12-ch 30-300 µL | 4662070 | | | | | | | | | | | | | | | | • | | | | | | |
| Novus 1-10 µL micro | 46200000 | | | | • | | | | | | | | | | | | | | | | | | |
| Novus 1-10 µL | 46200100 | | | | | | | | | | | | | | | | | | | | | | |
| Novus 5-50 µL micro | 46200200 | | | | | | | | | | | | | | | | | | | | | | |
| Novus 5-50 µL | 46200300 | | | | | | | | | | | | | | | | • | | | | | | |
| Novus 10-100 µL | 46200400 | | | | | | | | | | | | | | | • | • | | | | | | |
| Novus 30-300 µL | 46200500 | | | | | | | | | | | | | | | | • | | | | | | |
| Novus 100-1000 µL | 46200600 | | | | | | | | | | | | | | | | | • | | | | | • |
| Novus 0.5-5 mL | 46200700 | | | | | | | | | | | | | | | | | | | | | | |
| Novus 1-10 mL | 46200800 | | | | | | | | | | | | | | | | | | | | | | |
| Novus 8-ch 1-10 µL | 46300000 | • | | • | | | | | | | | | | | | | | | | | | | |
| Novus 12-ch 1-10 µL | 46300100 | • | | • | | | | | | | | | | | | | | | | | | | |
| Novus 8-ch 5-50 µL | 46300200 | | | | | | | | | | | | | | | | | | | | | | |
| Novus 12-ch 5-50 µL | 46300300 | | | | | | | | | | | | | | | | | | | | | | |
| Novus 16-ch 5-50 µL | 46300700 | | | | | | | | | | | | | | | | | | | | | | |
| Novus 8-ch 30-300 µL | 46300400 | | | | | | | | | | | | | | | | • | | | | | | |
| Novus 12-ch 30-300 µL | 46300500 | | | | | | | | | | | | | | | | • | | | | | | |
| Novus 8-ch 100-1200 µL | 46300800 | | | | | | | | | | | | | | | | | | | | | | |




Automated Liquid Handling

Delivering productivity with speed and precision

Thermo Scientific automated liquid handling systems are ideal for medium or high throughput applications, meeting multiple liquid handling requirements that include channels, microplate type or liquid volumes. Our automated liquid handling instruments reduce the time spent on repetitive pipetting tasks at any throughput level, delivering:

- Consistent results with increased efficiency for pipetting applications
- Improved productivity with higher liquid and microplate handling throughput
- Scalable and flexible options – including a broad array of accessories – to address a wide range of liquid handling requirements

Automated Liquid Handling



Thermo Scientific Multidrop Reagent Dispensers and Accessories

| | |
|--|----|
| Thermo Scientific Multidrop Combi nL Reagent Dispenser | 76 |
| Thermo Scientific Multidrop Combi Reagent Dispenser | 77 |
| Thermo Scientific Multidrop 384 Reagent Dispenser | 78 |
| Thermo Scientific Multidrop DW Reagent Dispenser | 79 |
| Thermo Scientific Multidrop Combi, 384 and DW Dispensing Cassettes | 80 |
| Standard tube dispensing cassette | 80 |
| Small Tube Dispensing Cassette | 80 |
| SMART dispensing cassettes | 80 |
| Thermo Scientific FILLit Software for Multidrop Reagent Dispensers | 82 |
| Multidrop Combi nL FILLit Software | 82 |

Thermo Scientific Matrix Dispensers and Accessories

| | |
|---|----|
| Thermo Scientific Matrix WellMate Microplate Dispenser | 83 |
| Thermo Scientific Matrix WellMate Disposable Tubing Cartridges | 84 |
| Thermo Scientific Matrix WellMate Stacker Base Unit | 85 |
| Thermo Scientific Matrix PlateMate Plus/WellMate Stacker Chimneys | 86 |
| Thermo Scientific Matrix Hydra II Liquid Handling System | 87 |
| Thermo Scientific Matrix Hydra DT Pipetting Workstation | 88 |

Thermo Scientific Versette Pipetting Workstations and Accessories

| | |
|--|----|
| Thermo Scientific Versette Pipetting Workstation | 89 |
| Thermo Scientific Versette Pipetting Heads | 90 |
| Thermo Scientific Versette Accessories | 91 |

Automation Tips and Accessories

| | |
|--|-----|
| Thermo Scientific Versette ClipTip Automation Tips | 92 |
| Thermo Scientific Matrix D.A.R.T.s Tips | 93 |
| Thermo Scientific Matrix Filtered D.A.R.T.s Tips | 94 |
| Thermo Scientific PocketTip D.A.R.T.s Automation Tips | 95 |
| Thermo Scientific Matrix D.A.R.T.s Tip Transfer Tool | 96 |
| Thermo Scientific Matrix Disposable Automation Reservoirs | 97 |
| Molecular BioProducts BioRobotix ART Filter Tips for Automated Workstations | 98 |
| Molecular BioProducts BioRobotix Pipet Tips, Standard, Black | 100 |
| Molecular BioProducts BioRobotix Pipette Tips: Standard, Natural, Nonsterile | 101 |

Microplate Stackers and Accessories

| | |
|--|-----|
| Thermo Scientific RapidStak Automated Microplate Stacker with Polara RS Software | 103 |
| Thermo Scientific RapidStak Accessories: Microplate Stacks | 104 |
| Thermo Scientific RapidStak Instrument Drivers | 104 |
| Thermo Scientific Polara RS | 105 |
| Thermo Scientific RapidStak DLL Programming Kits | 105 |

Microplate Movers


| | |
|---|-----|
| Thermo Scientific Orbitor RS Microplate Mover | 107 |
| Thermo Scientific CatalySt Express Microplate Handler | 108 |

Plastic Reservoirs • Storage Blocks

| | |
|---|-----|
| Thermo Scientific Nalgene Disposable Robotic Reservoirs | 109 |
| Thermo Scientific Nunc Disposable Plastic Reservoirs | 110 |
| Thermo Scientific Matrix Deepwell Storage Blocks | 111 |

Thermo Scientific Multidrop Combi nL Reagent Dispenser




 Increase throughput with a RapidStak Automated Microplate Stacker

The Multidrop Combi* nL bulk reagent dispenser for nanoliter-to-microliter volumes provides accurate, consistent, high throughput dispensing for laboratories across a volume range of 50 nL to 50 μ L.

The Multidrop Combi nL is a nano- to microvolume bulk reagent dispenser. It offers easy entry to precise, accurate and reliable low-volume dispensing for pharmaceutical and biotechnology laboratories.

Details

The Multidrop Combi nL reagent dispenser offers low volume dispensing for volumes 50 nL to 50 μ L into 96-, 384- and 1536- well plates with excellent precision and accuracy.

- Reliable low volume dispensing: For precise dispensing over a 50 nL to 50 μ L range, providing high quality assay data
- Increased flexibility: Dispenses repeatably into plates with variable height, allowing the user to fill rows, columns, or wells with the easy to use FILLit Software
- Effortless dispensing: The visual, icon-based onboard user interface makes all functions easy to setup, use, and maintain
- Reduced reagent costs: For low volume dispensing, a minimal dead volume and the backflush function minimize the use of expensive reagents
- High throughput: Fast dispensing combined with full robotic compatibility ensures increased assay throughput for laboratories requiring low volume assay formats

Warranty and Service Offering: One year

Recommended for: Assay development, primary and secondary screening, genomics and proteomics research, PCR setup, sequencing setup, cell based assays and bead based assays



| Specifications | |
|------------------------------------|---|
| Plate Types | 96-, 384- and -1536-well plates |
| Dispensing Volume Range | 50 nL to 50 μ L |
| Dispensing Precision | 50 nL: CV \leq 10% 0.5 μ L: CV \leq 5% 1-10 μ L: CV \leq 4% >10 μ L: CV \leq 2% |
| Dispensing Accuracy | <1 μ L: \pm 5% >1 μ L: \pm 2% |
| Dispensing Speed | 384-well plate: 50nL in 6 seconds, 1 μ L in 8 seconds 1536-well plate: 50nL in 21 seconds, 1 μ L in 27 seconds |
| Dispensing Increments | 1 nL increments 50-999 nL 10 nL increments 1.00-9.99 μ L 100 nL increments 10.0-50.0 μ L FILLit Software 1 nL increments |
| Dead Volume | <1.2 mL [†] |
| Interface | Serial RS-232, USB |
| Dimensions (W x D x H) | 14 x 14.8 x 8.6 in. (35.5 x 37.5 x 22 cm) |
| Weight | 9.6 kg (21.2 lbs.) |
| <i>† Reagent recovery possible</i> | |

| Cat. No. | Description |
|----------|---|
| 5840400 | Volume range 50 nL to 50 μ L; accommodates 96-, 384- and 1536-well microplates. |

Thermo Scientific Multidrop Combi Reagent Dispenser

The Multidrop Combi Bulk Reagent Dispenser offers unrivaled levels of flexibility and performance to meet the requirements of reagent dispensing in pharmaceutical and biotechnology laboratories.

The Multidrop Combi Reagent Dispenser is easy to use and combines the most versatile features for reagent dispensing with excellent performance. It's ideal for drug discovery, genomics and proteomics assays.

The Multidrop Combi bulk reagent dispenser uses an autoclavable 8-channel dispensing cassette. It supports a wide selection of plates and volume ranges, offering fast dispensing and high-throughput operation.



Details

- Provides precise dispensing over a 0.5 to 2500 μL range, ensuring reproducible assay data
- Accommodates microplates from 6 to 1536 wells and plate heights of 5 to 50 mm
- Visual icon-based graphic display makes it easy to use and program without formal training
- Minimal dead volume and back-flushing reduce reagent costs
- Easy to use FILLit Software provides increased flexibility and functionality
- Full robotic compatibility provides increased throughput
- Uses 8-channel detachable and autoclavable dispensing cassettes that are standard across the Multidrop range

Multidrop Combi SMART

- Adds an advanced built-in tracing system to the Multidrop Combi
- Provides improved reliability and cassette lifetime traceability, greatly enhancing the user's efficiency and reporting capability

Recommended for: Assay development, primary and secondary screening, compound storage, genomics and proteomics research, PCR setup, sequencing setup, cell based assays, bead based assays, and ELISA assays

Includes: One standard tube dispensing cassette and two small tube dispensing cassettes (Multidrop Combi cat. nos. 24072670, 24073290, 24073295 and Multidrop Combi SMART cat.nos. 24072675, 24073292 and 24073297)

Warranty: One year

| Specifications | |
|-------------------------|---|
| Plate Type | 6- to 1536-well plates |
| Dispensing Volume Range | 0.5 to 2500 μL |
| Dispensing Precision | Small tube dispensing cassette: 0.5 μL , CV \leq 10%; 2 μL , CV \leq 5%; 10 μL , CV \leq 3%; >10 μL , CV \leq 3% Standard tube dispensing cassette: 5 μL , CV \leq 10%; 20 μL , CV \leq 1.5%; 100 μL , CV \leq 1%; >100 μL , CV \leq 1% |
| Dispensing Accuracy | Small tube dispensing cassette: 2 μL , \pm 10% (typical); 10 μL , \pm 5% (typical); >10 μL , \pm 5% (typical) Standard tube dispensing cassette: 5 μL , \pm 3% (typical); 20 μL , \pm 2% (typical); 100 μL , \pm 1% (typical); >100 μL , \pm 1% (typical) |
| Dispensing Speed | 96-well plate: 10 μL in 3 sec.; 20 μL in 4 sec.; 100 μL in 10 sec. 384-well plate: 1 μL in 5 sec.; 5 μL in 5 sec.; 10 μL in 6 sec.; 20 μL in 9 sec. 1536-well plate: 1 μL in 14 sec.; 5 μL in 26 sec. |
| Dispensing Increments | Small tube dispensing cassette: 0.5 μL increments 0.5-50 μL Standard tube dispensing cassette: 5 μL increments 5-2500 μL |
| Dead Volume | <1 mL volume range 0.5-50 μL ; <7 mL volume range >50 μL [†] |
| Interface | Serial RS-232, USB |
| Dimensions (W x D x H) | 14 x 12.9 x 8.6 in. (35.5 x 33 x 22 cm) |
| Weight | 9.1 kg (20.1 lbs.) |

[†] Reagent recovery option

| Cat. No. | Description |
|----------|--------------------------------------|
| 5840300 | 100-240 V 50/60 Hz |
| 5840310 | With SMART option 100-240 V 50/60 Hz |

Increase throughput with a RapidStak Automated Microplate Stacker

Thermo Scientific Multidrop 384 Reagent Dispenser



The Multidrop 384 Bulk Reagent Dispenser provides reliable microvolume reagent dispensing for pharmaceutical and biotechnology laboratories.

The Multidrop 384 reagent dispenser provides continuous high-speed dispensing of various liquids with excellent precision and robot compatibility for higher productivity.

Details

Multidrop 384 instrument dispenses a wide range of solutions into 96 and 384 well plates across 5 to 395 μL volume range.

- Precise dispensing with reproducible dispensing results
- High-speed dispensing for high assay throughput
- Easy setup with quick plate and volume selection
- Easy to integrate with different robotic setups for higher throughput

Convenient, autoclavable dispensing cassette

- 8-channel detachable, autoclavable dispensing cassette to ensure sterile conditions and avoid cross contamination
- All reagent lines can be backflushed into the reagent bottle, minimizing the loss of expensive reagents

Includes: One standard tube dispensing cassette (Cat. No. 24072670)

Warranty and Service Offering: One year

Recommended for: Assay development, primary and secondary screening, compound storage, genomics and proteomics research, cell based assays and ELISA assays

▶ Increase throughput with a RapidStak Automated Microplate Stacker

| Specifications | |
|--------------------------------------|---|
| Plate Types | 96-, 384-well plates |
| Dispensing Volume Range | 5 to 395 μL |
| Dispensing Precision | 5 μL : CV \leq 10% (typical) 20 μL : CV \leq 1.5% (typical) 100 μL : CV \leq 1% (typical) |
| Dispensing Accuracy | 5 μL : \pm 3% (typical) 20 μL : \pm 2% (typical) 100 μL : \pm 1% (typical) |
| Dispensing Speed | 96-well plate: 20 μL in 5 seconds, 50 μL in 8 seconds 384-well plate: 20 μL in 20 seconds, 50 μL in 25 seconds |
| Dispensing Increments | 5 μL increments |
| Dead Volume | <7 mL [†] |
| Interface | RS-232 |
| Dimensions (W x D x H) | 12.2 x 12.6 x 6.1 in. (31 x 32 x 15.5 cm) |
| Weight | 6.2 kg (13.7 lbs.) |
| [†] Reagent recovery option | |

| Cat. No. | Description |
|----------|--------------------|
| 5840150 | 220-240 V 50/60 Hz |
| 5840157 | 100-120V 50/60Hz |

Thermo Scientific Multidrop DW Reagent Dispenser

The Multidrop DW Reagent Dispenser is a high-speed reagent dispenser designed for repetitive dispensing of large volumes in pharmaceutical and biotechnology laboratories.

The Multidrop DW Reagent Dispenser releases various liquids into 96 standard and DW plates across the 20 to 995 μL volume range.

This dispenser accelerates larger scale assays by providing superior precision and accuracy for repetitive dispensing of large volumes.

Details

- Precise dispensing ensures high-quality assay results and reproducible data
- Wide volume range provides flexibility and increased speed in assays requiring larger volumes
- With quick plate and volume selection combined with easy maintenance, the Multidrop DW makes assay setup quick and easy


Convenient, autoclavable dispensing cassette

- Uses an 8-channel detachable and autoclavable dispensing cassette to deliver sterile conditions and to avoid cross contamination
- All reagent lines can be backflushed into the reagent bottle, minimizing the loss of expensive reagents

Recommended for: Cell based assays, compound storage, genomics and proteomics research, ELISA assays and microbiological tests

Includes: One standard tube dispensing cassette (Cat. No. 24072670)

Warranty: One year



| Specifications | |
|-------------------------|--|
| Plate Types | 96-well plates; 96-deep well plates; 1.1 mL tubes in 96 tube rack; Micronic 96-place tube racks |
| Dispensing Volume Range | 20 to 995 μL |
| Dispensing Precision | 20 μL : CV \leq 1.5% (typical) 100 μL : CV \leq 1% (typical) 900 μL : CV \leq 0.5% (typical) |
| Dispensing Accuracy | 20 μL : \pm 2% (typical); 100 μL : \pm 1% (typical); 900 μL : \pm 1% (typical) |
| Dispensing Speed | 96-well plate: 20 μL in 5 seconds; 50 μL in 8 seconds; 300 μL in 26 seconds; 900 μL in 74 seconds |
| Dispensing Increments | 5 μL increments |
| Dead Volume | <7 mL [†] |
| Interface | RS-232 |
| Dimensions (W x D x H) | 12.2 x 12.6 x 7.1 in. (31 x 32 x 18 cm) |
| Weight | 6.2 kg (13.7 lbs.) |

[†] Reagent recovery option

| Cat. No. | Description |
|----------|-------------------|
| 5840177 | 100-120V 50/60 Hz |
| 5840170 | 200-240V 50/60 Hz |



Our comprehensive line of Multidrop Dispensing Cassettes offers optimal solutions for a wide volume range and a variety of reagents for the best performance and results.

The autoclavable Multidrop Dispensing Cassettes come with different tubing sizes, as well tip material and size. Choose between standard and small tube cassettes, or SMART versions that store dispensing lifetime data. Dispensing cassettes are available as single versions, 5-packs and 10-packs.

Standard tube dispensing cassette

Standard tube dispensing cassettes are designed to be used with Multidrop Combi, Multidrop 384 and Multidrop DW dispensers.

Details:

- Plastic cover structure, wide silicone tubing and polypropylene tip comb with 0.5 mm tip inner diameter
- Designed for volumes 5 μL and greater, and when used with Multidrop Combi Reagent Dispensers up to 2.5 mL
- Most suitable for dispensing volumes of 20 μL and greater
- Wide tubing allows high speed for dispensing of large volumes with the standard cassettes
- Standard tube cassettes are also available with longer tubing versions, ideal for automated systems where reagent vessels are remote from the instruments

Small Tube Dispensing Cassette

Small tube dispensing cassettes can be used with Multidrop Combi and Multidrop Micro reagent dispensers.

Details:

- Plastic cover structure, silicone tubing with small tubing diameter and polypropylene or metal tip with orifice 0.22 mm
- Tip structure and design optimized for efficient and reproducible dispensing of small volumes from 0.5 μL up to 50 μL
- Small tube dispensing cassettes are available with either plastic or metal tips to accommodate a variety of reagents
- Long small tube dispensing cassettes are ideally suited for use with automated systems. The eight channels are joined as a single tube, minimizing liquid resistance.

SMART dispensing cassettes

SMART dispensing cassettes are used with the Multidrop Combi SMART unit. These cassettes are equipped with a SMART information chip to store lifetime dispensing data for the cassette. When used with the Multidrop Combi SMART cassette version, the remaining usage time is displayed. SMART cassettes are available in standard and small tube versions, and are distinguished from normal cassette versions by an antenna tag.

| Cat. No. | Description |
|----------|--|
| 24072670 | Standard tube dispensing cassette |
| 24072671 | Standard tube dispensing cassette 5-pack |
| 24072672 | Standard tube dispensing cassette 10-pack |
| 24072677 | Long standard tube dispensing cassette (length of tubing 50-200 cm in 10 cm increments) |
| 24070300 | Dispensing cassette, aluminium covers (with 40 cm tubing set) |
| 24070307 | Special dispensing cassette, aluminium covers (length of tubing 50-200 cm in 10 cm increments) |
| 24073290 | Small tube plastic tips dispensing cassette |
| 24073291 | Small tube plastic tips dispensing cassette 5-pack |
| 24073295 | Small tube metal tip dispensing cassette |
| 24073296 | Small tube metal tip dispensing cassette 5-pack |
| 24073293 | Long small tube plastic tip dispensing cassette |
| 24073298 | Long small tube metal tip dispensing cassette |
| 24072675 | SMART Standard tube dispensing cassette (with 40 cm tubing set) |
| 24072676 | SMART Standard tube dispensing cassette 5-pack |
| 24072678 | SMART Long standard tube dispensing cassette |
| 24073292 | SMART Small tube plastic tip dispensing cassette |
| 24073001 | SMART Small tube plastic tip dispensing cassette 5-pack |
| 24073297 | SMART Small tube metal tip dispensing cassette |
| 24073002 | SMART Small tube metal tip dispensing cassette 5-pack |



Thermo Scientific FILLit Software for Multidrop Reagent Dispensers



FILLit* Software is an easy-to-use tool for controlling Multidrop reagent dispensers.

FILLit Software provides advanced control and increased flexibility for the use of Multidrop reagent dispensers. The protocols can be built and deployed quickly for execution either directly from a PC or transferred into the instrument's memory.

Details

- Extensive plate list containing plates from a variety of microplate manufacturers
- Select and fill discrete wells (Multidrop Combi nL), rows (Multidrop Combi nL) or columns
- User-created protocols can be stored in the database for ongoing use or modification
- Protocols can also be downloaded to the instrument for stand alone use, or transferred from the Multidrop to the PC for further modifications

Multidrop Combi nL FILLit Software

- Pre-calibrated liquid types stored in the database allow quick selection at startup
- Adding new liquids to the database is easy with the calibration wizard
- Individual valve control allows dispensing into discrete wells

| PC Requirements | |
|------------------|--|
| Operating System | Microsoft Windows* 2000 with Service Pack 4, Windows XP Professional with Service Pack 2 or Windows Vista 32-Bit Edition |
| RAM | 1 GB |
| Drive | CD-ROM drive (or CD compatible DVD-ROM drive) |
| Monitor | XVGA monitor with 1024 × 768 resolution |
| Port | One serial or USB port available |
| Free Disk Space | 0.5 GB |
| Browser | Microsoft Internet Explorer 6.0 (or greater) installed |

| Cat. No. | For Use with |
|----------|------------------------------------|
| 5188020 | Fillit Software Multidrop Combi nL |
| 5188010 | Fillit Software Multidrop Combi |

Thermo Scientific Matrix WellMate Microplate Dispenser

The WellMate dispenser brings superior flexibility and cost-efficiency to microplate dispensing.

With a volume range of 1 to 2000 μL , the WellMate dispenser accommodates a wide selection of labware and improves cost efficiency with disposable, pre-sterilized and pre-calibrated tubing cartridges. Programming flexibility allows any user to operate the instrument, and the optional plug-and-play WellMate Stacker provides walk-away capability to meet a range of throughput requirements.

Details

- Manually adjustable dispensing head height and dispensing locations allow one instrument to accommodate virtually any 6, 12, 24, 48, 96 or 384-well microplate
- Programming options to dispense volumes from 1-2000 μL in increments of 1.0 μL
- Fluid flow can be manually reversed to recover valuable reagent contained in the lines
- The pre-sterilized, pre-calibrated tubing cartridges are economical enough to be disposable, yet rugged enough to be autoclaved for long-term use
- Allows sterile plate filling of cell culture media into 6, 12, 24, 48, 96 or 384 tissue culture plates
- The WellMate unit can be placed in a dark room for dispensing light-sensitive compounds

Recommended for: Addition of bacteria, yeast, buffers, diluents, enzymes, substrates, and ligands, labeled compounds, including fluorophores and radio labeled compounds, microbeads for assays such as LOCI or loading combichem plates, viscous solutions including scintillation fluid or glycerol/sucrose solutions

Warranty: One year

| Specifications | |
|--------------------------------------|--|
| Plate type | 6- to 384- well plates |
| Dispensing Volume Range | 1 to 2000 μL |
| Dispensing Precision | $\pm 2.0\%$ or 1.0 μL (standard bore cartridge) $\pm 2.0\%$ or 0.25 μL (small bore cartridge) $\pm 2.0\%$ or 0.25 μL (stainless steel cartridge) |
| Dispensing Accuracy | $\pm 2.0\%$ or 1.0 μL (standard bore cartridge) $\pm 2.0\%$ or 0.25 μL (small bore cartridge) $\pm 2.0\%$ or 0.25 μL (stainless steel cartridge) |
| Dispensing Speed | 96-well plate: 10 μL : 4.5 s; 50 μL : 8.7 s; 100 μL : 13.2 s 384-well plate: 5 μL : 9 s; 10 μL : 11 s; 50 μL : 26.5 s |
| Interface | Serial RS-232, OCX |
| Dimensions (W \times D \times H) | 14.5 \times 12.0 \times 8.5 in. (36.8 \times 30.5 \times 21.6 cm) |
| Weight | 11 kg (24 lbs.) |

| Cat. No. | Description |
|----------|-------------------------------|
| 20110001 | WellMate Microplate Dispenser |



Thermo Scientific Matrix WellMate Disposable Tubing Cartridges



Disposable Tubing Cartridges for the WellMate dispenser offer a convenient solution for a variety of liquid types.

Select from a wide variety of tubing cartridges to dispense everything from sensitive cellular materials to high vapor pressure fluids.

Compatible with WellMate dispensers, these pre-sterilized tubing cartridges can be quickly and easily replaced, eliminating the need for service or calibration.

Details

- Disposable, pre-sterilized, pre-calibrated tubing
- Cartridges can be quickly and easily replaced, eliminating the need for service or calibration
- Standard-bore nozzle: For volumes above 20 μL and 6 to 96-well plate filling
- Small-bore nozzle: For lower volume (1.0 μL and above) applications
- Stainless steel, PTFE-coated: For use with viscous fluids, 6 to 384-well dispensing, and low-profile plates

Includes: Five cartridges per case

| Cat. No. | Description |
|--------------------------|------------------------------|
| Tubing Assemblies | |
| 201-30001 | Standard Bore |
| 201-30002 | Small Bore |
| 201-30003 | Stainless Steel, PTFE-Coated |
| Adapter | |
| 201-40002 | Tubing Cartridge Adapter |





Thermo Scientific Matrix WellMate Stacker Base Unit

The WellMate stacker base unit can dramatically increase the throughput of microplate dispensing.

This stacker is the perfect add-on to the WellMate base unit, increasing throughput, walk-away time and the efficiency of many microplate dispensing procedures.

Details

- When combined with short or tall chimneys, enables walk-away capability for virtually any throughput level
- Accommodates low profile, standard, and deep well height plates
- All WellMate stackers are field upgradable

Requires: Thermo Scientific Matrix WellMate Base Unit

Warranty: One year

Certifications: CE

| Cat. No. | Description |
|-----------|--|
| 201-20001 | WellMate Stacker |
| 501-90750 | Annual Service Contract for WellMate Stacker |



Base Unit

Thermo Scientific Matrix PlateMate Plus/WellMate Stacker Chimneys



PlateMate Plus/WellMate Stacker chimneys provide a complete solution with the WellMate Stacker Base Unit and PlateMate Plus.

Details

Tall Chimney Stackers

- Each chimney accommodates up to 50 standard height plates
- Compatible with deep well and low profile plates

Short Chimney Stackers

- Each chimney accommodates up to 25 standard height plates
- Compatible with deep well and low profile plates
- The WellMate dispenser, stacker and short chimney combinations fit into most standard-sized biological safety cabinets

Requires: PlateMate Plus Base Unit or WellMate Stacker Base Unit

| Cat. No. | Description |
|-----------|---------------------------|
| 501-30005 | Tall (50 plate capacity) |
| 501-30006 | Short (25 plate capacity) |



Thermo Scientific Matrix Hydra II Liquid Handling System

The Hydra II syringe, fixed tip-based automated benchtop liquid handling system operates with push-button ease.

The Hydra II benchtop 96-channel liquid handler is compact and capable of automating simple pipetting tasks, including low volume and viscous sample transfers with the single press of a button. It utilizes syringe-based technology to accurately and precisely dispense even the most sensitive solutions.

Details

- Syringe-based automatic pipetting enables 100 nL - 290 µL pipetting volume range
- Choice of two 96-format volume range syringe dispense heads allows unit to be used with any 96- or 384-microplate format, including standard, low profile, deep well and standard footprint tube racks
- Precise speed control: Aspirate and dispense speeds can be optimized and independently set based on fluid type or to facilitate optimal performance
- Compact size allows use on any benchtop or in a standard enclosure
- Wash pump accommodates two wash solutions for critical wash procedures

Programming Flexibility

- Can be programmed utilizing easy-to-read, four-line LCD display, or with intuitive ControlMate software
- ControlMate software, which is used consistently for all Thermo Scientific Matrix automated liquid handling platforms, is license-free and downloadable from controlmate.net
- Available ControlMate OLE package provides developers with an ActiveX automation interface for integrated use

Syringe/Needle Assemblies for Precise Low-Volume Dispensing

- Syringe barrels are made from borosilicate glass for excellent chemical capability
- Inner chambers are precision-machined and highly polished to provide leak-free operation
- Each syringe tip features three points of sealed contact for long-lasting performance

Includes: Hydra II base unit, single microplate stage, power cord, PC communication cable, RS232 to USB connection and user manual

Required Accessories: Tip wash station and pump module

Warranty and Service Offering: One year

Certifications: CE marked



Recommended for:

- Any 96- or 384-channel low-volume positive displacement pipetting procedure
- Dispensing of solvents
- Nucleic acid preps
- PCR setup
- Plate-to-plate transfers
- Reagent addition
- Reformatting into high-density plates

| Cat. No. | Description | Accuracy |
|----------|---|------------------|
| 109611 | 96 Base Unit, 100 µL, with Duraflex Needles | ±2.5% or 0.01 µL |
| 109621 | 96 Base Unit, 290 µL, with Duraflex Needles | ±2.5% or 0.05 µL |

Thermo Scientific Matrix Hydra DT Pipetting Workstation



The Hydra* DT Pipetting Workstation is a compact, automated liquid handler that is compatible with the complete line of disposable Thermo Scientific Matrix D.A.R.T.s tips.

The Hydra DT benchtop 96-channel liquid handler is compact and capable of automating routine pipetting tasks, including plate stamping, reformatting and sample transfers with the single press of a button. The Hydra DT is compatible with our extensive range of Matrix D.A.R.T.s* (Disposable Automation Research Tips), providing excellent application flexibility with easy tip changing. Its extensive volume range accommodates multiple pipetting tasks within one laboratory.

Details

- Disposable tip-based liquid handler with a variety of disposable tip options
- Combined volume range from 0.5 μL to 300 μL
- 96 simultaneous transfers with superior accuracy and precision, with zero carryover
- Compact design requires minimal bench or hood space
- Intuitive ControlMate software makes the system easy to use
- Touch pad programming option eliminates the need for a PC
- Simple tip clamp requires only 5 lb. of force and can be actuated with one hand
- Easy touch pad operation reduces repetitive motion strain
- Compatible with Matrix D.A.R.T.s tips
- Compatible with standard biological hoods; perform biological assays with ease and efficiency using sterile or sterile-filtered Matrix D.A.R.T.s tips
- Dynamic volume range with simple tip changes; the 100 μL units are compatible with 0.5 to 30 μL tips, and 1.0 to 100 μL tips on the same instrument
- Includes: Hydra DT base unit with single position stage, ControlMate software, power cable, user manual, PC communication cable and RS323 to USB adapter

Compatible with: Thermo Scientific Matrix D.A.R.T.s tips; ControlMate Software

Warranty and Service Offering: One year

Certifications: CE marked

Recommended for: Low throughput pipetting procedures; genomics, proteomics, cell based assays, IC50, EC50, microplate stamping, low volume microplate filling and microplate reformatting.

| Cat. No. | Description |
|-------------|---|
| 1096-DT-100 | Hydra DT Base Unit; 96 Channel; 0.5-100 μL |
| 1096-DT-300 | Hydra DT Base Unit; 96 Channel; 5-300 μL |

Thermo Scientific Versette Pipetting Workstation

The ultimate in versatility and scalability, our Versette* automated liquid handler is ideal for laboratories that want to move from handheld pipetting to automated handling.

The Versette automated liquid handling system is compact, user-scalable, and offers many pipetting options. It features 19 self-interchanging and quick-swap, single- through 384-channel pipetting heads, and a total volume range of 0.1 to 1250 μ L. With a compact modular design, the Versette liquid handler offers end users a choice between a two- or six-position stage, and two user-friendly programming options.

Details

- Compatible with 19 interchangeable pipetting heads, ranging from single- to 384-channels
- Pipetting heads are available in disposable and fixed tip formats
- Total volume range from 0.1 to 1250 μ L
- All pipetting heads include RFID tags to self-identify, track usage and store service information

Consumables Designed to Optimize Liquid Handling Performance

- Single-, 8-, and 12-channel pipetting tools utilize Thermo Scientific Versette ClipTip automation tips that securely seal to the pipetting head
- Versette ClipTip automation tips require low insertion and ejection force, extending the life of the instrument's components
- 96- and 384-channel pipetting heads utilize Thermo Scientific D.A.R.T.s* tips
- D.A.R.T.s tips utilize a unique surface seal design that delivers accurate and precise pipetting by forming an even seal across all pipetting channels

Options for Two- or Six-position Deck

- Two user-upgradeable labware capacity options are available for simple liquid handling procedures, medium throughput, and use as an integrated component

Space Matters

- Compact design allows it to fit on standard laboratory benches and into select biological hoods
- Two- and six-position units are equipped with a safety shield
- The six-position stage uses a dual-level design to further minimize space requirements

User-Friendly Programming Options

- Compatible with an onboard GUI interface or PC-based ControlMate* software
- Quickly and easily create or edit simple or complex pipetting procedures
- Ideal for procedures as diverse as automated serial dilution and hit picking to 96/384 plate replications



Recommended for: Can be used in many applications across academia, pharma and biotech, including drug discovery, cell-based, genomic and proteomic procedures. Ideal for life science researchers in low- to high-throughput laboratories who are currently performing handheld pipetting tasks but want to move to automation, or customers building integrated systems that require liquid handling.

Warranty: One year

| Cat. No. | Description | Dimensions |
|-------------|--|---|
| 65001BS | Versette base unit; requires pipetting head housing, stage, pipetting head(s) and accessories separately available | Included below with stage dimensions |
| 650-02-SMC | Single- and multi-head housing assembly for use with single-, 8-, 12-channel pipetting heads | Included below, individual weight 30 lb. (13.6 kg) |
| 650-02-NTC | 96- and 384- head housing assembly for use with 96- and 384-channel pipetting heads | Included below, individual weight 40 lb. (18.1 kg) |
| 650-03-TPS | Two-position stage, guarding included | 26.6" \times 21.4" \times 16.1" (67.6 \times 54.4 \times 40.9 cm); 20 lb. (9.1 kg) |
| 650-03-TPSR | Two-position robotic-friendly stage | 26.6" \times 21.4" \times 16.1" (67.6 \times 54.4 \times 40.9 cm); 15 lb. (6.8 kg) |
| 650-03-SPS | Six-position stage, guarding included (can also be removed for robotic use) | 26.6" \times 21.7" \times 26.7" (67.6 \times 55.1 \times 67.8 cm); 50 lb. (22.7 kg) |

Thermo Scientific Versette Pipetting Heads



Versette Pipetting Heads are available in disposable and fixed-tip formats with a total volume range of 0.1 µL to 1250 µL.

All pipetting heads include RFID tags to self-identify for tracking usage and service information.

Details

- Single- through 384-channel pipetting options
- Disposable and fixed tip options
- Pipetting heads are user interchangeable
- Single-, 8- and 12-channel models are self-interchanging and stored on instrument stage
- 96- and 384-channel models require manual swap that takes only minutes
- Thermo Scientific Matrix PlateMate pipetting heads can be upgraded to be compatible with Versette liquid handlers

Warranty: One year



| Cat. No. | Volume Range | Accuracy | Precision |
|---|----------------|------------------|------------------|
| Single Channel Air Displacement Pipetting Head | | | |
| 650-07-S12 | 0.5 to 12.5 µL | 1.0% or 0.42 µL | 0.6% or 0.023 µL |
| 650-07-S30 | 2 to 30 µL | 1.0% or 0.084 µL | 0.5% or 0.048 µL |
| 650-07-S125 | 10 to 125 µL | 1.0% or 0.334 µL | 0.4% or 0.14 µL |
| 650-07-S300 | 20 to 300 µL | 0.9% or 0.72 µL | 0.3% or 0.3 µL |
| 650-07-S1250 | 100 to 1250 µL | 0.8% or 2.5 µL | 0.2% or 0.8 µL |
| 8-Channel Air Displacement Pipetting Head | | | |
| 650-07-M812 | 0.5 to 12.5 µL | 2.0% or 0.084 µL | 1.2% or 0.063 µL |
| 650-07-M830 | 2.0 to 30 µL | 2.0% or 0.168 µL | 1.0% or 0.096 µL |
| 650-07-M8300 | 20 to 300 µL | 1.0% or 1.44 µL | 1.0% or 0.096 µL |
| 12-Channel Air Displacement Pipetting Head | | | |
| 650-07-M1212 | 0.5 to 12.5 µL | 2.0% or 0.084 µL | 1.2% or 0.063 µL |
| 650-07-M1230 | 2 to 30 µL | 2.0% or 0.168 µL | 1.0% or 0.096 µL |
| 650-07-M12300 | 20 to 300 µL | 1.0% or 1.44 µL | 1.0% or 0.096 µL |
| 96-Channel Air Displacement Pipetting Head | | | |
| 650069630 | 0.5 to 30 µL | 2.0% or 0.15 µL | 1.5% or 0.1 µL |
| 6500696300 | 5 to 300 µL | 2.0% or 1.0 µL | 1.5% or 0.75 µL |
| 384-Channel Air Displacement Pipetting Head | | | |
| 6500638430 | 0.5 to 30 µL | 2.0% or 0.15 µL | 1.5% or 0.1 µL |
| 65006384100 | 1 to 100 µL | 2.0% or 0.5 µL | 1.5% or 0.25 µL |
| 96-Channel Positive Displacement Piercing Pipetting Head | | | |
| 650069650SS [†] | 0.1 to 50 µL | 2.5% or 0.01 µL | 1.5% or 0.2 µL |
| 96-Channel Positive Displacement Pipetting Head | | | |
| 650069650DF ^{††} | 0.1 to 50 µL | 2.5% or 0.01 µL | 1.5% or 0.2 µL |
| 384-Channel Positive Displacement Pipetting Head | | | |
| 6500638450SS [‡] | 0.1 to 50 µL | 2.5% or 0.01 µL | 1.5% or 0.2 µL |
| 6500638450DF ^{††} | 0.1 to 50 µL | 2.5% or 0.01 µL | 1.5% or 0.2 µL |

[†] Stainless-steel Ceramic Coated
^{††} Duraflex
[‡] Stainless-steel PTFE Coated

Thermo Scientific Versette Accessories

Versette Accessories are for use with the Versette automated liquid handling system, enabling an effective solution.

Warranty: One year

| Cat. No. | Description |
|----------------|---|
| 650-05-96TTW | 96 Channel Tip Wash Station, Tall Height |
| 650-05-384TTW | 384 Channel Tip Wash Station, Tall Height |
| 650-08-VMJP | Vacuum Manifold Kit, 110V |
| 650-08-VMNA | Vacuum Manifold Kit, 100V |
| 650-08-VMEU | Vacuum Manifold Kit, 220V |
| 650-08-PRC | Piercing Manifold Kit |
| 650-08-1.5ML-L | 1.5 mL Lidded Tube and Vial Adapter |
| 650-08-2ML | 2.0 mL Tube and Vial Adapter |
| 650-08-8ML | 8.0 mL Tube and Vial Adapter |
| 650-08-BCR | Linear Barcode Reader |
| 650-08-13MM | Tube Adapter, 13 mm Vial |
| 650-08-WSTE | Versette ClipTip Automation Tips Disposal Chute and Container |
| 650-04-PUMP | Pump module 21.2" x 8.7" x 8.7" (54.6 x 22.1 x 22.1cm); 22 lb. (10 kg) |



Ask your representative about integrating the Thermo Scientific Orbiter RS Microplate Mover



Thermo Scientific Versette ClipTip Automation Tips



With a unique interface that locks the tip onto the Versette automated liquid handling system, Versette ClipTip Automation Tips provide stable, straight mounting every time. Each tip has a clasp or 'clip' that locks in place to ensure an effective seal, and will only release when mechanically ejected.

Details

- Unique design: On or Off only – tips lock in place, providing a secure fit that prevents loose, uneven or leaky tips
- Compatible with single-, 8-, and 12- channel Versette pipetting head options
- Range of versatile tips with three sizes for volumes: 1250 μL , 300/125 μL and 30/12.5 μL



| Tip Volume | Range | Non-Sterile | Sterile | Sterile, Filtered | Single Channel Air Displacement Pipetting Head | 8-Channel Air Displacement Pipetting Head | 12-Channel Air Displacement Pipetting Head |
|-----------------------|------------------------|-------------|---------|-------------------|--|---|--|
| 12.5/30 μL | 0.5-12.5 μL | 5606 | 5607 | n/a | 650-07-S12 | 650-07-M812 | 650-07-M1212 |
| 12.5 μL | 2-30 μL | 5606 | 5607 | 5608 | 650-07-S30 | 650-07-M830 | 650-07-M1230 |
| 30 μL | 2-30 μL | 5606 | 5607 | 5618 | 650-07-S30 | 650-07-M830 | 650-07-M1230 |
| 125/300 μL | 10-125 μL | 5706 | 5707 | 5708 | 650-07-S125 | 650-07-M8300 | 650-07-M12300 |
| 125/300 μL | 20-300 μL | 5706 | 5707 | 5708 | 650-07-S300 | 650-07-M8300 | 650-07-M12300 |
| 1250 μL | 100-1250 μL | 5806 | 5807 | 5808 | 650-07-S1250 | n/a | n/a |

All tips packed 10 racks of 96. Compatible with tubes, 96- and 384-well plates. 1250 μL tips packed 8 racks of 96.

Thermo Scientific Matrix D.A.R.T.s Tips

Matrix D.A.R.T.s (Disposable Automation Research Tips) have been designed to provide the ultimate performance in automated pipetting.

Unlike conventional automation tips, Matrix D.A.R.T.s tips are provided in disposable, rigid magazines that bond onto the pipetting heads by sealing against a silicone pad. This forms a definitive seal that minimizes tip-to-tip height variations and eliminates static interactions, delivering consistent dispensing into the bottom of microplate wells.

By creating an extremely consistent seal, Matrix D.A.R.T.s tips offer unmatched consistency when aspirating and dispensing into the bottom of microplate wells.

Details

- Compatible with the Versette, Hydra DT, PlateMate 2x3 and PlateMate Plus Automated Pipetting Systems
- Stepped-body design minimizes volume displacement and provides tip straightness for better well access
- Highly polished molds and cores ensure low liquid retention
- Small internal tip orifice provides better small volume droplet formation, resulting in optimal accuracy and precision
- Sterile, sterile filtered, wide-bore and extended-length tip options available, offering the flexibility to choose tips for any application
- Thorough inspection for proper formation as well as straightness eliminates defects
- Unique face seal design provides consistent tip-to-tip height for precise low volume dispensing
- Entire assembly can be quickly removed and replaced, eliminating extra handling of the tips that can cause damage or contamination

Includes: 20 magazines per case



| Cat. No. | Description | Volume [Tip] | Sterile | Packaging |
|----------|-------------------------------------|--------------|---------|------------------------------|
| 5586 | 96 Tips | 30 µL | No | 96 Tips/Rack, 20 Racks/Case |
| 5587 | 96 Tips | 30 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5506 | 96 Tips, Extended Length | 30 µL | No | 96 Tips/Rack, 20 Racks/Case |
| 5507 | 96 Tips, Extended Length | 30 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5526 | 96 Tips | 100 µL | No | 96 Tips/Rack, 20 Racks/Case |
| 5527 | 96 Tips | 100 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5516 | 96 Tips | 300 µL | No | 96 Tips/Rack, 20 Racks/Case |
| 5517 | 96 Tips | 300 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5546 | 96 Tips, Extended Length, Wide Bore | 300 µL | No | 96 Tips/Rack, 20 Racks/Case |
| 5547 | 96 Tips, Extended Length, Wide Bore | 300 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5536 | 96 Tips, Extended Length | 300 µL | No | 96 Tips/Rack, 20 Racks/Case |
| 5537 | 96 Tips, Extended Length | 300 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5301 | 384 Tips | 12.5 µL | No | 384 Tips/Rack, 10 Racks/Case |
| 5302 | 384 Tips | 12.5 µL | Yes | 384 Tips/Rack, 10 Racks/Case |
| 5416 | 384 Tips, Extended Length | 30 µL | No | 384 Tips/Rack, 20 Racks/Case |
| 5417 | 384 Tips, Extended Length | 30 µL | Yes | 384 Tips/Rack, 20 Racks/Case |
| 5316 | 384 Tips | 30 µL | No | 384 Tips/Rack, 20 Racks/Case |
| 5317 | 384 Tips | 30 µL | Yes | 384 Tips/Rack, 20 Racks/Case |
| 5326 | 384 Tips | 100 µL | No | 384 Tips/Rack, 20 Racks/Case |
| 5327 | 384 Tips | 100 µL | Yes | 384 Tips/Rack, 20 Racks/Case |
| 5321 | 384 Tips | 100 µL | No | 384 Tips/Rack, 10 Racks/Case |
| 5322 | 384 Tips | 100 µL | Yes | 384 Tips/Rack, 10 Racks/Case |

Thermo Scientific Matrix Filtered D.A.R.T.s Tips



Matrix Filtered D.A.R.T.s (Disposable Automation Research Tips) feature a stepped body design that minimizes volume displacement and provides better well access.

Unlike conventional automation tips, our Matrix filtered D.A.R.T.s tips are provided in disposable rigid magazines and bond onto the pipetting heads by being squeezed against a silicone pad. This forms a definitive seal that minimizes tip-to-tip height variations, delivering consistent dispensing into the bottom of microplate wells.

Details

- Compatible with Versette, Hydra DT, PlateMate 2x3, and PlateMate Plus systems
- Stepped body design minimizes volume displacement and ensures tip straightness for better well access
- Small internal tip orifice provides better small-volume droplet formation, resulting in optimal accuracy and precision
- Filter, wide-bore, and extended-length tip options available; offers the flexibility to choose tips for any application
- 100% inspection ensures the tips are free of defects; inspected for proper formation as well as straightness
- Unique face seal design provides consistent tip-to-tip height for precise low volume dispensing
- Disposable magazine eliminates the need to flip tips out of a traditional box and into a metal carrier
- Entire assembly can be quickly removed and replaced; eliminates extra handling of the tips that can cause damage or contamination


| Cat. No. | Description | Volume [Tip] | Sterile | Packaging |
|----------|-------------------------------------|--------------|---------|------------------------------|
| 5508 | 96 Tips, Extended Length | 30 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5588 | 96 Tips | 30 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5526 | 96 Tips | 100 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5518 | 96 Tips | 300 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5538 | 96 Tips, Extended Length | 300 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5548 | 96 Tips, Extended Length, Wide Bore | 300 µL | Yes | 96 Tips/Rack, 20 Racks/Case |
| 5318 | 384 Tips | 30 µL | Yes | 384 Tips/Rack, 20 Racks/Case |
| 5418 | 384 Tips, Extended Length | 30 µL | Yes | 384 Tips/Rack, 20 Racks/Case |
| 5328 | 384 Tips | 100 µL | Yes | 384 Tips/Rack, 20 Racks/Case |

PocketTip D.A.R.T.s (Disposable Automation Research Tips) offer low-volume transfer capability with existing Thermo Scientific automated liquid handling equipment, ideal for compound screening procedures.

PocketTip D.A.R.T.s enable the transfer between 50-250 nL of sample volume using Versette, PlateMate and Hydra DT liquid handling instruments. These low-volume transfers reduce the cost of reagents, labware and sample preparation while improving data quality by eliminating intermediate pipetting steps and providing soft in-tip mixing.

Details

- Available in 96- and 384- channel format with 50, 100 and 250 nL pocket sizes
- Sterile format is available for sensitive liquid handling procedures or operations that require sterility
- Enables existing Thermo Scientific automated liquid handlers to accommodate low-volume dispensing for compound screening procedures
- Requires only minor procedural adjustments to Versette, PlateMate and Hydra DT instrumentation
- By transferring nanoliter volumes directly to assay plates, timely aqueous dilution steps that can cause precipitation and result in unreliable data are eliminated
- Sample and assay materials are gently combined with soft in-tip mixing, removing the need for offline procedures and improving data quality with thorough mixing that will not damage sensitive material
- Consistent accuracy and precision are tightly controlled during the molding process
- Purity and straightness validation is also performed on compatible instrumentation to ensure only the highest quality
- With instrument compatibility range starting from the Versette, Hydra DT on-the-fly dose response curves or assay development can be set up with simple Hydra DT programming
- Higher throughput needs can be met with the Versette, PlateMate 2x3 and CatalySt Express, or PlateMate Plus instruments
- Automation sales specialists are available to assist with instrument programming and modifications, as well as with instrument validation regarding any Thermo Scientific automated liquid handler



| Cat. No. | Volume [Tip] | Sterile | Packaging |
|----------|--------------|---------|-----------------|
| 5696-100 | 100 nL | No | Racks, 20 × 96 |
| 5696-250 | 250 nL | No | Racks, 20 × 96 |
| 5696-50 | 50 nL | No | Racks, 20 × 96 |
| 5697-100 | 100 nL | Yes | Racks, 20 × 96 |
| 5697-250 | 250 nL | Yes | Racks, 20 × 96 |
| 5697-50 | 50 nL | Yes | Racks, 20 × 96 |
| 5756-100 | 100 nL | No | Racks, 20 × 384 |
| 5756-250 | 250 nL | No | Racks, 20 × 384 |
| 5756-50 | 50 nL | No | Racks, 20 × 384 |
| 5757-100 | 100 nL | Yes | Racks, 20 × 384 |
| 5757-250 | 250 nL | Yes | Racks, 20 × 384 |
| 5757-50 | 50 nL | Yes | Racks, 20 × 384 |

Thermo Scientific Matrix D.A.R.T.s Tip Transfer Tool



The Matrix D.A.R.T.s Tip Transfer Tool simplifies operations and increases efficiency when loading Matrix D.A.R.T.s tips into serial dilution or custom configuration magazines used with the Thermo Scientific Versette, PlateMate Plus, PlateMate 2 x 3, and Hydra DT.

Using this tool to transfer Matrix D.A.R.T.s tips from a magazine into serial dilution or custom magazines is fast and efficient. An entire row or column can be transferred in one step.

The tool supports a variety of tip configurations on Thermo Scientific automated liquid handling instruments. It simplifies and increases efficiency when loading Matrix D.A.R.T.s into serial dilution or custom configuration magazines, and can be easily modified to work with either 96- or 384-format Matrix D.A.R.T.s.

Details

- Compatible with 96- and 384-format Matrix D.A.R.T.s tips, the tool is easily modified to work with either one
- Supports a variety of tip configurations on Thermo Scientific Versette, PlateMate and Hydra DT instruments
- Simply arrange the 96- or 384-tip fittings into desired configuration, such as 4.5 mm or 9 mm spacing
- Easy assembly offers fast swapping between 96- and 384-format tip fittings
- Operates like a pipetter by picking up tips on tip fittings, and using the plunger to eject tips into serial dilution or custom-configured magazine
- The tool is lightweight and requires low force ejection for easy tip transfers

Includes: Transfer tool, tip fittings for 96- and 384-format, and two stripper plates

| Part Number | Description |
|-------------|---|
| 501-30045 | Tip Transfer Tool; Matrix DARTs |
| 501-30046 | Tip Transfer Tool Fittings; 96 channel |
| 501-30047 | Tip Transfer Tool Fittings; 384 channel |





Thermo Scientific Matrix Disposable Automation Reservoirs

Matrix Disposable Automation Reservoirs feature the same shape and dimensions as standard microplates, and are designed to reduce reagent waste.

Now you have the flexibility to handle reagents on virtually any automated liquid handler. Hydrophilic surface treatment results in a lessened meniscus, allowing the liquid to spread out evenly along the reservoir bottom. The result: Less liquid is needed to cover the bottom and allows access to the liquid by all pipette tips.

With the same dimensions as standard microplates, these disposable reservoirs can be used in plate handlers or stackers. The removable splash baffles prevent reagents from splashing out when being moved on robotic decks. They are compatible with virtually any automated liquid handler.

Details

- Economical and disposable, they eliminate the risk of cross-contamination and the inconvenience of washing and autoclaving
- Standard microplate dimensions; can be used in plate handlers or stackers
- Built-in grooves in reservoir bottom reduce reagent waste and allow consistent aspiration of minimal amounts of reagent
- Hydrophilic surface allows the liquid to spread out evenly, so minimal reagent is needed to cover reservoir bottom
- Polypropylene construction offers excellent chemical resistance
- Individually wrapped reservoirs stay sterile until ready for use; no need to autoclave
- Built-in pour spout allows easy emptying of reagents

Removable Splash Baffles

- Prevent reagents from splashing out when being moved on robotic decks
- Enables reservoirs to be moved faster, and to be filled closer to the top with reagent
- Baffles can be removed to increase reservoir's capacity



| Cat. No. | Description | Mfr. No. |
|--------------------------------------|----------------------|----------|
| Nonsterile | | |
| 1064058 | 96 Channel | 1064058 |
| 1064056 | 96 Channel, Deepwell | 1064056 |
| 1064057 | 384 channel | 1064057 |
| Sterile, Individually Wrapped | | |
| 1064158 | 96 Channel | 1064158 |
| 1064156 | 96 Channel, Deepwell | 1064156 |
| 1064157 | 384 Channel | 1064157 |

Molecular BioProducts BioRobotix ART Filter Tips for Automated Workstations



Molecular BioProducts BioRobotix ART Filter Tips include a filter that prevents contaminating liquids and aerosols from being drawn into your workstation.

These sterile tips prevent contaminating liquids and aerosols from being drawn into your workstation, helping to maintain integrity of samples and reducing the risk of cross-contamination.


Designed for maximum accuracy and repeatability, BioRobotix ART filter pipette tips are available in two colors: Natural (clear) and black (liquid sensing). Constructed of premium-grade virgin polypropylene, all tips are rigorously tested on the workstations for which they were designed.

Details

- Compatible with a wide range of automated workstations
- Designed for maximum accuracy and repeatability
- Available in two colors: Natural (clear) and black (liquid sensing)
- Constructed of premium-grade, virgin polypropylene
- Tips are rigorously tested on the workstations for which they were designed

| Cat. No. | Description – Validated Workstations |
|--------------|--|
| Black | |
| 901-011 | Tecan 10 μ L Conductive Sterile Filtered |
| 902-011 | Tecan 50 μ L Conductive Sterile Filtered |
| 903-011 | Tecan 175 μ L Conductive Sterile Filtered |
| 903-011G | Tecan 175 μ L Conductive Sterile Filtered Wide Bore |
| 903-033 | Tecan 175 μ L Conductive Filtered Econopack |
| 904-011 | Tecan 1000 μ L Conductive Sterile Filtered |
| 904-011G | Tecan 1000 μ L Conductive Sterile Filtered Wide Bore |
| 904-013 | Tecan 1000 μ L Conductive Sterile Filtered Econo Pack |
| 918-011 | Beckman Biomek FX SPAN-8 20 μ L Conductive Sterile Filtered |
| 919-011 | Beckman Biomek FX SPAN-8 130 μ L Conductive Sterile Filtered |
| 951-011 | Qiagen 250 μ L Conductive Sterile Filtered |
| 953-011 | Qiagen 1100 μ L Conductive Sterile Filtered |
| 953-011G | Qiagen 1100 μ L Conductive Sterile Filtered Wide Bore |

| Cat. No. | Description – Validated Workstations |
|----------------|--|
| Natural | |
| 901-021 | Tecan 10 µL Clear Sterile Filtered |
| 902-021 | Tecan 50 µL Clear Sterile Filtered |
| 903-021 | Tecan 175 µL Clear Sterile Filtered |
| 904-021 | Tecan 1000 µL Clear Sterile Filtered |
| 906-021 | Tecan Te-MO 50 µL Clear Sterile Filtered |
| 907-021 | Tecan Te-MO 100 µL Clear Sterile Filtered |
| 908-021 | Tecan Te-MO 150 µL Clear Sterile Filtered |
| 912-021 | Beckman Biomek FX 10 µL Clear Sterile Filtered |
| 915-021 | Beckman Biomek 20 µL Clear Sterile Filtered |
| 915-021-05 | Beckman Biomek 20 µL Clear Sterile Filtered Low Retention |
| 916-021 | Beckman FX Extended Length 40 µL Clear Sterile Filtered |
| 916-021-05 | Beckman FX Extended Length 40 µL Clear Sterile Filtered Low Retention |
| 917-021 | Beckman Biomek 130 µL Clear Sterile Filtered |
| 917-021-05 | Beckman Biomek 130 µL Clear Low Retention |
| 917-021G | Beckman Biomek 130 µL Clear Wide Bore |
| 918-021 | Beckman Biomek FX 20 µL Clear Sterile Filtered |
| 918-021-05 | Beckman Biomek FX SPAN-8 20 µL Clear Sterile Filtered Low Retention |
| 919-021 | Beckman Biomek FX 130 µL Clear Sterile Filtered |
| 919-021-05 | Beckman Biomek FX 130 µL Clear Sterile Filtered Low Retention |
| 919-021G | Beckman Biomek FX 130 µL Clear Sterile Filtered Wide Bore |
| 920-021 | Caliper Zymark 50 µL Clear Sterile Filtered |
| 921-021 | Caliper Zymark 80 µL Clear Sterile Filtered |
| 923-021 | Caliper Zymark 100 µL Clear Sterile Filtered |
| 925-021 | PerkinElmer Plate Trak 20 µL Clear Sterile Filtered |
| 925-021-05 | PerkinElmer Plate Trak 20 µL Clear Sterile Filtered |
| 927-021 | PerkinElmer Plate Trak 50 µL Clear Sterile Filtered |
| 927-021-05 | PerkinElmer Plate Trak 50 µL Clear Sterile Filtered |
| 929-021 | PerkinElmer Plate Trak 235 µL Clear Sterile Filtered |
| 929-021-05 | PerkinElmer Plate Trak 235 µL Clear Sterile Filtered |
| 931-021 | PerkinElmer CCS/FLIPR Type A 384 25 µL Clear Sterile Filtered |
| 935-021 | Beckman 384 25 µL Clear Sterile Filtered |
| 935-021-05 | Beckman 384 25 µL Clear Sterile Filtered Low Retention |
| 938-021 | Beckman 384 Extended Length 40 µL Clear Sterile Filtered |
| 938-021-05 | Beckman 384 Extended Length 40 µL Clear Sterile Filtered Low Retention |
| 951-021 | Qiagen 250 µL Clear Sterile Filtered |
| 953-021 | Qiagen 1100 µL Clear Sterile Filtered |
| 984-021 | Biotek Precision SBS Format 50 µL Clear Sterile Filtered |
| 985-021 | Biotek Precision SBS Format 100 µL Clear Sterile Filtered |
| 986-021 | Biotek Precision SBS Format 150 µL Clear Sterile Filtered |



Molecular BioProducts BioRobotix Pipet Tips, Standard, Black



Molecular BioProducts BioRobotix standard black pipette tips are constructed of high quality carbon-impregnated polypropylene for liquid sensing.

Designed for maximum accuracy, repeatability and the ability to detect liquids, BioRobotix black polypropylene tips are compatible with a variety of liquid sensing automated workstations. Constructed of premium-grade virgin polypropylene with embedded carbon, all tips are rigorously tested on the systems for which they were designed.

Details

- Compatible with a wide range of liquid-sensing automated workstations from multiple vendors
- Black (liquid-sensing) tips are designed for maximum accuracy, repeatability and liquid-sensing capability
- Constructed of premium-grade, carbon-embedded virgin polypropylene
- Tips are rigorously tested on the systems for which they were designed

| Cat. No. | Description – Validated Workstations |
|----------|--|
| 901-251 | Tecan 20 μ L Conductive Sterile |
| 901-252 | Tecan 20 μ L Conductive Non Sterile |
| 901-253 | Tecan 20 μ L Conductive Econo Pack |
| 902-251 | Tecan 50 μ L Conductive Sterile |
| 902-252 | Tecan 50 μ L Conductive Non Sterile |
| 902-253 | Tecan 50 μ L Conductive Econo Pack |
| 903-033 | Tecan 175 μ L Conductive Filtered Econopack |
| 903-251 | Tecan 200 μ L Conductive Sterile |
| 903-251G | Tecan 200 μ L Conductive Sterile Wide Bore |
| 903-252 | Tecan 200 μ L Conductive Non Sterile |
| 903-252G | Tecan 200 μ L Conductive Wide Bore Non Sterile |
| 903-253 | Tecan 200 μ L Conductive Non Sterile Econo Pack |
| 904-251 | Tecan 1000 μ L Conductive Sterile |
| 904-251G | Tecan 1000 μ L Conductive Sterile Wide Bore |
| 904-252 | Tecan 1000 μ L Conductive Non Sterile |
| 904-252G | Tecan 1000 μ L Conductive Wide Bore |
| 904-253 | Tecan 1000 μ L Conductive Econo Pack |
| 916-251 | Beckman FX Extended Length 50 μ L Conductive Sterile |
| 916-252 | Beckman FX Extended Length 50 μ L Conductive Non Sterile |
| 918-251 | Beckman Biomek FX SPAN-8 20 μ L Conductive Sterile |
| 918-252 | Beckman Biomek FX SPAN-8 20 μ L Conductive Non Sterile |
| 919-251 | Beckman Biomek FX SPAN-8 250 μ L Conductive Sterile |
| 919-252 | Beckman FX SPAN-8 250 μ L Conductive Non Sterile |
| 951-251 | Qiagen 300 μ L Conductive Sterile |
| 951-252 | Qiagen 300 μ L Conductive Non Sterile |
| 951-252G | Qiagen 300 μ L Conductive Wide Bore |
| 951-253 | Qiagen 300 μ L Conductive Econo Pack |
| 953-251 | Qiagen 1100 μ L Conductive Sterile |
| 953-251G | Qiagen 1100 μ L Conductive Sterile Wide Bore |
| 953-252 | Qiagen 1100 μ L Conductive Non Sterile |
| 953-253 | Qiagen 1100 μ L Conductive Econo Pack |

Molecular BioProducts BioRobotix Pipette Tips: Standard, Natural, Nonsterile

Molecular BioProducts BioRobotix clear polypropylene pipette tips are designed for maximum accuracy and repeatability.

Designed for maximum accuracy and repeatability, BioRobotix pipette tips are compatible with a variety of automated workstations from multiple vendors. Constructed of premium-grade virgin polypropylene, all tips are rigorously tested on the systems for which they were designed.

| Cat. No. | Description – Validated Workstations |
|------------|--|
| 901-261 | Tecan 20 µL Clear Sterile |
| 901-262 | Tecan 20 µL Clear Non Sterile |
| 901-263 | Tecan 20 µL Clear Econo Pack |
| 902-261 | Tecan 50 µL Clear Sterile |
| 902-261XP | Tecan 50 µL Clear Sterile PATCH XPRESS |
| 902-262 | Tecan 50 µL Clear Non Sterile |
| 902-262XP | Tecan 50 µL Clear PATCH XPRESS |
| 902-263 | Tecan 50 µL Clear Econo Pack |
| 903-261 | Tecan 200 µL Clear Sterile |
| 903-262 | Tecan 200 µL Clear Non Sterile |
| 903-263 | Tecan 200 µL Clear Non Sterile Econo Pack |
| 904-261 | Tecan 1000 µL Clear Sterile |
| 904-262 | Tecan 1000 µL Clear Non Sterile |
| 904-263 | Tecan 1000 µL Clear Non Sterile Econo Pack |
| 906-261 | Tecan Te-MO 50 µL Clear Sterile |
| 906-262 | Tecan Te-MO 50 µL Clear Non Sterile |
| 907-261 | Tecan Te-MO 100 µL Clear Sterile |
| 907-262 | Tecan Te-MO 100 µL Clear |
| 908-261 | Tecan Te-MO 200 µL Clear Sterile |
| 908-262 | Tecan Te-MO 200 µL Clear Non Sterile |
| 912-261 | Beckman Biomek FX 10 µL Clear Sterile |
| 912-262 | Beckman Biomek FX 10 µL Clear Non Sterile |
| 915-261 | Beckman Biomek 20 µL Clear Sterile |
| 915-261-05 | Beckman Biomek 20 µL Clear Sterile Low Retention |
| 915-262 | Beckman Biomek 20 µL Clear Non Sterile |
| 915-262-05 | Beckman Biomek 20 µL Clear Low Retention |
| 916-261 | Beckman FX Extended Length 50 µL Clear Sterile |
| 916-261-05 | Beckman FX Extended Length 50 µL Clear Sterile Low Retention |
| 916-262 | Beckman FX Extended Length 50 µL Clear Non Sterile |
| 916-262-05 | Beckman FX Extended Length 50 µL Clear Low Retention |
| 917-261 | Beckman Biomek 250 µL Clear Sterile |
| 917-261-05 | Beckman Biomek 200 µL Clear Sterile Low Retention |
| 917-261G | Beckman Biomek 200 µL Clear Sterile Wide Bore |
| 917-262 | Beckman Biomek 250 µL Clear Non Sterile |
| 917-262-05 | Beckman Biomek 200 µL Clear Low Retention |
| 918-261 | Beckman Multimek 20 µL Clear Sterile |
| 918-261-05 | Beckman Biomek FX SPAN-8 20 µL Clear Sterile Low Retention |
| 918-262 | Beckman Multimek 20 µL Clear Non Sterile |
| 918-262-05 | Beckman Biomek FX SPAN-8 20 µL Clear Low Retention |
| 919-261 | Beckman Biomek 250 µL Clear Sterile |
| 919-261-05 | Beckman Biomek FX 250 µL Clear Sterile Low Retention |
| 919-261G | Beckman Biomek FX 250 µL Clear Sterile Wide Bore |



Details

- Compatible with a wide range of automated workstations
- Constructed of premium-grade virgin polypropylene
- Tips are rigorously tested on the systems for which they were designed
- Packaged in trays

*Molecular BioProducts BioRobotix Pipette Tips
>> continued on next page*

Molecular BioProducts BioRobotix Pipette Tips: Standard, Natural, Nonsterile (continued)

| Cat. No. | Description – Validated Workstations |
|------------|---|
| 919-262 | Beckman Multimek 250 µL Clear Non Sterile |
| 919-262-05 | Beckman Biomek FX 250 µL Clear Low Retention |
| 919-262G | Beckman Biomek FX 250 µL Clear Wide Bore |
| 920-261 | Zymark 50 µL Clear Sterile |
| 920-262 | Zymark 50 µL Clear Non Sterile |
| 921-261 | Zymark 100 µL Clear Sterile |
| 921-262 | Zymark 100 µL Clear Non-Sterile |
| 923-261 | Zymark 200 µL Clear Sterile |
| 923-262 | Zymark 200 µL Clear Non Sterile |
| 925-261 | Perkin-Elmer Plate Trak 20 µL Clear Sterile |
| 925-261-05 | Perkin-Elmer Plate Trak 20 µL Clear Sterile |
| 925-262 | Perkin-Elmer Plate Trak 20 µL Clear |
| 925-262-05 | Perkin-Elmer Plate Trak 20 µL Clear |
| 927-261 | Perkin-Elmer Plate Trak 50 µL Clear Sterile |
| 927-261-05 | Perkin-Elmer Plate Trak 50 µL Clear Sterile |
| 927-262 | Perkin-Elmer Plate Trak 5 µL Clear |
| 927-262-05 | Perkin-Elmer Plate Trak 5 µL Clear |
| 929-261 | Perkin-Elmer Plate Trak 235 µL Clear Sterile |
| 929-261-05 | Perkin-Elmer Plate Trak 235 µL Clear Sterile |
| 929-262 | Perkin-Elmer Plate Trak 235 µL Clear |
| 929-262-05 | Perkin-Elmer Plate Trak 235 µL Clear |
| 931-261 | CCS/FLIPR Type A 384 30 µL Clear Sterile |
| 931-262 | CCS/FLIPR Type A 384 30 µL Clear Non Sterile |
| 935-261 | Beckman 384 30 µL Clear Sterile |
| 935-261-05 | Beckman 384 30 µL Clear Sterile Low Retention |
| 935-262 | Beckman 384 30 µL Clear Non Sterile |
| 935-262-05 | Beckman 384 30 µL Clear Low Retention |
| 937-261 | FLIPR Liberty Type B 384 30 µL Clear Sterile |
| 937-262 | FLIPR Liberty Type B 384 30 µL Clear Non Sterile |
| 938-261 | Beckman 384 Extended Length 50 µL Clear Sterile |
| 938-261-05 | Beckman 384 Extended Length 50 µL Clear Sterile Low Retention |
| 938-262 | Beckman 384 Extended Length 50 µL Clear |
| 938-262-05 | Beckman 384 Extended Length 50 µL Clear Low Retention |
| 946-261 | Velocity11 384 10 µL Clear Sterile |
| 946-262 | Velocity11 384 10 µL Clear Non Sterile |
| 947-261 | Velocity11 384 30 µL Clear Sterile |
| 947-262 | Velocity11 384 30 µL Clear Non Sterile |
| 948-261 | Velocity11 384 50 µL Clear Sterile |
| 948-262 | Velocity11 384 50 µL Clear Non Sterile |
| 949-261 | Velocity11 384 70 µL Clear Sterile |
| 949-262 | Velocity11 384 70 µL Clear Non Sterile |
| 951-261 | Qiagen 300 µL Clear Sterile |
| 951-262 | Qiagen 300 µL Clear Non Sterile |
| 953-261 | Qiagen 1100 µL Clear Sterile |
| 953-262 | Qiagen 1100 µL Clear Non Sterile |
| 984-261 | Biotek Precision SBS Format 50 µL Clear Sterile |
| 984-262 | Biotek Precision SBS Format 50 µL Clear |
| 985-261 | Biotek Precision SBS Format 100 µL Clear Sterile |
| 985-262 | Biotek Precision SBS Format 100 µL Clear |
| 986-261 | Biotek Precision SBS Format 200 µL Clear Sterile |
| 986-262 | Biotek Precision SBS Format 200 µL Clear |





Thermo Scientific RapidStak Automated Microplate Stacker with Polara RS Software

The RapidStak Automated Microplate Stacker can be set up in minutes, providing walkaway time and relief from monotonous tasks.

Fast and reliable, the RapidStak Automated Microplate Stacker features dual-instrument loading for optimum versatility. Plates can be quickly delivered and processed using two instruments simultaneously. The RapidStak stacker provides valuable walkaway time and relief from monotonous tasks. Setup takes only minutes with most Thermo Scientific and third party instruments, even in the field. It's the perfect partner for your Multidrop bulk reagent dispenser product line, or Thermo Scientific ALPS line of automated sealers.

Details

- Rapid, high throughput performance
- High speed is made possible via a unique buffer-nest technology that enables parallel processing of stacking and plate loading
- Unique choreography and bidirectional telescoping arm design allow it to deliver plates to the instrument while performing the next stack/destack operation simultaneously
- Multitasking yields a load/unload time of less than six seconds, for up to 300 plates per hour throughput

Quick and Easy Setup

- Rapid, easy setup in the field for the least amount of instrument downtime – no tools required for setup
- No PC required for automation of Multidrop dispensers, Wellwash Versa microplate washers or ALPS Sealers
- No instrument modifications or additional options needed
- Does not interfere with the standalone use of the instruments

Flexible, Expandable

- Accommodates six to 22 mm SBS conforming plates
- Standard model holds two Staks for a capacity of 50 plates (15 mm size; one stack remains empty); can be upgraded to RapidStak 2x units
- RapidStak 2x stackers hold four Staks, handling up to 150 plates for even greater walkaway time
- 30- and 50-plate Staks are easily interchangeable and can be used with all RapidStak models

Automate Two Instruments with Polara RS Software

- Scheduling software transforms any RapidStak stacker into a complete, optimized benchtop assay system
- Lets you program the stacker to perform almost any plate-loading application
- Provides even greater versatility by automating virtually any two instruments simultaneously

Compliance: CSA approved

Certifications: CE marked UL listed

| Cat. No. | Description |
|----------|-------------------------|
| F01350 | RapidStak |
| F01351 | RapidStak 2x |
| F01489 | RapidStak, Shortened |
| F01490 | RapidStak 2x, Shortened |
| F01364 | Upgrade Kit |

| Specifications | |
|-----------------------------------|---|
| Compatible Plate Types | 6 to 22 mm SBS-conforming plates |
| Maximum Plate Capacity (10/15 mm) | Up to 50 |
| Plate Delivery/Retrieval Time | 6 seconds |
| Throughput | Up to 300 plates per hour |
| Interface | Serial RS-232C |
| W × D | Standard unit: 38.1 × 57.5 cm Short unit: 38.1 × 49.5 cm |
| Stack Height | 25 stack: 52.5 cm 30 stack: 60.2 cm 50 stack: 87.9 cm |



Thermo Scientific RapidStak Accessories: Microplate Stacks



RapidStak Accessories: Microplate Stacks are compatible with the RapidStak and RapidStak 2x units.

| Cat. No. | Model |
|----------|----------------|
| F01492 | 25-plate Stack |
| F01362 | 30-plate Stack |
| F01363 | 50-plate Stack |

Thermo Scientific RapidStak Instrument Drivers

RapidStak Instrument Drivers connect an ever-increasing number of laboratory instruments to Polara automated lab systems.

Polara RS Instrument Drivers enable you to easily integrate a wide variety of lab instruments and devices into a Polara RS automated lab system. They support instrument pooling, enabling you to add instruments for higher throughput without having to change your methods. You never have to specify robot operations in your method; Polara RS Instrument Drivers determine those for you.

| Cat. No. | For Use with Thermo Scientific instruments |
|----------|---|
| RSI0012 | ALPS 3000 |
| RSI0013 | ALPS 300 |
| RSI0014 | ALPS 100 |
| RSI0004 | Fluoroskan Ascent, Fluoroskan Ascent FL and Luminoskan Ascent |
| RSI0010 | Microscan MS-3 Barcode Reader |
| RSI0002 | Multidrop 384/Micro/DW |
| RSI0003 | Multidrop Combi |
| RSI0017 | Multidrop Combi nL |
| RSI0001 | Polara RS new instrument interface |
| RSI0011 | Variomag Microplate Shaker |
| RSI0006 | Varioskan Flash (only works with 2-stak RapidStak) |
| RSI0020 | Multiskan FC |
| RSI0022 | Wellwash Versa |
| RSI0015 | Multiskan Spectrum |



Thermo Scientific Polara RS

Polara RS Software is a proprietary software package that transforms any RapidStak Microplate Stacker into a complete, optimized benchtop assay solution. Automate two instruments quickly and easily.


| Cat. No. | Model |
|----------|--------------------------------------|
| F01517 | Polara RS, Single Instrument License |
| F01518 | Polara RS, Dual Instrument License |

Thermo Scientific RapidStak DLL Programming Kits

RapidStak DLL Programming Kits are designed to enable the connection of the RapidStak Microplate Stacker with other instruments.

This tool kit is for software engineers who want to enable any programmable device, such as a laboratory instrument, to control the RapidStak stacker via a serial communications connection.

| Cat. No. | Model |
|----------|-------------|
| F01410 | Single Side |
| F01411 | Dual Side |





***Depend on
the full range
of superior
Thermo Scientific
solutions...***

Advanced design features combined with quality construction, offering years of dependable service

Cold Storage

Researchers worldwide protect more than two billion samples inside Thermo Scientific cold storage equipment. With +4°C refrigerators to -196°C cryogenic freezers and the proven performance of Thermo Scientific Nalgene and Nunc Cryo vials, you're free to concentrate on your work without worrying about your valuable samples.

Thermo Scientific Orbitor RS Microplate Mover

The Orbitor* RS Microplate Mover provides industrial-sized performance within a compact, benchtop mover.

Developed using proven technology from our Thermo Scientific RapidStak and Thermo Scientific Dimension4 product lines, the Orbitor RS high-speed microplate mover offers reliable performance with totally flexible plate handling. Extensive vertical reach allows multiple stacked or high density instruments to be loaded in a small footprint. A bidirectional telescoping arm provides superior reach, improved user safety, and unlimited base rotations within a 360° workspace. The Orbitor RS expands the existing range of Thermo Scientific automated technologies for drug discovery and addresses the demand for increased throughput, storage capacity and operational flexibility.

Details

Reliable

- Collision detection and recovery
- Plate sensing in the gripper
- Servo gripper does not drop a plate when the power goes out

Fast

- Delivers a plate in as little as four seconds
- Bidirectional telescoping arm
- Unlimited 360° rotations
- Motion blending for very smooth movements

Superior Reach and Small Footprint

- Integrate up to four instruments
- Reach radius of ±0.4 meters
- Additional storage capacity (up to nine hotels or stacks)
- Fits in most biosafety cabinets

Flexible Plate Handling

- Supports random or sequential access storage
- Landscape or portrait orientation within same stack or hotel
- Adjustable gripper
- De-lid plates via gripper
- Internal re-grip station

Easy to teach

- Manual teaching
- Easy-to-use interface for teaching
- SmartMove* technology teaches paths around obstacles
- Teach wizards are available for teaching stacks, hotels and instruments

Software

- Programmer's Tool Kit available for OEM opportunities
- Operates with Momentum workflow software
- Error detection and recovery
- Dynamic scheduler ensures efficient instrument use

Safe

- Collision detection automatically stops if the mover senses contact



| Specifications | |
|----------------------|---|
| Nominal Payload | 320 g (11.3 oz) |
| Input Voltage | 110V, 115V or 230V |
| Power Line Frequency | 50/60Hz |
| Weight | 25 kg (55 lb.) |
| Temperature Range | 4 to 40°C (39 to 104°F) |
| Humidity | 80%, noncondensing |
| Storage | 24 to 96 random access; 80 to 440 stacked storage |
| Interface | RS232 |

| Cat. No. | Model |
|--|-------|
| Please contact your Thermo Scientific sales representative | |

Thermo Scientific CataLyst Express Microplate Handler



The CataLyst Express Microplate Handler provides reliable, precise automation of microplate handling for a wide range of instruments and containers.

The Thermo Scientific CataLyst Express handler enables researchers to easily automate small volume assays, such as selectivity or specificity tests. It is the first truly industrial robot to be packaged as a ready-to-use, no-tools-required product for laboratory use. Controlled by Thermo Scientific Momentum workflow software, the CataLyst Express unit transfers microplates between instruments and the included storage units. The only setup and teaching required is physically installing an instrument relative to the CataLyst Express handler and teaching the path to the instrument.

The CataLyst Express hardware consists of a Thermo Scientific CataLyst-5 robotic arm, three microplate storage units and an integrated onboard re-grip station. A wide variety of microplate types is supported. A re-grip station enables the CataLyst Express handler to rotate microplates to either landscape or portrait orientation.

| Specifications | |
|----------------------|--|
| Nominal Payload | 1 kg (2.21 lb.) |
| Input Voltage | 110V, 115V or 230V |
| Power Line Frequency | 50/60Hz |
| Weight | 76 kg (167 lb.) |
| Temperature Range | 10 to 40°C (50 to 104°F) |
| Humidity | 80%, noncondensing |
| Storage | 45 to 285 random access; 80 to 760 stacked storage |
| Interface | RS232 |

| Cat. No. | Model |
|----------|----------------------------|
| F01243 | Standard Configuration |
| F01244 | User-defined Configuration |

Details

High throughput, Extreme flexibility

- Continuous, error-free operation for maximum walk-away time and throughput
- Easily automates virtually any application, including the loading and unloading of complex instruments
- Up to three times faster than competing movers for quicker results
- Random and sequential storage handles a wide range of containers
- Can be easily expanded to accommodate up to 760 stacked plates or 285 random-access plates by using the Thermo Scientific microplate carousel
- Momentum software automates multiple applications and schedules them simultaneously
- Mover Teach user interface makes setup and configuration effortless

Advanced Safety Features

- Remote Emergency Stop control
- Interlocked safety guarding that halts the arm when opened

Industrial-Strength Robotics

- Dexterous articulated robot provides five degrees of freedom, allowing for highly accurate positioning and plate placement
- System's 360° base rotation covers a maximum work area and allows very flexible instrument placement

Ordering Information: Pre-configured PC has six serial ports, 15 digital inputs, 12 digital outputs, one analog input, four relay outputs and two serial ports

Compliance: Meets OSHA requirements; CSA approved

Certifications: CE marked

Thermo Scientific Nalgene Disposable Robotic Reservoirs

Nalgene* Disposable Robotic Reservoirs, featuring a standard microplate format in a 300 mL volume, increase walk-away time and minimize dead volume.

Details

- Polypropylene construction withstands robotic liquid handling, making it easier to run automated assays and increase walk-away time
- Compatible with most robotic arms, plate stackers (hotels), plate cranes and platform configurations
- Reservoir volume is 300 mL (345 mL brim capacity), with molded-in graduations every 100 mL
- Nestable, for easy automated stacking and handling
- Accepts bar-coded labels on two sides
- Available sterile and nonsterile with flat bottoms or with convoluted bottoms that minimize dead volume
- Flat-bottom configuration can be used with single-, 6-, 12-, 96-, 384- and 1536-well formats
- Convoluted-bottom design is most effective with a 96-well format

Ordering Information: Dimensions: 3.4" L x 5 W" x 1.75" H (86 x 128 x 43 mm)

Compliance: Conform to the ANSI* standard microplate footprint



| Cat. No. | Description |
|-----------|---------------------|
| 1200-1300 | Flat |
| 1200-1301 | Flat, Sterile |
| 1200-2300 | Convoluted |
| 1200-2301 | Convoluted, Sterile |



Thermo Scientific Nunc Disposable Plastic Reservoirs



Nunc Disposable Plastic Reservoirs are disposable, batch traceable, and are fully compatible with automated instrumentation.

Featuring a standard microplate format, the revolutionary convoluted design minimizes dead volume and has baffles to reduce splashing. The flat bottom version is designed to meet a variety of assay needs in both manual and automated environments. Our reservoirs feature excellent chemical resistance and low adsorption of organic molecules or cells.

Made from the highest quality plastics, these reservoirs are ideal for use in biotechnology, pharmaceutical and research laboratories, as well as in the production of vaccines and diagnostic kits.

Details

- Made of low-binding polypropylene material
- Highest usable capacity on the market: 300 mL capacity with graduations (345 mL brim)
- Fully compatible with automated instrumentation
- Dimensional compliance – standard microplate footprint
- Sterile version, RNase/DNase-free available
- Polypropylene for low binding
- With or without 96 indentations to reduce minimum required working volume (required volume: 2.1 mL for indented format)
- Baffles to reduce splashing
- Reservoirs are stackable
- Excellent chemical resistance and low adsorption of organic molecules or cells

| Cat. No. | Surface | Sterile |
|----------|----------|---------|
| 370905 | Flat | No |
| 370906 | Flat | Yes |
| 370907 | Indented | No |
| 370908 | Indented | Yes |

Thermo Scientific Matrix Deepwell Storage Blocks

Matrix Deepwell Storage Blocks are constructed of chemical-resistant polypropylene and meet all standard microplate specifications.

Details


- Alphanumeric notations simplify the identification of wells
- Industry standard footprint ensures compatibility with sealers, stackers and liquid handlers
- Chimney well design reduces risk of carryover when filling wells
- 1 and 2 mL blocks have rounded bottoms and are ideal for use as mother plates when creating compound libraries using automated liquid handling instruments
- Chimney well design of the 1 mL block reduces risk of cross-contamination between wells during pipetting, while the 2 mL block has cut-away sides for easy handling and placement on robotic systems
- Industry standard footprint offers compatibility with sealers, stackers and liquid handlers
- Participating in our customized barcode service enables truly accurate tracking of plates

Unmatched Packaging

- Matrix microplates are available in quick tear bags for easy access
- Inverted plate packaging prevents potential contamination associated with reaching across a plate and facilitates loading of stackers in the correct orientation right out of the bag

Optional Barcoding

- Program flexibility includes choice of code types, color-codes and barcode positions
- Stringent quality control provides barcodes that are accurate, readable and applied correctly
- Each code is verified against a database to guarantee no duplicates at your laboratory site



| Cat. No. | Description |
|----------|---|
| 4211MTX | 1 mL deep-well; 40 blocks/case |
| 4212MTX | 1 mL deep-well; Sterile; 40 blocks/case |
| 4221MTX | 2 mL deep-well; 40 blocks/case |
| 4222MTX | 2 mL deep-well; Sterile; 40 blocks/case |
| 4325MTX | 384-well; 80 blocks/case |
| 4326MTX | 384-well; Sterile; 80 blocks/case |



Microplate Instrumentation

The most complete solutions for all of your laboratory needs.

Thermo Scientific microplate instruments offer precise results and efficient performance, making them ideal for cell biology, molecular biology and immunology applications in the fields of cancer research, drug development, proteomics and genomics.

- Flexible, easy-to-use microplate washers support routine and varied research applications
- A choice of dedicated and multi-technology microplate readers offers a wide variety of applications
- Superior microplate incubator shakers deliver accurate temperature control and efficient orbital shaking





Microplate Instrumentation



| | |
|---|-----|
| Incubator/Shakers | |
| Thermo Scientific iEMS Incubator/Shaker | 114 |
| Thermo Scientific iEMS Incubator/Shaker HT | 115 |
| Multimode Readers | |
| Thermo Scientific Appliskan Multimode Reader | 116 |
| Thermo Scientific Varioskan Flash Multimode Reader | 118 |
| Microplate Fluorometer and Luminometers | |
| Thermo Scientific Fluoroskan Ascent FL Microplate Fluorometer and Luminometer | 120 |
| Thermo Scientific Fluoroskan Ascent Microplate Fluorometer | 122 |
| Thermo Scientific Luminoskan Ascent Microplate Luminometer | 124 |
| Microplate Photometers | |
| Thermo Scientific Multiskan EX Microplate Photometer | 126 |
| Thermo Scientific Multiskan FC Microplate Photometer | 127 |
| Thermo Scientific Multiskan GO Microplate Spectrophotometer | 128 |
| Thermo Scientific Multiskan Spectrum Microplate Spectrophotometer | 130 |
| Microplate Washers | |
| Thermo Scientific Wellwash 4 Mk2 Microplate Washer | 132 |
| Thermo Scientific Wellwash Microplate Washer | 133 |
| Thermo Scientific Wellwash Versa Microplate Washer | 134 |
| Software | |
| Thermo Scientific Ascent Software | 136 |
| Thermo Scientific Skanlt Software | 137 |



Thermo Scientific iEMS Incubator/Shaker



The iEMS* Incubator/Shaker is a high-performance microplate incubator and shaker designed for ELISA applications.

The iEMS Incubator/Shaker provides superior temperature control and efficient orbital shaking to dramatically increase sensitivity and specificity of ELISA assays. It reduces incubation times, providing high performance and productivity that meet even the highest assay demands. Up to nine 96-well plates can be processed simultaneously, delivering constant incubation temperatures up to 40°C.

Details

- Unique design of the iEMS thermal microplate holder eliminates the edge effect, ensuring even microplate heating from all sides
- Accurate temperature uniformity across the plate provides high reliability for assays requiring elevated temperatures
- The unit incorporates a powerful variable-speed orbital shaker, offering efficient mixing for even very viscous liquids
- Shaking motion enhances the reaction in wells and reduces incubation times, increasing throughput and productivity

Includes: Nine iEMS thermal holders (Cat. No. 5921200)

Warranty: One year

| Specifications | |
|--------------------------------|---|
| Programmable Temperature Range | 14°C to 40°C |
| Controlled Incubation Range | Ambient 3°C to 40°C |
| Resolution | 0.1°C |
| Programmable Incubation Time | Up to 48 h in step of 1 s |
| Inaccuracy | ±0.3°C |
| Uniformity | <0.3°C across whole plate |
| Evaporation | ≤ 300 mg/plate after 1 hour (at 37°C, without plate sealer) |
| Shaking Frequency | 400 to 1400 rpm in step of 250 rpm (5 speeds) |
| Diameter | 1 mm (radius 0.5 mm) |
| Programmable Shaking Time | Up to 48 h in step of 1 second |
| Programmable Interval Time | Up to 48 h in step of 1 second |

| Cat. No. | Description |
|----------|--|
| 5112200 | iEMS Incubator/Shaker 220-240V, 50/60 Hz |
| 5112207 | iEMS Incubator/Shaker 100-120V 50/60 Hz |
| 5921200 | iEMS Thermal Holder |

Thermo Scientific iEMS Incubator/Shaker HT

The iEMS Incubator/Shaker HT is a high-performance microplate incubator and orbital shaker for assays requiring temperatures up to 69°C.

The iEMS Incubator/Shaker HT offers accurate temperature control and orbital shaking for increased productivity. It accommodates three 96-well plates in a small footprint, making it perfect for medium-sized laboratories.

Designed for extended temperatures, the iEMS Incubator/Shaker HT is ideal for applications, such as DNA hybridization and primer extension assays.

Details

- Offers an extended temperature range up to 69°C for wider range of applications
- Unique design of the iEMS thermal holder HT eliminates the edge effect and ensures microplates are heated evenly from all sides
- Accurate temperature uniformity offers high reliability for assays requiring elevated temperatures
- Incorporates a powerful variable-speed orbital shaker, providing efficient mixing for even very viscous liquids
- Shaking motion enhances the reaction in wells and reduces incubation times, increasing throughput and productivity

Includes: Three iEMS thermal holders HT (Cat. No. 5921210)

Warranty: One year

| Specifications | |
|--------------------------------|--|
| Programmable Temperature Range | 14°C to 69°C |
| Controlled Incubation Range | Ambient 3°C to 69°C |
| Resolution | 0.1°C |
| Programmable Incubation Time | Up to 48 h in step of 1 s |
| Warming Speed | 35 min. from 24° to 65°C |
| Inaccuracy | ±0.5°C |
| Uniformity | <0.6°C across whole plate |
| Evaporation | No evaporation as film on the plate (All specs are with film on the plate) |
| Shaking Frequency | 400 to 1400 rpm in step of 250 rpm (5 speeds) |
| Diameter | 1 mm (radius 0.5 mm) |
| Programmable Shaking Time | Up to 48 h in step of 1 second |
| Programmable Interval Time | Up to 48 h in step of 1 second |

| Cat. No. | Description |
|----------|---|
| 5112250 | iEMS Incubator/Shaker HT 220-240V, 50/60 Hz |
| 5112257 | iEMS Incubator/Shaker HT 100-120V 50/60 Hz |
| 5921210 | iEMS thermal holder HT |



Thermo Scientific Appliskan Multimode Reader



The Appliskan* Multimode Reader is a filter-based multitechnology microplate reader for photometric, luminometric and fluorometric research applications.

The Appliskan reader offers all detection technologies in one compact and robust instrument for a wide range of application needs. It features an onboard shaker, an incubator and up to two onboard dispensers, and accommodates 6- to 384-well plates.

Details

- All technologies included – fluorometry, luminometry, absorbance, Time-Resolved Fluorescence (TRF) and FP
- Excellent sensitivity with luminometric and TRF applications
- Small footprint fits into even the most crowded laboratories
- Versatile Thermo Scientific SkanIt Software provides easy assay design, flexible data handling and convenient report formatting, even with challenging applications
- Numerous plate formats from 6- to 384-well
- Up to two dispensers for exact follow-up of kinetic reactions
- Fluorometric wavelength range up to 820 nm provides a high sensitivity and a wide linear dynamic range in all fluorometric applications, even with red-sensitive fluorochromes and demanding TR-FRET applications

Recommended for: Cellular assays, cytotoxicity and cell proliferation assays, DNA quantitation, ELISA/FIA/LIA assays, enzyme kinetic studies, europium assays, FRET and BRET assays, GPCR assays, Ion channel assays, kinase assays, multi-label assays, phagocytosis, protein assays, reporter gene assays and signal transduction

Includes: SkanIt Software, Excitation 485 nm and Emission 535 nm filters. Other filters available upon request.

Warranty: One year



| Specifications | |
|---|---|
| Excitation Wavelength Range (Fluorometry) | 200 nm to 1000 nm |
| Emission Wavelength Range (Fluorometry) | 360 nm to 820 nm |
| Sensitivity (Fluorometry) | Fluorescence intensity: <2 fmol fluorescein/well, 384-well plate; Time-resolved fluorescence: <20 amol Europium/well, 384-well plate |
| Precision | Fluorescence Polarization: <10mP 1nM Fluorescein, 96-well plate |
| Dynamic Range (Fluorometry) | >5 decades |
| Wavelength Range (Luminometry) | Standard mode: 360 nm to 820 nm; High-sensitive mode: 300 nm to 630 nm |
| Sensitivity (Luminometry) | High-sensitive mode: <10 amol ATP/well, 384-well plate; Standard mode: <200 amol ATP/well, 384-well plate |
| Dynamic Range (Luminometry) | >5 decades |
| Wavelength Range (Photometry) | 200 nm to 1000 nm |
| Measurement Range (Photometry) | 0 to 4.0 Abs |
| Linearity (Photometry) | 0 to 2.5 Abs (96-well plates) at 450 nm, $\pm 2\%$, 0 to 2 Abs (384-well plates) at 450 nm, $\pm 2\%$ |
| No. of Dispensers | Up to two |
| Dispensing Volume | 5 to 500 μL with 1 μL increments |
| Incubator Temperature Range | From ambient +4°C to 45°C, at ambient 25°C, at ambient 25°C |
| Shaking | Linear shaking |
| Plate Types | 6 to 384-well plates |
| Wavelength Selection | Filters (Ex/Abs 12.5 mm, Em 25.4 mm) |
| Measurement Types | Fluorescence intensity, time-resolved fluorescence, fluorescence polarization, absorbance and luminescence |
| Light Source | Xenon flash lamp |
| User Interface | Requires, but does not include, a personal computer |
| Computer Interface | Serial RS-232-C port |
| W \times D \times H | 14.8 \times 19.5 \times 13.5 in. (37.5 \times 49.5 \times 34 cm) |
| Weight | 27 kg (60 lbs.) |

| Cat. No. | Description |
|---------------------|-------------------------------|
| Instruments | |
| 5230000 | Appliskan |
| 5230010 | Appliskan with one dispenser |
| 5230020 | Appliskan with two dispensers |
| Upgrade Kits | |
| 460SP400 | First Dispenser Kit |
| 460SP420 | Second Dispenser Kit |

Thermo Scientific Varioskan Flash Multimode Reader



The Varioskan* Flash spectral scanning multimode reader offers optimal performance with unlimited wavelength selection for your most demanding research assays.

The Varioskan Flash reader provides flexibility for a variety of different applications, with unlimited wavelength selection, up to three onboard dispensers, unparalleled optical performance and the advanced Thermo Scientific SkanIt Software.

Varioskan Flash spectral scanning multimode reader includes fluorescence intensity, time-resolved fluorescence (TRF), photometric and optional luminometric detection technologies. It is designed for optimization and analysis of binding assays, ADMETox, molecular biology assays, enzyme kinetic studies, ion-channel and cell signaling assays, and other assays.

Details


- Unlimited wavelength selection for spectral analysis and measurement at any single wavelength: The optimal measurement wavelengths can be identified and easily selected for any assay at any time
- Extremely easy measurement setup with automated internal functions, offering maximum sensitivity and full dynamic range for every assay
- Onboard dispensers for exact follow-up of kinetic reactions: Essential for flash luminescence assays, Ca^{2+} flux studies and other rapid kinetic applications
- Increased assay throughput: Reads up to 1536-well plates and can easily be integrated with automated systems
- High performance incubator for temperature-critical assays
- Easy assay setup, flexible data handling and convenient report formatting with Thermo Scientific SkanIt Software



Increase throughput with a RapidStak Automated Microplate Stacker

Recommended for: Apoptosis assays; Ca^{2+} flux assays, cell proliferation, cellular assays, cytotoxicity and ADMETox, direct DNA, RNA and protein quantitation, ELISA/ FIA/TRF-ELISA assays, enzyme kinetic studies, europium assays, FRET assays, TR-FRET assays, BRET assays, GPCR assays, Ion channel assays, kinase assays, multilabel assays, reporter gene assays, signal transduction, tryptophan and tyrosine UV fluorescence

Warranty: One year



| Specifications | |
|--|---|
| Fluorescence Intensity/Time-Resolved Fluorescence | |
| Plate types | 6 - 1536-well plates |
| Wavelength selection | Double excitation and double emission monochromators |
| Excitation wavelength range | 200 - 1000 nm |
| Emission wavelength range | 270 - 840 nm |
| Excitation/emission bandwidth | 5 nm and 12 nm/12 nm |
| Sensitivity/dynamic range | Fluorescence intensity, top reading: < 0.4 fmol fluorescein/well, > 6 decades, 384-well plate Fluorescence intensity, bottom reading: < 4 fmol fluorescein/well, > 5.5 decades, 384-well plate Time-resolved fluorescence, top reading: < 120 amol europium/well, > 6 decades, 384-well plate |
| Luminometry | |
| Plate types | 6 - 1536-well plates, spectral scanning 6 - 384-well plates |
| Wavelength selection | All wavelengths, filters and double monochromators |
| Wavelength range | 360 - 670 nm, spectral scanning 270 - 840 nm |
| Sensitivity/dynamic range | < 7 amol ATP/well, > 7 decades, flash ATP reaction, 384-well plate |
| Photometry | |
| Plate types | 6 - 384-well plates |
| Wavelength selection | Double monochromators |
| Wavelength range | 200 - 1000 nm |
| Bandwidth | 5 nm |
| Linearity | 0 - 4 Abs (96-well plate) at 450 nm, $\pm 2\%$ 0 - 3 Abs (384-well plate) at 450 nm, $\pm 2\%$ |
| Accuracy | $\pm 2\%$ or 0.003 Abs, whichever is greater, at 200 - 399 nm (0 - 2 Abs) $\pm 1\%$ or 0.003 Abs, whichever is greater, at 400 - 1000 nm (0 - 3 Abs) |
| Precision | SD < 0.001 Abs or CV < 0.5%, whichever is greater, at 450 nm (0 - 3 Abs) |
| Dispensing | |
| No. of dispensers | Up to 3, automatic dispensing position control |
| Plate types | 6 - 384-well plates |
| Syringe size | 1 mL (standard), 5 mL (on request) |
| Dispensing volume | 1 - 10 000 μL , with 1 μL increments (1 mL syringe) Automatic safety control based on maximum well volume |
| Dispensing speed | 30 s, 96-well plate 80 s, 384-well plate (5 μL /well, 1 mL syringe, 0.40 mm tip) |
| Incubator | |
| | From ambient + 4°C to 45°C, at ambient 25°C |
| Shaker | |
| | Orbital with adjustable speed and diameter |
| General Features | |
| Measurement speed | 96-well plate in 15 s, 384-well plate in 45 s, and 1536-well plate in 135 s (minimum kinetic interval time from A1 back to A1) |
| Spectral scanning speed | < 2 s/well, 400 - 500 nm, 1 flash, 2 nm steps |
| Measurement types | Fluorescence intensity, time-resolved fluorescence, photometry and optional luminometry, all with spectral scanning |
| Light source | Xenon flash lamp |
| Dimensions (H x W x D) | 500 x 540 x 580 mm 19.7 x 21.3 x 22.8 in. |
| Weight | 55 - 67 kg (121 - 148 lbs.), depending on the configuration |

| Cat. No. | Model |
|----------|---|
| 5250030 | Varioskan Flash, top reading |
| 5250040 | Varioskan Flash, top and bottom reading |
| 5250500 | Varioskan LumiSens option, factory fitted (also enables luminometric spectral scanning) |
| 5250510 | Dispenser option, with 1 mL syringe, factory fitted |

Thermo Scientific Fluoroskan Ascent FL Microplate Fluorometer and Luminometer



The Fluoroskan Ascent* FL Microplate Fluorometer and Luminometer excels at rapid kinetic assays, cytotoxicity and related applications.

The Thermo Scientific Fluoroskan Ascent FL Microplate Reader is a compact and robust instrument with excellent optical performance for a variety of fluorometric and luminometric research applications. It offers versatile plate formats, fast reading speeds, up to three dispensers and top/bottom reading of plates.

Details

Fast reading speed

- Reads a 96-well plate in just 15 seconds, essential for kinetic applications
- Great for a wide variety of applications, including cytotoxicity, ion channel, FRET and BRET
- The combination of fluorometric GFP measurement and a luminometric luciferase measurement from the same well is fast due to the versatility of assay programming

Onboard dispensers for exact follow-up of kinetic reactions

- Equipped with up to three reagent dispensers for rapid kinetic assays, such as Ca^{2+} flux
- Supports simultaneous dispensing and reading, enabling monitoring of kinetic measurements from the very start of the reaction
- Very low dead volume and the backflush capability allow conservation of expensive reagents

High sensitivity for both top and bottom reading

- Fiberless direct illumination optics for both top and bottom reading
- High sensitivity, wide dynamic range and low cross talk ensure accurate and precise results for 1- to 384-well plates

Specially designed for automation

- Reads up to 384-well plates for increased assay throughput
- Easily integrated with automated systems

Recommended for: Ca^{2+} flux assays, cell proliferation, cytotoxicity, cell adhesion, DNA quantitation, reporter gene assays, hybridization assays, quantitation of PCR products, FRET assays, BRET assays, molecular beacon assays, immunoassays, enzyme activity, bacterial quantitation, phagocytosis, and oligonucleotide assays

Includes: PC software and filter pairs: Excitation: 355 nm/Emission: 460 nm, Excitation: 485 nm/ Emission: 538 nm

Warranty: One year



Increase throughput with a RapidStak Automated Microplate Stacker

| Specifications | |
|----------------------|---|
| Wavelength Range | Excitation: 320 nm - 700 nm; Emission: 360 nm - 670 nm |
| Excitation Filters | Up to eight filters in the excitation filter wheel. 355 nm and 485 nm filters included as standard. Other filters available upon request. |
| Emission Filters | Up to six filters in the emission filter wheel. 460 nm and 538 nm filters included as standard. Other filters available upon request. |
| Sensitivity | Fluorometry: 2 fmol fluorescein/well in a black 96-well plate Luminometry: 40 amol ATP/well using flash reaction, white 384-plate |
| Dynamic Range | Fluorometry: >6 decades Luminometry: >9 decades over whole gain setting area |
| No. of dispensers | Up to 3 |
| Dispensing volume | 1 to 1000 μ L in 1 μ L increments |
| Dispensing Speed | 25 seconds, 96-well plate, 5 μ L/well |
| Plate Types | 1- to 384-well plates |
| Measurement Speed | 15 seconds, 96-well plate |
| Wavelength selection | Filters |
| Light Source | Quartz-halogen Lamp |
| Detector | Photomultiplier Tube |
| Incubator | From ambient +3°C to 45°C, at ambient 25°C |
| Shaking | Orbital shaker |
| Computer Interface | Serial RS-232C port |
| User Interface | Requires a computer (not included) |
| Dimensions (HxWxD) | 13.4 \times 16.5 \times 16.5 in. (34 \times 42 \times 42 cm) |
| Weight | Basic unit: 21 kg (46 lbs.); three optional dispensers add 3.5 kg |

| Cat. No. | Description |
|---|--|
| 5210450 | Fluoroskan Ascent FL |
| 5210460 | Fluoroskan Ascent FL with one dispenser |
| 5210462 | Fluoroskan Ascent FL with two dispensers |
| 5210463 | Fluoroskan Ascent FL with three dispensers |
| <i>Includes PC software and filter pairs: Ex 355 nm/ Em 460 nm, Ex 485 nm/ Em 538 nm; Additional filters available upon request</i> | |

Thermo Scientific Fluoroskan Ascent Microplate Fluorometer



The Fluoroskan Ascent Microplate Fluorometer is a compact and robust instrument with excellent optical performance for a variety of research applications, including quantitation of double-stranded DNA in a solution.

The Fluoroskan Ascent Microplate Fluorometer is ideal for life science research applications, including fluorometric protein and enzyme studies, molecular interactions, nucleic acid quantification, reporter gene, fluorometric kinase, immuno and cell based assays.

Details

Fast reading speed

- Enables the most rapid throughput
- Ideal for a wide variety of applications, such as cytotoxicity, ion channel and FRET

Onboard dispensers for exact follow-up of kinetic reactions

- Equipped with up to three reagent dispensers for fast kinetic assays, such as Ca^{2+} flux
- Supports simultaneous dispensing and reading, enabling monitoring of fast kinetic measurements from the very start of the reaction
- Very low dead volume and the backflush capability allow conservation of expensive reagents

High sensitivity for both top and bottom reading

- Fiberless direct illumination optics for both top and bottom reading
- Ensures a high sensitivity, a wide dynamic range, and accurate and precise results for 1- to 384-well plates

Specially designed for automation

- The Fluoroskan Ascent fluorometer can read up to 384-well plates for increased assay throughput
- Easily integrated with automated systems

Recommended for: Ca^{2+} flux assays, cell proliferation, cytotoxicity, multi-drug resistance, cell adhesion, DNA quantitation, reporter gene assays, hybridization assays, quantitation of PCR products, FRET assays, molecular beacon assays, immunoassays, enzyme activity, neonatology, bacterial quantitation, phagocytosis, oligonucleotide assays and ADMETox

Includes: PC software and filter pairs: Excitation: 355 nm/Emission: 460 nm, Excitation: 485 nm/ Emission: 538 nm

Warranty: One year



Increase throughput with a RapidStak Automated Microplate Stacker

| Specifications | |
|---|---|
| Excitation Wavelength Range (Fluorometry) | 320 nm to 700 nm |
| Emission Wavelength Range (Fluorometry) | 360 nm to 800 nm |
| Excitation Filters (Fluorometry) | Up to eight filters in the excitation filter wheel. 355 nm and 485 nm filters included as standard. Other filters available upon request. |
| Emission Filters (Fluorometry) | Up to eight filters in the emission filter wheel. 460 nm and 538 nm filters included as standard. Other filters available upon request. |
| Sensitivity (Fluorometry) | 2 fmol fluorescein/well in a black 96-well plate |
| Dynamic Range (Fluorometry) | >6 decades |
| No. of dispensers | Up to 3 |
| Dispensing volume | 1 to 1000 μ L in 1 μ L increments |
| Dispensing Speed | 25 seconds, 96-well plate, 5 μ L/well |
| Plate Types | 1- to 384-well plates |
| Measurement Speed | 15 seconds, 96-well plate |
| Wavelength selection | Filters |
| Light Source | Quartz-halogen lamp |
| Detector | Photomultiplier Tube |
| Incubator | From ambient +3°C to 45°C, at ambient 25°C |
| Shaking | Orbital shaker |
| Computer Interface | Serial RS-232C port |
| User Interface | Requires a computer (not included) |
| Dimensions (H x W x D) | 13.4 x 16.5 x 16.5 in. (34 x 42 x 42 cm) |
| Weight | Basic unit: 21 kg (46 lbs.); three optional dispensers add 3.5 kg |

| Cat. No. | Description |
|----------|---|
| 5210470 | Fluoroskan Ascent |
| 5210480 | Fluoroskan Ascent with one dispenser |
| 5210482 | Fluoroskan Ascent with two dispensers |
| 5210483 | Fluoroskan Ascent with three dispensers |

Includes PC software and filter pairs: Ex 355 nm/ Em 460 nm, Ex 485 nm/ Em 538 nm; Additional filters available upon request



Thermo Scientific Luminoskan Ascent Microplate Luminometer



The Luminoskan Ascent Microplate Luminometer excels at luminometric applications, including reporter gene, immuno and cell based assays, enzyme studies, molecular interactions and more.

The Luminoskan Ascent Microplate Luminometer is a compact and robust instrument with excellent optical performance for a variety of luminometric research applications, offering versatile plate formats, fast reading speeds and up to three dispensers.

Details

Fast reading speed

- Reads a 96-well plate in just 15 seconds
- Essential for kinetic applications, such as enzyme kinetics and phagocytosis assays

Onboard dispensers for exact follow-up of kinetic reactions

- Equipped with up to three reagent dispensers for fast kinetic assays, such as ATP assays
- Supports simultaneous dispensing and reading, enabling monitoring of fast kinetic measurements from the very start of the reaction
- Very low dead volume and the backflush capability enable conservation of expensive reagents

High sensitivity

- DLReady* certified for excellent sensitivity
- High sensitivity, a wide dynamic range and low crosstalk ensure accurate and precise results in all luminometric assays


Specially designed for automation

- Reads up to 384-well plates for increased assay throughput
- Easily integrates with automated systems

Recommended for: Reporter gene assays, immunoassays with luminescent substrates, cell proliferation and cytotoxicity assays, intracellular Ca²⁺ assays, ATP assays, phagocytosis assays, reactive oxygen assays, microbiological assays, enzyme assays, BRET assays and ADMETox

Includes: PC software

Warranty: One year



Increase throughput with a RapidStak Automated Microplate Stacker



| Specifications | |
|--------------------------------|---|
| Wavelength Range (Luminometry) | 270 nm to 670 nm |
| Sensitivity (Luminometry) | 10 amol ATP/well using flash reaction, white 384-well plate |
| Dynamic Range (Luminometry) | >9 decades over whole gain setting area |
| No. of dispensers | Up to 3 |
| Dispensing volume | 1 to 1000 μ L in 1 μ L increments |
| Dispensing Speed | 25 seconds, 96-well plate, 5 μ L/well |
| Plate Types | 1 to 384-well plates |
| Measurement Speed | 15 seconds, 96-well plate |
| Wavelength selection | Filters |
| Detector | Photomultiplier Tube |
| Incubator | From ambient +3°C to 45°C, at ambient 25°C |
| Shaking | Orbital shaker |
| User Interface | Requires a computer (not included) |
| Computer Interface | Serial RS-232C port |
| Dimensions (H x W x D) | 13.4 x 16.5 x 16.5 in. (34 x 42 x 42 cm) |
| Weight | 46 lb. (21 kg) |

| Cat. No. | Description |
|----------|---|
| 5300160 | Luminoskan Ascent |
| 5300170 | Luminoskan Ascent with one dispenser |
| 5300172 | Luminoskan Ascent with two dispensers |
| 5300173 | Luminoskan Ascent with three dispensers |

*Includes PC software and filter pairs: Ex 355 nm/ Em 460 nm, Ex 485 nm/ Em 538 nm;
Additional filters available upon request*

Thermo Scientific Multiskan EX Microplate Photometer



The Multiskan* EX Microplate Photometer combines ease of use and proven reliability – an ideal benchtop microplate photometer for basic ELISA applications.

The Multiskan EX photometer is a basic reader for endpoint and kinetic assays, with PC software offering a wide range of data reductions and report generation.

Details

- Easy to use on-board software to run ELISA based applications
- Versatile Ascent Software provides flexible assay programming and data handling for research applications
- Wavelength range from 400 nm up to 750 nm
- Fast reading of 96 well plate
- Exceptional reliability and extended three year warranty
- Multiskan EX photometer with onboard software complies with the IVD (in vitro diagnostics) Directive 98/79/EC (Note: the EU Directive is not valid in North America)

Recommended for: Immunoassays, protein assays, growth curve and hormone assays, endotoxins, food diagnostics, HIV assays, hybridization assays, mini-sequencing assays, cytotoxicity, cell adhesion, signal transduction, endotoxins, antioxidants and food diagnostics

Warranty: Three years

| Specifications | |
|------------------------|---|
| Wavelength Selection | Filters |
| Wavelength Range | 400 nm to 750 nm |
| Filters | 8-position filter wheel; Standard filters: 405 nm, 450 nm, 620 nm; Other filters available upon request |
| Read-out Range | 0 to 3.5 Abs |
| Linearity (Photometry) | 0 to 2 Abs, $\pm 2\%$ at 405 nm |
| Accuracy (Photometry) | $\pm 2\%$ or 0.007 Abs, whichever is greater, typical value $\pm 1\%$ (0 to 2 Abs) at 405 nm |
| Precision | CV < 0.5% (0.3 to 1.5 Abs) at 405 nm CV < 1% (1.5 to 2 Abs) at 405 nm |
| Resolution | 0.001 Abs |
| Plate Types | 96-well plate |
| Measurement Speed | 5 seconds, 96-well plate |
| Light Source | Quartz tungsten halogen lamp |
| Detector | Eight (8) silicon photodetectors |
| Shaking | Linear shaking, three speeds |
| User Interface | On-board or PC control |
| Computer Interface | RS-232C serial interface |
| W × D × H | 16.5 × 12.6 × 5.5 in. (42 × 32 × 14 cm) |
| Weight | 11 kg |

| Cat. No. | Description |
|----------|---------------------------------|
| 51118177 | Multiskan EX 100-120V 50/60Hz |
| 51118170 | Multiskan EX 200-240V, 50/60 Hz |



Thermo Scientific Multiskan FC Microplate Photometer

The Multiskan FC Microplate Photometer is a filter-based microplate photometer with a 340 - 850 nm wavelength range, ideal for a range of applications from enzyme kinetic studies to Lowry assays.

The Multiskan FC photometer incorporates our 30-plus years of experience and leadership in microplate photometry. It delivers proven performance and reliable results for a wide variety of applications, with features including built-in self diagnostics, IQ/OQ/PQ and verification tools, and an easy to use visual user interface with multiple language availability.

Details

- Wide wavelength range of 340 nm to 850 nm
- Fast reading of both 96- and 384-well plates
- Shaking and incubation up to 50°C for temperature critical assays
- Ease of use through the large color screen
- Variety of languages (English, Chinese, French, German, Japanese, Portuguese, Russian, Spanish)
- Superior usability and logical workflow with Thermo Scientific SkanIt Software
- USB port for easy data transfer



Recommended for: Immunoassays (ELISA), protein assays, endotoxins, cytotoxicity and proliferation assays, enzyme assays, and growth curves

Warranty: One year

Download a variety of Multiskan FC ready-made software protocols from thermoscientific.com

Increase throughput with a RapidStak Automated Microplate Stacker

| Specifications | Description |
|----------------------|---|
| W x D x H | 11.4 x 15.7 x 8.3 in. (29 x 40 x 21 cm) |
| Wavelength Range | 340 nm to 850 nm |
| Wavelength Selection | Filters |
| Linearity (405 nm) | 96-well plate: 0 to 4 Abs, ±2%, normal mode; 0 to 3 Abs, ±2%, fast mode 384-well plate: 0 to 3 Abs, ±2%, normal mode; 0 to 2.5 Abs, ±2%, fast mode |
| User Interface | On-board or PC control |
| Readout range | 0 to 6 Abs |
| Accuracy (405 nm) | ±1% (0.3 to 3 Abs), ±2% (3 to 4 Abs) |
| Precision (405 nm) | CV ≤0.2% (0.3 to 3 Abs), CV ≤1.0% (3 to 4 Abs) |
| Light Source | Quartz-halogen lamp |
| Detector | Photodetector |
| Display | High contrast color display with 480 x 272 dots |
| Plate Types | 96- and 384-well plates |
| Measurement Speed | 6 seconds, 96-well plate, fast mode; 12 seconds, 96-well plate, normal mode; 11 seconds, 384-well plate, fast mode; 33 seconds, 384-well plate, normal mode |
| Resolution | 0.001 Abs |
| Shaking | Linear shaking with three modes: Slow, medium, fast |
| Filters | 8-position filter wheel; filters installed: 405 nm, 450 nm, and 620 nm. Additional filters can be ordered separately. |
| Optional Incubator | Temperature range from ambient +4°C up to 50°C |

| Cat. No. | Description |
|----------|-----------------------------|
| 51119000 | Multiskan FC |
| 51119100 | Multiskan FC with incubator |

Thermo Scientific Multiskan GO Microplate Spectrophotometer



The Multiskan GO UV/Vis Microplate Spectrophotometer offers freely-selectable wavelengths for 96- and 384-well plates and various types of cuvettes.

The Multiskan GO spectrophotometer features a broad wavelength range, including UV area and path length correction, as well as fast reading speed. It is the ideal tool for virtually any photometric research application: DNA, RNA, protein analysis and more. It offers exceptional usability through its intuitive, user-friendly internal software and comprehensive Thermo Scientific SkanIt Software, both available with multiple language options. An automatic power save function reduces energy consumption. Additionally, the Multiskan GO unit meets the RoHS (Restriction of Hazardous Substances) directive.

The Multiskan GO spectrophotometer has been designed to deliver first-rate performance and high quality results with minimal user effort. At instrument start up, extensive automatic self diagnostics verifies all major instrument functionalities. Continuous runtime control of optical and mechanical functions also ensures stable day-to-day and year-after-year performance and reliability.

Details

- Freely-selectable wavelengths from 200 to 1000 nm meets the demands of various assays
- Both microplate and cuvette reading for any throughput requirements
- Very fast plate measurements and a full spectrum of a sample in less than 10 seconds
- High quality data guaranteed by extensive self diagnostics
- Energy consumption of the instrument is reduced more than 70% when the power save function is activated.
- Visual internal software on a large color screen for quick measurements
- Easy and logical assay setup for demanding assays with the powerful SkanIt Software
- A selection of eight languages (English, Chinese, French, German, Japanese, Portuguese, Russian, Spanish)

Recommended for: DNA and RNA quantitation and purity, protein assays, enzyme assays, kinetic assays, immunoassays, cell proliferation and cytotoxicity assays, apoptosis assays, reporter gene assays and GPCR assays

Warranty: One year

Download a variety of ready-made software protocols for the Multiskan GO spectrophotometer from thermoscientific.com



| Specifications | |
|--------------------------------------|---|
| Wavelength Selection | Monochromator |
| Light Source | Xenon flash lamp |
| Wavelength Range | 200-1000 nm with 1 nm steps |
| Read-out Range | Up to 4 Abs |
| Bandwidth | <2.5 nm |
| Linearity | 0-2.5 Abs, 2% at 450 nm (96-well plate, cuvette with beam window ≥ 2 mm) |
| Accuracy | 1.0% + 0.003 Abs (0-2.0 Abs) 2% (2.0-2.5 Abs) at 450 nm |
| Precision | SD<0.003 Abs or CV<0.5% at 450 nm (Precision mode) SD<0.003 Abs or CV<1% at 450 nm (Fast mode) |
| Plate Types | 96- and 384-well plates |
| Cuvettes | 12.5 (W) \times 12.5 (D) \times 40-58 (H) mm |
| Measurement Speed | 6 s with 96-well plate 10 s with 384-well plate (from A1 back to A1) |
| Plate Shaking | Linear |
| Spectral Scanning Speed | 10 s from 200 to 1000 nm with 1nm steps |
| Incubation | From ambient +2°C to 45°C |
| User Interfaces | Stand alone use: 4.5-inch color display and keypad PC-control: SkanIt Software |
| USB Connections | PC Memory device port for data export External printer (HP PCL5 compatible) |
| Mains Input | 100-240V, 50/60Hz |
| Max Power Consumption | <110 W |
| Power Save Consumption | <2.5 W |
| Dimensions (W \times D \times H) | 11.2 \times 16.9 \times 10.2 in. (28.5 \times 43 \times 26 cm) |
| Weight | 10.8 kg (23.8 lbs.) |

| Cat. No. | Description |
|----------|------------------------------|
| 51119200 | Multiskan GO without cuvette |
| 51119300 | Multiskan GO with cuvette |

Thermo Scientific Multiskan Spectrum Microplate Spectrophotometer



The Multiskan Spectrum Spectrophotometer is a high performance instrument for endpoint, kinetic and spectral scanning applications, reading 6-to 384-well microplates and cuvettes.

The Multiskan Spectrum Spectrophotometer combines cuvette reading capabilities and advanced microplate reading in one convenient, compact instrument. It features unlimited wavelength selection, measurement at low UV to visible wavelengths, and flexible sample formats ranging from cuvettes to screening assays in 384 well plate formats.

The Multiskan Spectrum unit excels at nucleic acid and protein analysis, enzyme assays, cytotoxicity and cell proliferation assays and apoptosis assays.

Details

- Unlimited wavelength selection for both spectral analysis and for easy setup when changing assays
- Supports assays requiring measurement from low UV to visible wavelengths
- Reads both cuvettes and 6- to 384-well microplates, offering flexibility from single sample measurements to screening assays
- The Multiskan Spectrum spectrophotometer can easily be integrated with automated systems for increased assay throughput
- Onboard incubator for temperature-sensitive assays and onboard shaker with selectable speed
- Easy assay setup, flexible data handling and convenient report formatting with Thermo Scientific SkanIt Software

Recommended for: Nucleic acid analysis, enzyme assays, protein analysis, kinetic assays, cytotoxicity and proliferation assays, apoptosis assays and ELISA assays

Warranty: One year



Increase throughput with a RapidStak Automated Microplate Stacker



| Specifications | |
|-------------------------|---|
| Wavelength selection | Monochromator |
| Wavelength range | 200 to 1000 nm |
| Read-out range | 0 to 4 Abs |
| Linearity | 0 to 3 Abs, $\pm 2\%$ at 450 nm |
| Accuracy | $\pm 1.0\%$ or ± 0.005 Abs (0 to 2 Abs) $\pm 2.0\%$ (2 to 3 Abs) |
| Precision | SD < 0.005 Abs or CV $< 1\%$ (0 to 2 Abs) CV $< 2\%$ (2 to 3 Abs) |
| Bandwidth | 2 nm |
| Wavelength accuracy | ± 1 nm |
| Stray light | $< 0.02\%$ at 230 nm |
| Plate types | 6- to 384-well plates |
| Cuvette types | Two cuvettes, for sample and reference From standard to ultramicro; glass, plastic or quartz |
| Measurement speed | 96-well plate: 20 s typical 384-well plate: 60 s typical |
| Light source | Xenon flash lamp |
| Detector | Photodiodes |
| Incubation | From ambient $+4^{\circ}\text{C}$ to 45°C at ambient 25°C |
| Shaking | Plate: linear Cuvette: magnetic stirrer |
| W \times D \times H | 15.6 \times 17.7 \times 8.5 in. (39.7 \times 45 \times 21.7 cm) |
| Weight | 12.5 kg (27.5 lbs.) |

| Cat. No. | Description |
|----------|--|
| 51118600 | Multiskan Spectrum with cuvette, SkanIt Software Drug Discovery Edition |
| 51118650 | Multiskan Spectrum with cuvette, SkanIt Software Research Edition |
| 51118700 | Multiskan Spectrum without cuvette, SkanIt Software Drug Discovery Edition |
| 51118750 | Multiskan Spectrum without cuvette, SkanIt Software Research Edition |

Thermo Scientific Wellwash 4 Mk2 Microplate Washer



The Wellwash* 4 Mk2 Microplate Washer combines easy operation and reliable washing performance for excellent finishing results.

The Wellwash 4 Mk2 Microplate Washer is the ideal strip washer for routine ELISA applications.

The uncompromised performance combined with flexibility and easy-to-use wash cards make the Wellwash 4 Mk2 ideal for routine ELISA applications in clinical, veterinary, food and agricultural laboratories.

Details

- Unique co-axial wash heads provide high washing efficiency and contamination-free results needed for routine ELISA applications
- Pre-programmed wash cards offer a quick and easy start for routine users
- Programmable program cards and interchangeable 8- and 12-way wash heads deliver flexibility for research work
- Supplied with convenient, quick release bottle caps for easy buffer changes
- Pump shuts down automatically after 30 seconds to reduce wear and noise levels
- The Wellwash 4 Mk2 washer complies with the IVD (in vitro diagnostics) Directive 98/79/EC (Note: the EU Directive is not valid in North America)

Includes: One 8-way and one 12-way wash head, 2 L wash and 2 L waste bottles

Warranty: One year

| Specifications | |
|-----------------------------------|-----------------------------------|
| Wash/Waste Bottle Capacity | 2 L wash bottle, 2 L waste bottle |
| Priming Consumption | 15 to 20 mL |
| Residual Volume | <5 µL per well |
| Precision | 5% |
| Operating Cycles | Continuous |
| Wash Heads | 8- and 12-way |
| Operating Pressure (self-limited) | 0.5 bar |
| Weight | 10 kg |

| Cat. No. | Description |
|----------|------------------------------|
| 5160772 | Wellwash 4 Mk2 110-120V 60Hz |
| 5160770 | Wellwash 4 Mk2 220-240V 50Hz |

Thermo Scientific Wellwash Microplate Washer

The Wellwash Microplate Washer is easy to use and convenient for routine ELISA applications.

The Wellwash washer offers ease of use and convenience through an easy-to-view graphical user interface, local language versions and a USB port. Equipped with one wash bottle and 1 x 8 way wash head, it is designed for washing 96-well plates in routine ELISA applications.

Details

Easy and convenient to use

- Large color display and keyboard for easy operation
- Display is clear, showing everything needed in one place
- User interface is intuitive for quick protocol development
- Easy-to-use software requires minimal training for fast startup
- Unpressurized bottles are safe and secure to use
- Liquid level sensors are present in both wash and waste bottles for security
- Automatic prime feature provides safe performance
- Aerosol cover prevents aerosolization of infectious diseases
- Plate sensor recognizes presence of plate
- The Wellwash unit is designed to efficiently remove excess liquid from wells
- Low residual volumes deliver optimal washing performance and reliable assay results



Recommended for: ELISA-based assays

Includes: 1 x 2 L wash bottle, 1 x 2 L waste bottle, 1 x 8 wash head, aerosol cover

Warranty: One year

| Specifications | |
|------------------------|---|
| Number of Wash Bottles | 1 |
| Wash Heads (Optional) | 1 x 8, 1 x 12 |
| Plate Types | 96-well plates |
| Data Connection | USB |
| Cell Washing | No |
| Bottles | Non-pressurized |
| Display | LCD Color Display, 4.3 in. |
| Wash Volume | 50 to 1000 µL |
| Prime Volume | 5 to 100 mL |
| Rinse Volume | 50 to 100 mL |
| Dispense Volume | 50 to 400 µL |
| Residual Volume | <1.5 µL |
| Dimensions | 34.5 x 38.5 x 24 cm (13.6 x 15.2 x 9.4 in.) |
| Weight | 9 kg (20 lbs.) |

| Cat. No. | Description |
|----------|-------------|
| 5165000 | Wellwash |





Thermo Scientific Wellwash Versa Microplate Washer



With its superior versatility, the Wellwash Versa Microplate Washer is easy to use for ELISA and cell washing applications in research and routine laboratories.

The Thermo Scientific Wellwash Versa washer offers ease of use and convenience through a graphical user interface, local language versions and a USB port. Designed for high performance and versatility, it provides reliable and secure washing of 96- or 384-well plates, as well as sensitive cell washing.

The Wellwash Versa washer is ideal for research work. It includes two wash bottles and one rinse bottle, and the ability to select single or double 8- or 12-way wash heads, a double 8-way cell wash head, or a single 16-way wash head.

Details

Easy and convenient to use

- Large color display and keyboard for easy operation
- Clear display shows everything on one screen
- Intuitive user interface supports fast protocol development
- Easy-to-use software requires minimal training for quick startup

Reliable and secure performance

- Unpressurized bottles are safe and secure to use
- Liquid level sensors are present in both wash and waste bottles
- Automatic rinse and prime features provide safe performance
- An aerosol cover prevents aerosolization of infectious diseases
- Plate sensor recognizes presence of plate

Optimal washing performance

- Sweep mode ensures an extremely low residual volume in the well, resulting in reliable assay results
- Washing performance is further optimized by adjustable parameters, such as dispense and aspiration height and aspiration speed

Added versatility

- USB flash memory stick transfers protocols between two or more washers, and facilitates easy internal software updates and downloading of log files and reports
- Specially-designed wash is available for cell washing
- Choose from double 2 x 8 and 2 x 12 wash heads, 1 x 16 wash head, and two wash bottles and one rinse bottle for application flexibility

Recommended for: ELISA-based assays, cell washing

Includes: 2 x 2 L wash bottles, 1 x 2 L rinse bottle, 1 x 4 L waste bottle, 2 x 8 wash head, aerosol cover

Warranty: One year



Find best practices in washing at
www.thermoscientific.com/wellwash


| Specifications | |
|------------------------|--|
| Number of Wash Bottles | Up to 3 |
| Wash Heads (Optional) | 1 × 8, 1 × 12, 2 × 8, 2 × 12, 1 × 16, 2 × 8 cell |
| Plate Types | 96-well, 384-well plates |
| Data Connection | USB, RS-232-C |
| Cell Washing | Yes |
| Bottles | Non-pressurized |
| Display | LCD Color Display, 4.3 in. |
| Wash Volume | 50 to 1000 μ L (96), 20 to 300 μ L (384) |
| Prime Volume | 5 to 100 mL |
| Rinse Volume | 5 to 100 mL |
| Dispense Volume | 50 to 400 μ L (96), 20 to 120 μ L (384) |
| Residual Volume | <1.5 μ L |
| Dimensions | 34.5 × 38.5 × 24 cm (13.6 × 15.2 × 9.4 in.) |
| Weight | 9 kg (20 lb.) |

| Cat. No. | Description |
|----------|----------------|
| 5165010 | Wellwash Versa |



Increase throughput with a RapidStak Automated Microplate Stacker





Ascent* Software offers a flexible, easy-to-use way to control all models of the Ascent microplate reader family.

Ascent Software can be used with the following instruments:

- Fluoroskan Ascent
- Fluoroskan Ascent FL
- Luminoskan Ascent
- Multiskan EX

Ascent Software has been designed to operate on all microplate instruments belonging to the Ascent family, including fluorometers, luminometers and photometers. It features a clear and easy-to-follow approach with highly visual assay protocol setup, user-friendly data handling, and effective integration to automation and LIMS systems.

Details

Ascent Software is divided into two desktops: The Procedure Desktop controls the instrument actions and allows visual, easy setup of even the most complicated assays. The Results Desktop enables data to be reduced, calculated and reported according to assay requirements.

- Almost any assay can easily be created by selecting and adding the desired assay steps to the steplist
- Each step, such as Measure, Incubate or Shake, corresponds to an instrument action
- There can be up to 99 steps in one single assay
- Comprehensive in-built calculations, including Blank Subtraction, Curve Fit and Kinetic Data Reduction, are available
- The function tool can be used to create any customized calculation
- The Remote Control Interface enables easy integration with robotics and HIS/LIMS systems

| Cat. No. | Description | Port |
|-----------|--|-------------|
| 5185410CD | Fluoroskan Ascent and Fluoroskan Ascent FL | Serial port |
| 5185430CD | Multiskan EX | Serial port |
| 5185450CD | Luminoskan Ascent | Serial port |





Thermo Scientific SkanIt Software

SkanIt* Software is the ultimate tool for both microplate reader control and data handling.

Two editions of SkanIt Software are available: A Research Edition for scientists working in life science research, and a Drug Discovery Edition, offering features needed for compliance with the FDA's 21 CFR Part 11, for the drug discovery industry.

SkanIt Software is compatible with the following Thermo Scientific Microplate Readers: Varioskan Flash spectral scanning multimode reader, Multiskan GO spectrophotometer, Multiskan FC filter-based photometer, Multiskan Spectrum spectrophotometer and Appliskan filter-based multimode reader. This powerful software supports optimal use of the instrument's features and provides visual workflow, plus effortless data reduction with built-in calculations.

Details

- Create new protocols by adding steps to the steplist; each step corresponds to an instrument action (e.g. dispense, measure)
- For multi-label assays, several wavelengths can be measured almost simultaneously
- Follow real-time spectral scanning and kinetic curves during the measurement
- Select from a variety of built-in calculations, or create custom calculations
- Export data easily either manually or automatically
- Create comprehensive reports
- Remote control interface provides easy integration with LIMS and automated systems



| Cat. No. | Description | Port |
|----------|--|--------------------|
| 5187100 | SkanIt Software for Multiskan FC, Research Edition | USB |
| 5187080 | SkanIt Software for Varioskan Flash, Research Edition | Serial port or USB |
| 5187030 | SkanIt Software for Multiskan Spectrum, Research Edition | Serial port |
| 5187060 | SkanIt Software for Appliskan, Research Edition | Serial port |
| 5187119 | SkanIt Software for Multiskan GO, Research Edition | USB |
| 5187090 | SkanIt Software for Varioskan Flash, Drug Discovery Edition | Serial port or USB |
| 5187040 | SkanIt Software for Multiskan Spectrum, Drug Discovery Edition | Serial port |





Nucleic Acid Purification and Electroporation

Flexible KingFisher systems deliver high purity yields with quality results

Thermo Scientific KingFisher systems offer a rapid, automated and reproducible purification workflow by using magnetic rods to move particles through the purification phases of binding, washing and elution. After sample lysis, nucleic acids bind efficiently with coated magnetic beads in a suitable buffer. Washing then eliminates contaminants and elutes high purity DNA or RNA into the adjustable volume of elution buffer.

By eliminating laborious manual steps and purification problems, such as filter clogging and risk of contamination, the KingFisher system produces consistent, high purity yields over multiple experiments.





Nucleic Acid Purification



Thermo Scientific KingFisher Systems and Accessories

- Thermo Scientific KingFisher Flex Magnetic Particle Processors 140
- Thermo Scientific Consumables for KingFisher Flex Systems 141
- Thermo Scientific KingFisher mL Magnetic Particle Processors 142
- Thermo Scientific Consumables for KingFisher mL Systems 143
- Thermo Scientific KingFisher Magnetic Particle Processors 144
- Consumables for KingFisher Systems 144
- Thermo Scientific BindIt Software for KingFisher Instruments 145
- Thermo Scientific KingFisher Kits 146

Electroporation Cuvettes

- Molecular BioProducts Electroporation Cuvettes 149



Thermo Scientific KingFisher Flex Magnetic Particle Processors

The KingFisher* Flex System offers highly versatile, automated magnetic particle processing for DNA/RNA, protein or cell purification from virtually any source.

The KingFisher Flex Magnetic Particle Processor is a truly flexible solution for various sample processing requirements. With the 24 configuration, customers can raise the processing volume up to 5 mL, producing an increased yield of the purified product. For higher throughput needs, 96 samples can be processed with the 96 magnet head.

Using our revolutionary magnetic particle separation technology, the KingFisher Flex processor provides the fastest and easiest method for sample preparation, with excellent reproducibility and quality from a variety of sample materials.

Thermo Scientific KingFisher Kits enhance our unique nucleic acid purification workflow, providing an optimized high-throughput method for outstanding flexibility.



Recommended for: Genomics and proteomics, target identification, veterinary assays, biomarker discovery and quality control

Warranty: One year

Details

- Fully automated, high-speed purification of nucleic acids, proteins and cells
- High throughput: Up to 96 samples can be processed in under 20 minutes
- Expanded sample volume with 24 configuration increases the yield by ten-fold
- Open and flexible system allows use of any magnetic particle based kit to meet application demands
- Easy-to-use BindIt software provides instrument control, protocol creation and modification
- Optimized purification protocols and instructions for KingFisher Kits are available at www.thermoscientific.com/kingfisher

| Specifications | |
|--------------------------------------|--|
| Processing Volumes | 24: 200 μ L to 5 mL; 96: 20 to 1000 μ L |
| Capacity | 96 or 24 samples/run |
| Collection Efficiency of Particles | \geq 95% |
| Particle Size | ca. $>$ 1 μ m |
| Magnetic Rods | 24 or 96 in one frame; Interchangeable magnet heads |
| Plate Types | 24 or 96-well plates: KingFisher 24 deepwell plate (200 to 5000 μ L) Microtiter 96 deepwell plate (50 to 1000 μ L) KingFisher 96 plate (20 to 200 μ L) PCR plate (20 to 100 μ L) |
| Heating Temperature | From 4°C above ambient temperature to 96°C |
| Keyboard Display | START / PAUSE / STOP / OK / TURNTABLE ROTATING CLOCKWISE / TURNTABLE ROTATING COUNTERCLOCKWISE / four cursor keys / LCD |
| Dimensions (H \times W \times D) | 26.8 \times 23.6 \times 15 in. (68 \times 60 \times 38 cm) |
| Weight | 28 kg (62 lb.) |

| Cat. No. | Description |
|----------|---------------------------------------|
| 5400610 | KingFisher Flex with 96 PCR head |
| 5400620 | KingFisher Flex with 96 KF head |
| 5400630 | KingFisher Flex with 96 deepwell head |
| 5400640 | KingFisher Flex with 24 deepwell head |

Thermo Scientific Consumables for KingFisher Flex Systems

The Consumables for KingFisher Flex Systems are designed specifically to maximize the value of these instruments.


The KingFisher Flex system is compatible with KingFisher 24 Deepwell plates, Microtiter 96 Deepwell plates, KingFisher 96 plates and fully skirted rigid PCR plates. Specially designed tip combs that protect the magnets during the process are available for various plate types and applications.

KingFisher Flex instruments utilize disposable plastics made of polypropylene. The plastics – tip combs and microplates – are ideal for magnetic particle processing due to their low binding affinity for biomolecules. The special design enables excellent recovery of magnetic beads.



See a complete line of KingFisher Kits in this catalog for:

- Blood DNA
- Total RNA
- Cell and Tissue DNA
- Viral NA
- Plant DNA



| Cat. No. | Description | For Use with |
|----------|--|---|
| 95040450 | Microtiter Deepwell 96 plate | KingFisher Flex with 96 deepwell head |
| 97002514 | 96-tip comb for PCR magnet | KingFisher Flex with 96 PCR head |
| 97002524 | 96-tip comb for KF magnets | KingFisher Flex with 96 KF head |
| 97002534 | 96-tip comb for DW magnets | KingFisher Flex with 96 deep well head |
| 97002610 | KingFisher Flex 24 deepwell tip comb and plate | KingFisher Flex with 24 deep well head |
| 97002540 | 96 KF microplate (200 µL) | KingFisher Flex with 96 KF head, KingFisher Flex with 96 deep well head |
| 95040470 | KingFisher Flex 24 deepwell plate | KingFisher Flex with 24 deep well head |

Thermo Scientific KingFisher mL Magnetic Particle Processors



The KingFisher mL System introduces automated, low-throughput sample preparation into your laboratory workflow.

KingFisher mL magnetic particle processor is the best choice for higher processing volumes, up to 1 mL. It allows all purification steps to be carried out in a single strip of five tubes and 15 samples processed per run. With the capability of releasing target molecules as low as 50 μL , samples like DNA or RNA from larger starting volumes can be isolated and concentrated simultaneously.

Thermo Scientific KingFisher Kits complete the unique nucleic acid purification workflow, providing an optimized high-throughput method for exceptional flexibility.

Details

Based on the innovative technology of transferring magnetic particles instead of liquids, the KingFisher mL processor offers rapid and reproducible purification of high-quality DNA, RNA, proteins and cells for various types of downstream applications.

- High-speed purification of nucleic acids, proteins and cells
- Open and flexible system allows the use of any magnetic particle based kit to meet their application demands
- Easy-to-use BindIt Software provides instrument control, protocol creation and modification
- Optimized purification protocols and instructions for KingFisher Kits are available at www.thermoscientific.com/kingfisher

Warranty and Service Offering: One year

Recommended for: Genomics and proteomics, drug discovery, biomarker discovery, quality control and veterinary assays

| Specifications | |
|--------------------------------------|--|
| Processing Volume | 50 to 1000 μL |
| Capacity | 15 samples/run |
| Collection Efficiency of Particles | $\geq 95\%$ |
| Magnetic Particle Size | ca. $>1\mu\text{m}$ |
| Magnetic Rods | 3 \times 5 format |
| Vessel Type | Special tube strip, 1 \times 5 tubes |
| Tip Combo | Special design, 1 \times 5 format |
| Keyboard Display | START/STOP/TWO CURSOR KEYS/LCD |
| Dimensions (W \times H \times D) | 11.4 \times 11.4 \times 12.2 in. (29 \times 29 \times 31 cm) |
| Weight | 10 kg (23 lbs.) |


| Cat. No. | Description |
|----------|----------------------------------|
| 5400050 | KingFisher mL, 100-240V, 50/60Hz |

Thermo Scientific Consumables for KingFisher mL Systems


The Consumables for KingFisher mL Systems are made specifically to maximize the value of your KingFisher mL instruments.

The KingFisher mL instrument accommodates a maximum of 15 tube strips, which are compatible with the specially-designed tip comb. Each sample processing operation uses one tube strip containing five tubes; one tip comb with five tips is used for processing five samples at a time.

The KingFisher mL processor utilizes disposable plastics made of high quality polypropylene. The plastics, tip combs and tube strips are ideal for magnetic particle processing due to their low binding affinity for biomolecules. The special design also permits an effective recovery of magnetic beads.



| Cat. No. | Description |
|----------|--|
| 97002111 | KingFisher mL Tip comb, 800 pieces |
| 97002121 | KingFisher mL Tube, 20 × 45 pieces |
| 97002141 | KingFisher mL Combi 240, tubes and tip combs for 240 samples |



Thermo Scientific KingFisher Magnetic Particle Processors



The KingFisher System excels at purifying small-scale samples using magnetic particle processing.

Based on the innovative technology of transferring magnetic particles instead of liquids, the KingFisher Magnetic Particle Processor offers rapid and reproducible purification of high-quality DNA, RNA, proteins and cells for various types of downstream applications.

It is designed to automate the time-consuming sample preparation process of nucleic acids, proteins and cells from virtually any source. All purification steps are carried out in microplates with simple push button operation. Up to 24 samples can be processed per run.

Thermo Scientific KingFisher Kits complete the unique nucleic acid purification workflow, providing an optimized high-throughput method for superior flexibility.

Details

- Ideal for high speed purification of nucleic acids, proteins and cells
- Open and flexible system allows customers to use any magnetic particle based kit to meet their application demands
- Easy-to-use BindIt Software provides instrument control, protocol creation and modification

Recommended for: Genomics and proteomics, drug discovery, biomarker discovery, quality control, and veterinary assays

Warranty: One year

| Cat. No. | Description |
|----------|-------------------------------|
| 5400000 | KingFisher, 100-240V, 50/60Hz |

| Specifications | |
|--|--|
| Processing Volume | 20 to 200 μ L |
| Capacity | 24 samples/run |
| Collection Efficiency of the Particles | \geq 99% |
| Magnetic Particle Size | ca. $>1\mu$ m |
| Magnet Rods | 2 \times 12 format |
| Plate Type | Special design, 8 \times 12 grid format |
| Tip Comb | Special design, 2 \times 12 format |
| Keyboard Display | START/STOP/two cursor keys/LCD |
| Dimensions (W \times D \times H) | 11.4 \times 11.4 \times 12.2 in. (29 \times 29 \times 31 cm) |
| Weight | 10 kg (23 lbs.) |



Consumables for KingFisher Systems

Designed specifically for use with KingFisher instruments.

The KingFisher system uses disposable plastics made of high quality polypropylene. The plastics, tip combs and tube strips are ideal for magnetic particle processing due to their low binding affinity for biomolecules. The special design also permits an effective recovery of magnetic beads.

| Cat. No. | Description |
|----------|------------------------------|
| 97002070 | KingFisher tip comb |
| 97002080 | KingFisher plate 100 μ L |
| 97002084 | KingFisher plate 200 μ L |


Thermo Scientific BindIt Software for KingFisher Instruments

The BindIt Software Version 3.1 is a versatile tool for protocol creation, modification and control of KingFisher instruments.

Thermo Scientific BindIt Software has been designed to enable the development of custom-made protocols for custom applications. BindIt Software provides an easy-to-use graphical user interface that accelerates user productivity.

Details

- Protocols are created and stored in a PC database using BindIt Software
- Once a protocol has been created, the protocol can either be transferred to the KingFisher instrument memory or executed directly from the software
- Protocols run directly from the software are not stored in instrument memory
- Based on a steplist, the parameters for the active step are shown on the screen
- All steps have default parameters that can be changed according to the demands of the application
- Plates and reagents used are defined in the plate layout
- BindIt Software generates a status report, including the run log, plate layout and step parameters for QC purposes
- It enables the KingFisher Flex particle processor to interface with liquid handling, robotics and plate stacking instruments, providing a fully automated solution and the highest possible throughput



| Specifications | |
|------------------|--|
| Operating System | Windows XP Professional with Service Pack 2, Windows Vista; 32-bit Edition; Business Edition |
| RAM | 1 GB RAM |
| Drive | CD-ROM drive |
| Monitor | XVGA monitor with 1024 by 768 resolution |
| Port | One serial or USB port available |
| Free Disk Space | 0.5 GB |
| Browser | Microsoft Internet Explorer 6.0 (or greater) installed |

| Cat. No. | Description | Includes |
|----------|--|----------|
| 5189009 | BindIt 3.1 Software for KingFisher instruments | CD only |

Thermo Scientific KingFisher Kits



Blood DNA



Total RNA



Cell and Tissue DNA



Viral NA



Plant DNA

The KingFisher Kits complete the unique nucleic acid purification workflow, providing an optimized high throughput method for superior flexibility.

The KingFisher Nucleic Acid Purification Kits, together with the KingFisher magnetic particle processor and consumables, enable the rapid and reproducible purification of high quality DNA or RNA from a wide variety of materials. The superior performance of KingFisher Kits produces DNA or RNA that is free of inhibitors and is ready for use in many downstream applications. Easy-to-use protocols facilitate safe sample handling as well as minimal user input, maximizing walk-away time.

Details:

- Excellent performance and consistent results: Optimized to produce high yields of quality DNA or high integrity RNA
- Suitable for wide range of sample materials
- Large volume sample input possible with KingFisher Blood DNA Kit: Up to 3 mL of blood per sample
- Flexible throughput from one sample to up to 96 samples per run
- Walk-away solution for all throughputs
- Minimized hands-on time increases overall efficiency
- Compatible with most common downstream analysis and applications

All kits contain:

- KingFisher magnetic beads
- Lysis buffer, binding buffer, 2-3 wash buffers and elution buffer
- Kit-specific reagents to enhance performance



| Cat. No. | Description | Sample type | Instrument model | Format | Sample input | Samples per run | Reagents sufficient for | Elution volume | Run time/minutes |
|----------|---|---|------------------|----------------|--|-----------------|-------------------------|----------------|------------------|
| 97010060 | KingFisher Blood DNA kit 60 preps | Fresh or frozen whole blood, treated with EDTA or citrate | KingFisher mL | Tube strip | 250 µL | 15 | 60 samples | 150 µL* | 45 |
| 97010196 | KingFisher Blood DNA kit 1 × 96 preps | Fresh or frozen whole blood, treated with EDTA or citrate | KingFisher Flex | 96 magnet head | 250 µL | 96 | 96 samples | 150 µL* | < 35 |
| | | | | 24 magnet head | 3 mL | 24 | 9 samples | 500 µL* | 60 |
| 97020060 | KingFisher Total RNA kit 60 preps | Cells and tissues and samples stored in RNA stabilization reagent. | KingFisher mL | Tube strip | Up to 20 mg of tissue, up to 2 × 10 ⁶ cells | 15 | 60 samples | 150 µL* | 60 |
| 97020196 | KingFisher Total RNA kit 1 × 96 preps | Cells and tissues and samples stored in RNA stabilization reagent. | KingFisher Flex | 96 magnet head | Up to 20 mg of tissue, up to 2 × 10 ⁶ cells | 96 | 96 samples | 150 µL* | < 60 |
| 97030060 | KingFisher Cell and Tissue DNA kit 60 preps | Tissue, cells, cultured bacteria | KingFisher mL | Tube strip | Up to 20 mg of tissue, up to 1 × 10 ⁷ cells, up to 1 mL overnight culture of bacteria | 15 | 60 samples | 150 µL* | 25 |
| 97030196 | KingFisher Cell and Tissue DNA kit 1 × 96 preps | Tissue, cells, cultured bacteria | KingFisher Flex | 96 magnet head | Up to 20 mg of tissue, up to 1 × 10 ⁷ cells, up to 1 mL overnight culture of bacteria | 96 | 96 samples | 150 µL* | 25 |
| 97040060 | KingFisher Viral NA kit 60 preps | Cell-free body fluids such as serum, plasma, urine. Buccal swabs | KingFisher mL | Tube strip | 200 µL | 15 | 60 samples | 100 µL* | 50 |
| 97040196 | KingFisher Viral NA kit 1 × 96 preps | Cell-free body fluids such as serum, plasma, urine. Buccal swabs | KingFisher Flex | 96 magnet head | 200 µL | 96 | 96 samples | 100 µL* | 40 |
| 97050196 | KingFisher Plant DNA kit 1 × 96 preps | Fresh, frozen or lyophilized plant tissue and cells, such as leaves, seeds, grains. | KingFisher Flex | 96 magnet head | Up to 50 mg of fresh tissue or 10 mg of dried tissue | 96 | 96 samples | 150 µL* | < 30 |

DNA/RNA purified with all KingFisher Kits is suitable for PCR, Real Time PCR and other enzymatic downstream analysis methods
 * Recommended elution volume, but user-adjustable

| Cat. No | Description |
|----------|---|
| 97010060 | KingFisher Blood DNA kit 60 preps for KingFisher mL |
| 97010196 | KingFisher Blood DNA kit 1 × 96 preps for KingFisher Flex |
| 97020060 | KingFisher Total RNA kit 60 preps for KingFisher mL |
| 97020196 | KingFisher Total RNA kit 1 × 96 preps for KingFisher Flex |
| 97030060 | KingFisher Cell and Tissue DNA kit 60 preps for KingFisher mL |
| 97030196 | KingFisher Cell and Tissue DNA kit 1 × 96 preps for KingFisher Flex |
| 97040060 | KingFisher Viral NA kit 60 preps for KingFisher mL |
| 97040196 | KingFisher Viral NA kit 1 × 96 preps for KingFisher Flex |
| 97050196 | KingFisher Plant DNA kit 1 × 96 preps for KingFisher Flex |



*Depend on
the full range
of superior
Thermo Scientific
solutions...*

▶ Cell Culture

Essential products for your cell culture laboratory

Thermo Scientific Nunc cell culture products have been used by researchers worldwide for more than 50 years. We take pride in manufacturing products with consistent high quality to ensure you get the most reproducible and reliable results in your research.

With surface areas ranging from 0.013 cm² to 25,280 cm², Nunc cell culture products allow for easy scale up when expanding cultures, and our surfaces address most applications.



Molecular BioProducts Electroporation Cuvettes

The durable polycarbonate construction of Molecular BioProducts Electroporation Cuvettes enables them to withstand pulses required by many commonly used electroporators.

Molecular BioProducts electroporation cuvettes are thoroughly cleansed to eliminate trace inhibitor materials, and then sterilized with electronic beam irradiation. Gap tolerances are precisely set for reproducible field strength delivery and provide high transformation efficiency.

Details:

- Ideal for bacterial transformation and transfection
- Seamless plastic molding to prevent sample leaking
- Parallel aluminum plates for uniform sample treatment
- Frosted caps for easy labeling
- 1, 2 and 4 mm gap width available with color-coded caps
- Stringent manufacturing process to ensure reproducible results
- Individually wrapped, single-use packaging
- Pre-sterilized transfer pipetter included in the package to facilitate sample transfer

| Cat. No. | Description | Unit Packaging |
|----------|--------------------------------------|----------------|
| 5510 | 1 mm gap (Pre-sterilized, white cap) | 50 cuvettes |
| 5520 | 2 mm gap (Pre-sterilized, blue cap) | 50 cuvettes |
| 5540 | 4 mm gap (Pre-sterilized, green cap) | 50 cuvettes |





PCR and Real-Time PCR Instruments and Consumables

For genomics and related research, dependable thermal cyclers and high-quality consumables are the keys to productivity and success.

Thermo Scientific PCR products provide a range of PCR solutions for superior results.

Thermo Scientific Arktik and Piko thermal cyclers offer unparalleled performance and reliability for end-point PCR, whether working with individual reactions or high-throughput projects.

The PikoReal Real-Time PCR system is a highly compatible platform for personal bench top use, combining high performance and outstanding features.

Our ultra-thin wall tubes and plates represent the fast generation of PCR consumables, bringing significantly improved performance in fast PCR and realtime PCR assays. In combination with high quality UTW consumables, PCR and qPCR reactions are performed with high consistency, reduced reaction time, and less consumption of plastics and reagents.





PCR and Real-Time PCR Instruments and Consumables



Thermal Cyclers
 Thermo Scientific Arktik Thermal Cycler152
 Thermo Scientific Piko Thermal Cycler153

Real-Time PCR System
 Thermo Scientific Piko Real Real-Time PCR System154

PCR Tubes, Plates and Caps
 Thermo Scientific Piko PCR Plates156
 Thermo Scientific Piko Plate Illuminator157
 Sealing Films for Piko157
 Molecular BioProducts PCR Tubes159
 Molecular BioProducts PCR Plates and Caps160
 Molecular BioProducts PCR Strip Tubes and Caps160

Surface Decontamiant
 Molecular BioProducts DNA AWAY Surface Decontaminant161
 Molecular Bioproducts RNase AWAY Surface Decontaminant161
 Molecular BioProducts EasyStart PCR Mix-in-a-Tube162

Storage Reaction Tubes
 Molecular BioProducts HotStart Storage Reaction Tubes163



Thermo Scientific Arktik Thermal Cycler



The Arktik* Thermal Cycler is a reliable and flexible PCR instrument for the everyday requirements of any laboratory.


The Arktik Thermal Cycler is suitable for a wide range of applications, from individual reactions to high-throughput projects. It offers flexibility with three different interchangeable blocks: 96-well block, 384-well block and a dual block with two 48-well units. A broad temperature gradient simplifies temperature optimization. Ease-of-use is the key feature of the Arktik Thermal Cycler's user interface. The protocol is displayed graphically, and the programming is simple with intuitive menus.

Details

- **Interchangeable blocks** with excellent thermal precision
- **Dual 48-well block** for multi-user option
- Broad (up to 30°C) and accurate temperature gradient†
- **Reliable and easy to use**, with low noise level
- **Over-tightening protection** system in the heated lid
- Accepts virtually all standard PCR plastics

Ordering Alerts: Purchase of a heat block is required for operation of the Arktik Thermal Cycler.

Recommended for: PCR



| Specifications | |
|----------------------|--|
| Size (W × D × H) | 29 × 38 × 29 cm |
| Weight | 10.5 kg |
| Power Requirements | 100-240V, 50-60Hz |
| Block Configurations | 96 × 0.2 mL, 384 × 0.03 mL, 2 × 48 × 0.2 mL (dual block) Interchangeable |
| Max. Ramp Rate | Up to 3°C per second |
| Thermal Uniformity | ±0.4°C at 90°C |
| Thermal Accuracy | ±0.3°C |
| Thermal Range | 4 to 99.9°C |
| Gradient Range† | Max. 30°C |
| User Interface | Semi-graphical |
| Number of Programs | 4950 |
| Heated Lid | Manually adjustable, over-tightening protection system |
| Warranty | Two years |

| Cat. No. | Description |
|----------|--|
| TCA0001 | Arktik Thermal Cycler base unit with gradient† |
| TCA0002 | Arktik Thermal Cycler base unit without gradient |
| TCA0096 | 96-well block |
| TCA4848 | Dual block (2 × 48 wells) |
| TCA0384 | 384-well block |

† Gradient feature not available in the U.S. and Germany.

This product is licensed under U.S. Patent Nos. 5,552,580 and 5,496,517.

Notice: Purchase of this instrument conveys a limited, non-transferable immunity from suit for the purchaser's own internal research and development and applied fields other than human in vitro diagnostics under one or more of Canadian Patent 1,339,653, U.S. Patent 5,475,610 (claims 160-163 only) or corresponding claims in their unexpired non-U.S. counterparts owned by Applied Biosystems.

Thermo Scientific Piko Thermal Cycler

The Piko* Thermal Cycler delivers top performance with twice the speed and just half the size of regular thermal cyclers.

The Piko Thermal Cycler delivers high performance in a compact package. At just half the size of other PCR instruments, the Piko Thermal Cycler meets the highest criteria in thermal performance and can complete a PCR protocol in less than 15 minutes. This is achieved using unique technical advances that allow significant reduction in PCR run times and overall size of the instrument. The Piko Thermal Cycler is an ideal solution for both conventional and fast PCR applications. Piko Thermal Cyclers are available with two different block configurations: 24-well and 96-well. The 24-well cyclers accept all standard low profile single tubes, 8-tube strips and 24-well Piko PCR Plates. The 96-well instrument utilizes 96-well Piko PCR Plates.

Details

- **Superior thermal performance:** Consistent results from well-to-well due to extremely short settling times and high temperature uniformity across the block
- **Time saving:** PCR in as little as 15 minutes enabled by fast ramp rates and quick settling times
- **Reagent saving:** Repeatable yields from as little as 5 µL when using UTW vessels and heat sealers
- **Space saving:** One of the smallest footprints available, fits into a tiny bench space
- **Energy saving:** Short protocols and low wattage require only 25% of the power consumption of typical thermal cyclers
- **Repeatable sealing prevents sample evaporation:** Automatic heated lid provides consistent and tight sealing from run to run

| Specifications | |
|----------------------------|---|
| Size (W × D × H) | 16 × 17 × 23 cm |
| Weight | 4 kg (with external power supply) |
| Power Requirements | 100-240V, 50-60Hz |
| Block Configurations | 24-well (well volume 0.225 mL, sample volume max. 50 µL) 96-well (well volume 0.05 mL, sample volume max. 20 µL) |
| Max. Ramp Rate | >5°C heating and >4.5°C cooling |
| Typical Thermal Uniformity | ±0.3°C |
| Thermal Accuracy | ±0.2°C |
| Thermal Range | 4-99.9°C |
| User Interface | Semi-graphical and list mode programming |
| Number of Programs | Over 1000 |
| Heated Lid | Fully automated and motorized function |
| Warranty | Two years |

| Cat. No. | Description |
|----------|------------------------------|
| TCP0024 | Piko Thermal Cycler, 24-well |
| TCP0096 | Piko Thermal Cycler, 96-well |



Recommended for: PCR, High Performance PCR

Significant savings compared to conventional cyclers:

- 81% less energy consumed
- 71% less plastics used

This product is licensed under U.S. Patent Nos. 5,552,580 and 5,496,517.

Notice: Purchase of this instrument conveys a limited, non-transferable immunity from suit for the purchaser's own internal research and development and applied fields other than human in vitro diagnostics under one or more of Canadian Patent 1,339,653, U.S. Patent 5,475,610 (claims 160-163 only) or corresponding claims in their unexpired non-U.S. counterparts owned by Applied Biosystems.

Thermo Scientific PikoReal Real-Time PCR System



The PikoReal* Real-Time PCR System offers excellent performance in a small footprint.

The PikoReal Real-Time PCR System is a highly compatible gene quantification and genotyping platform in 24- and 96-well block formats.

The system is designed for a minimal footprint, making it ideal for personal benchtop use and field applications.

The PikoReal Real-Time PCR System incorporates innovative technologies to achieve fast performance while meeting the highest thermal requirements – all with reduced energy, plastics, and reagent consumption for real cost savings.

It also uses Ultra Thin Wall (UTW) microwell plates that are 25 percent the size of conventional microwell plates, yet compatible with multichannel pipetters, reagent dispensers and automated liquid handling systems.

Details

- **Five detection channels** enabling up to four-target multiplexing
- **High temperature uniformity** over the wells throughout the temperature range
- **Minimal footprint** for personal benchtop use
- **Piko format design** for significant savings in plate and reagent consumption
- **Automatic heated lid** for optimal temperature and sealing pressure
- **Remote control and monitoring** over Ethernet connection or stand-alone control using USB flash drives

Recommended for: Gene quantification and genotyping

| Specifications | |
|--|--|
| Thermal Block Formats | 24-well, 96-well (not interchangeable) |
| Sample Volume, Thermal Block | 10 to 50 μ L (PikoReal 24), 5 to 20 μ L (PikoReal 96) |
| Consumables, Thermal Block | 24-well or 96-well Piko PCR Plate; for 24-well block ; also 8-well strips and 0,2 mL single tubes |
| Max. Heating Rate, Thermal Block | >5°C/sec |
| Max. Cooling Rate | 4.5°C/sec |
| Temperature Range, Thermal Block | 4 to 99.9°C |
| Temperature Accuracy, Thermal Block | \pm 0.2°C |
| Temperature Uniformity, Thermal Block | \pm 0.3°C at 95°C |
| Temperature Range, Heated Lid | 30 to 110°C |
| Control, Heated Lid | Automatic temperature and pressure setting |
| Excitation | 5 LEDs |
| Excitation Range | 475 to 640 nm |
| Pre-Calibrated Dyes | FAM, SYBR Green, HEX, Yakima Yellow, ROX, Texas Red, Cy 5, SYBR Green for HRM (during 2011) |
| Multiplex | Up to 4 targets |
| Dynamic Range | 10 orders of magnitude |
| Sensitivity | 1 copy (theoretical) |
| Scan Time for Four Multiplexing Channels | <10 sec |
| Software Analysis Modes | Absolute quantification; relative quantification; melt curve analysis; allelic discrimination, high resolution melting (during 2011) |
| Operating Systems | Windows XP, Windows 7 |
| Communication | Ethernet (up to 10 instruments can be operated from a single PC) or stand-alone control with USB flash drive |
| Power Usage | 200W maximum |
| Dimensions (W x D x H) | 300 x 230 x 310 mm |
| Weight | 10 kg |

| Cat. No. | Description |
|----------|---------------------------------------|
| TCR0024 | PikoReal 24-well Real-Time PCR System |
| TCR0096 | PikoReal 96-well Real-Time PCR System |

Not available in U.S. before May 2011

Not available in Canada, Brazil, U.K., Germany, Austria, Switzerland, Italy, Spain, France, Belgium, The Netherlands, Luxemburg, Denmark and Sweden before May 2012

This product is licensed under U.S. Patent Nos. 5,552,580 and 5,496,517.

Notice: Purchase of this instrument conveys a limited, non-transferable immunity from suit for the purchaser's own internal research and development and applied fields other than human in vitro diagnostics under one or more of Canadian Patent 1,339,653, U.S. Patent 5,475,610 (claims 160-163 only) or corresponding claims in their unexpired non-U.S. counterparts owned by Applied Biosystems. The purchase of this product includes a limited, non-transferable license under specific claims of U.S. Patent Nos. 6,174,670, 6,569,627 and 5,871,908, owned by the University of Utah Research Foundation or Evotec Biosystems GmbH and licensed to Idaho Technology, Inc. and Roche Diagnostics GmbH.



The Ultra Thin Wall Technology in these tubes and plates represents the fast generation of PCR consumables, delivering significantly improved performance in fast PCR and real-time PCR assays. Each well wall is approximately 50 percent thinner than standard thin-walled tubes and plates. This further reduces the thermal barrier to heat flow into and out of the PCR sample, resulting in faster and more robust reactions.

Individual Tubes with Flat Caps

- Caps form a secure seal, yet are easy to open and close
- Ultra Thin Wall low profile available for fast PCR applications, such as with the Piko Thermal Cycler

Low Profile Strip Tubes with Attached Ultra-Clear Flat Caps†

- PCR strip tubes with separate caps for easy sample access
- Compatible with Piko thermal cyclers
- Ultra Thin Wall for fast PCR applications
- Low profile reduces dead space and increases PCR efficiency
- Ultra-clear flat caps ideal for use in qPCR assays

Piko PCR plates

- Modular Ultra Thin Wall PCR plates for Piko Thermal Cyclers
- Four plates, each the size of a microscope slide, can be reversibly snapped into a rigid frame, producing the equivalent of a standard 96-or 384-well microplate

Details

- Wall thickness is 50 percent of conventional thin-walled tubes
- Flexible format saves cost and reduces waste
- Compatible with standard multichannel pipettors and liquid handlers
- Well spacing and footprint correspond to industry (ANSI) dimensions
- Available in clear and white

Recommended for: PCR, Fast PCR, High Performance PCR

| Cat. No. | Color | Reaction Volume | Quantity |
|--|-------|-----------------|------------|
| Individual Tubes with Flat Caps | | | |
| TUC0010 | Clear | | 960 tubes |
| TUC0011 | White | | 960 tubes |
| Low Profile Strip Tubes with Attached Ultra-Clear Flat Caps | | | |
| TUC0080 | Clear | | 250 strips |
| TUC0081 | White | | 250 strips |
| 24-Well PCR Plate | | | |
| SPL0240 | Clear | 50 µL | 200 plates |
| SPL0241 | White | 50 µL | 200 plates |
| 96-Well PCR Plate | | | |
| SPL0960 | Clear | 20 µL | 200 plates |
| SPL0961 | White | 20 µL | 200 plates |
| Plate Frame | | | |
| SFR0241 | White | n/a | 50 plates |
| SFR0961 | White | n/a | 50 plates |

† Not available in U.S.


Thermo Scientific Piko Plate Illuminator

The Piko Plate Illuminator simplifies sample loading.

The Piko Plate Illuminator plate rack is state-of-the-art and provides a simple way to track the loading of reactants by illuminating the target well(s) from below with white LED lights.

- **Compatibility:** works with all standard single channel, 8-channel and 16-channel pipetters
- **Ease of use:** compact and easy-to-use instrument with pre-programmed loading patterns
- **Two-in-one design:** both 24-well and 96-well Piko PCR Plates fit into the rack

Recommended for: Sample handling, sample preparation



| Specifications | |
|------------------------|--|
| Dimensions (W x D x H) | 16 x 9 x 3 cm |
| Weight | 0.1 kg |
| Vessel Types | 24-well and 96-well Piko PCR Plates |
| Programs | 10 preset and user-selectable loading patterns |
| Power Requirements | 100-240V, 50-60Hz |

| Cat. No. | Description |
|----------|------------------------|
| PIP2496 | Piko Plate Illuminator |

Sealing Films for Piko

The Thermo Scientific optical adhesive sealing film is the perfect fit for Piko PCR Plates.

Details

- Optical clarity with minimal autofluorescence
- Tight sealing for sample volumes as low as 5 μ L
- Easy application without costly sealing devices
- Cryo-compatible adhesive is effective to -70°C

Recommended for: PCR, real-time PCR, qPCR, high performance PCR

| Cat. No. | Description | Quantity |
|----------|-------------------|--------------|
| ASF0020 | Adhesive film | 400 seals/Cs |
| TCS1080 | 8 Caps per strip | 120 strips |
| HSF0021† | Heat sealing film | 400 seals/Cs |

† For Piko only.

Biological Safety Cabinets

With unique technologies that create a very safe and comfortable working environment, our biological safety cabinets provide best-in-class safety, ergonomics and energy efficiency for today's most demanding laboratory applications.

Depend on the full range of superior Thermo Scientific solutions...



CO₂ Incubators

Unmatched choice. Advanced technology. Proven results. Our CO₂ incubators deliver innovative market-leading designs and application-based solutions – surrounding your samples with an environment you can trust to achieve your cell culturing goals.

Molecular BioProducts PCR Tubes

Molecular BioProducts PCR Tubes are engineered for consistent and reliable amplifications. Our PCR tubes are designed to provide uniform well contact for selected thermal cyclers. All products are free of RNase, DNase, DNA, ATP and endotoxins.

Details

- Secure seal lips prevent evaporation
- Thin-wall design for optimal heat transfer
- 100 percent virgin polypropylene for high clarity
- Packaged in puncture-proof bags for extra protection
- ZipLock* closure for easy resealing
- 10 individually sealed bags of 100 tubes to minimize cross contamination
- Available in both 0.2 and 0.5 mL sizes and assorted colors
- Special color packs available
- Specify desired color by adding the corresponding letter to the catalog number. Assorted colors are randomly selected. B-Blue, G-Green, O-Orange, P-Purple, R-Red, Y-Yellow, A-Assorted. Example: 3414B is 0.2 mL tubes in blue.



3412A

3414A

3430A

3432

| Cat. No. | Description | Volume | Unit Packaging |
|----------|--|--------|-----------------------|
| 3412 | Thin-wall tubes with flat caps (Natural) | 0.2 mL | 10 bags, 10 tubes/bag |
| 3412A | Thin-wall tubes with flat caps (Assorted colors) | 0.2 mL | 10 bags, 10 tubes/bag |
| 3414 | Thin-wall tubes with dome caps (Natural) | 0.2 mL | 10 bags, 10 tubes/bag |
| 3414A | Thin-wall tubes with dome caps (Assorted colors) | 0.2 mL | 10 bags, 10 tubes/bag |
| 3430 | Thin-wall tubes with flat caps (Natural) | 0.5 mL | 10 bags, 10 tubes/bag |
| 3430A | Thin-wall tubes with flat caps (Assorted colors) | 0.5 mL | 10 bags, 10 tubes/bag |
| 3432 | Thin-wall tubes with dome caps (Natural) | 0.5 mL | 10 bags, 10 tubes/bag |
| 3432A | Thin-wall tubes with dome caps (Assorted colors) | 0.5 mL | 10 bags, 10 tubes/bag |

Molecular BioProducts PCR Plates and Caps




Molecular BioProducts PCR Plates are compatible with all major thermal cyclers. They are molded under clean room conditions and individually wrapped to prevent contamination.

Details

- Fits tightly inside thermal cyclers
- Alphanumeric grid for sample tracking
- Elevated rims for sturdy non-flexing performance
- Plates are individually wrapped to prevent cross contamination

| Cat. No. | Description | Unit Packaging | Volume |
|----------|---------------------------------------|------------------------|--------|
| 3416 | 96-well PCR thin-wall plate (Natural) | 25 plates | 0.2 mL |
| 3418C | 8-PCR strip dome caps (Natural) | 10 bags, 12 strips/bag | n/a |
| 3417C | 12-PCR strip dome caps (Natural) | 10 bags, 8 strips/bag | n/a |
| 3481 | Adhesive plate sealing tape (Clear) | 100 sheets | n/a |



Molecular BioProducts PCR Strip Tubes and Caps

These PCR Strip Tubes and Caps are designed for easy reaction setup and reliable PCR amplifications. Under our stringent injection molding conditions, each cap perfectly aligns and fits to the corresponding tube on the strip.

Details

- Caps and tubes align and fit tightly in either direction
- Thin-wall allows effective heat transfer
- Available in natural and assorted colors
- Compatible with major thermal cyclers
- Specify desired color by adding the corresponding letter to the catalog number.
- Assorted colors are randomly chosen. B-Blue, G-Green, O-Orange, P-Purple, R-Red, Y-Yellow, A-Assorted. Example: 3418B is 8-PCR strip tubes in blue.

| Cat. No. | Description | Unit Packaging | Volume |
|----------|-------------------------------------|------------------------|--------|
| 3418 | 8-PCR strip tubes (Natural) | 10 bags, 12 strips/bag | 0.2 mL |
| 3418A | 8-PCR strip tubes (Assorted colors) | 10 bags, 12 strips/bag | 0.2 mL |
| 3418C | 8-PCR strip dome caps (Natural) | 10 bags, 12 strips/bag | n/a |
| 3417 | 12-PCR strip tubes (Natural) | 10 bags, 8 strips/bag | 0.2 mL |
| 3417C | 12-PCR strip dome caps (Natural) | 10 bags, 8 strips/bag | n/a |

Molecular BioProducts DNA AWAY Surface Decontaminant

Molecular BioProducts DNA AWAY* surface decontaminant eliminates unwanted DNA and DNase from glassware and plasticware without affecting subsequent DNA samples. It is ideal for removing DNA contamination from pipetters, PCR tube racks and bench surface prior to PCR set-up.

Details

- Degrades DNA more quickly and effectively than autoclaving
- Suitable for gel boxes, pipetters, benchtops, thermal cyclers and other equipment
- Use before performing PCR for a DNA-free work area
- Simply apply the ready-to-use formula to the surface to decontaminate, then wipe dry or rinse clean with water

| Cat. No. | Description |
|----------|--------------------------------------|
| 7010 | 250 mL bottle |
| 7008 | Canister of 25 wipes |
| 7009 | Box of 35 individually wrapped wipes |



Molecular Bioproducts RNase AWAY Surface Decontaminant

Molecular BioProducts RNase AWAY* formulation and wipes eliminate RNase from laboratory surfaces and prevent degradation of your precious RNA samples.

RNase AWAY decontaminant is ideal for decontaminating equipment and instruments, bench tops, glassware and plasticware. It reduces the dependency on carcinogenic DEPC treatments and saves time needed to bake glassware.

Details

- Use on pipetters, gel boxes or RNA prep areas
- Leaves no residue to interfere with gel polymerization or staining
- Chemically stable and nonabrasive
- Does not contain strong acids
- Apply to surface to decontaminate, then rinse with RNase-free water

| Cat. No. | Description |
|----------|--------------------------------------|
| 7002 | 475 mL spray bottle |
| 7003 | 1 L bottle |
| 7005 | 4 L bottles |
| 7000 | 250 mL bottle |
| 7006 | Canister of 25 wipes |
| 7007 | Box of 35 individually wrapped wipes |



Molecular BioProducts EasyStart PCR Mix-in-a-Tube



Molecular BioProducts EasyStart* mix reduces tedious PCR setup protocols to a single step.

EasyStart mix consists of a non-specific PCR reagent mixture hermetically sealed inside a thin-walled reaction tube, reducing the risk of contamination and providing optimum heat transfer. An aqueous solution of enzyme, template DNA and primers is added above the wax layer prior to thermal cycling; as the temperature rises, the wax melts, combining the reagents at the proper cycling temperatures.

Details

- Delivers better specificity and yield than conventionally prepared mixtures
- Compatible with a wide range of enzymes
- Reagent mix contains buffer, dNTPs, MgCl₂ and dH₂O in perfect proportion to ensure exceptional batch consistency and reproducible results
- Wax layer provides safe, easy storage of reagent mixture, prevents oxidation and evaporation and eliminates wasted reagents

Storage Conditions:

- Tubes can be stored at ambient temperatures for extended periods
- To remove PCR product, simply pierce the wax layer with a pipette tip

Includes: Choose 0.2 to 0.5 mL volumes; use 0.2 mL tubes for 0.2 mL heating blocks, and 0.5 mL tubes for 0.5 or 0.6 mL heating blocks

| Cat. No. | Description | Tube Volume | Reaction Volume |
|----------|------------------------------------|-------------|-----------------|
| 6028 | EasyStart Micro 20 (96 reactions) | 0.2 mL | 20 µL |
| 6228 | EasyStart Micro 20 (96 reactions) | 0.2 mL | 20 µL |
| 6020 | EasyStart Micro 50 (96 reactions) | 0.2 mL | 50 µL |
| 6024 | EasyStart Micro 100 (96 reactions) | 0.2 mL | 100 µL |
| 6022 | EasyStart 50 (96 reactions) | 0.5 mL | 50 µL |
| 6025 | EasyStart 100 (96 reactions) | 0.5 mL | 100 µL |




Molecular BioProducts HotStart Storage Reaction Tubes

Molecular BioProducts HotStart* Storage Reaction Tubes are thin-walled with pre-positioned wax bead to eliminate tedious oil or wax overlays.

HotStart Storage Reaction Tubes are ideal for storing pre-aliquoted master mixes in a freezer. To prepare the PCR cocktail, simply melt the wax bead and let it harden over reagents, forming a barrier to keep reagents separate until the annealing temperature is reached in the cyclor. The result: significant time savings while achieving dependable reproducibility.

Details

- Bead is stable at room temperature
- Provides synchronous reaction startup, thereby reducing mis-primers, primer dimers, and premature annealing
- Particularly recommended for complex DNA or cDNA, very low copy number targets, more than 30 thermal cycles, multiple PCR or RT-PCR
- After thermal cycling reactions, PCR product can remain in tube for later analysis; simply pierce wax layer with pipette tube to remove product
- Hinged caps provide tight seal
- Pre-sterilized; RNase- and DNase-free



| Cat. No. | Description | Tube Volume | Reaction Volume |
|----------------------|--|-------------|-------------------|
| Tubes | | | |
| 6002 | HotStart 50 (96 reactions) | 0.5 mL | 25 to 50 μ L |
| 6005 | HotStart 100 (96 reactions) | 0.5 mL | 60 to 100 μ L |
| 6008 | HotStart Micro 20 (96 reactions) | 0.2 mL | 15 to 25 μ L |
| 6010 | HotStart Micro 50 (96 reactions) | 0.2 mL | 25 to 50 μ L |
| 6014 | HotStart Micro 100 (96 reactions) | 0.2 mL | 60 to 100 μ L |
| 8-Tube Strips | | | |
| 6208 | HotStart Micro 20 Strips (96 reactions with dome caps) | 0.2 mL | 15 to 25 μ L |
| 6210 | HotStart Micro 50 Strips (96 reactions with dome caps) | 0.2 mL | 25 to 50 μ L |
| Bulk-Packed | | | |
| 6302 | HotStart 50 Bulk (480 reactions) | 0.5 mL | 25 to 50 μ L |
| 6305 | HotStart 100 Bulk (480 reactions) | 0.5 mL | 60 to 100 μ L |
| 6308 | HotStart Micro 20 Bulk (480 reactions) | 0.2 mL | 15 to 25 μ L |
| 6310 | HotStart Micro 50 Bulk (480 reactions) | 0.2 mL | 25 to 50 μ L |



Microcentrifuge Tubes and Racks

Ideal for all of your sample preparation applications.


Sample preparation requires proven microcentrifuge tubes and accessories that can be trusted in virtually any application.

We offer a broad range of tubes in a wide variety of sizes and shapes to meet your needs.

Rated for use up to 30,000 xg and built rugged with the highest quality materials, our products are designed for ease of use and are autoclavable, boil-proof and freezable to -80°C.



Microcentrifuge Tubes and Racks



Centrifuge and Microcentrifuge Tubes



| | |
|---|-----|
| Snap-Cap Centrifuge Tubes | 166 |
| Locking Lid Microcentrifuge Tubes | 166 |
| Capless Microcentrifuge Tubes | 167 |
| Screw Cap Microcentrifuge Tubes – Conical | 167 |
| Screw Cap Microcentrifuge Tubes – Free Standing | 168 |
| Screw Caps for Microcentrifuge Tubes | 168 |

Microcentrifuge and Microtiter Tubes


| | |
|---------------------------------|-----|
| Microtiter Tubes | 169 |
| Specialty Microcentrifuge Tubes | 169 |

Flipper Racks • Cryogenic and Microtube Racks

| | |
|---|-----|
| Molecular BioProducts 4-Way Flipper Racks | 170 |
| Molecular BioProducts 81-Well Cryogenic Rack with Lid | 170 |
| Molecular BioProducts FlipStrip Microtube Racks with Lids | 171 |
| Molecular BioProducts Reversible Microtube Racks with Lids | 172 |
| Molecular BioProducts 96-Well Flipper Microtube Racks with Lids | 172 |




Snap-Cap Centrifuge Tubes



Snap-cap tubes are easy to open and close, offering the most convenient choice for routine laboratory procedures including sample storage, centrifugation and preparation. Autoclavable, boil proof and freezable to -80°C. Made from virgin polypropylene.

| Cat. No. | Description | Color | Rated at |
|---|--------------------------------|---------|-----------|
| 3446 | 0.6 mL, Graduated, Non-Sterile | Natural | 30,000 xg |
| 3449 | 0.6 mL, Graduated, Sterile | Natural | 30,000 xg |
| 3448 | 1.5 mL, Graduated, Non-Sterile | Natural | 26,000 xg |
| 3451 | 1.5 mL, Graduated, Sterile | Natural | 26,000 xg |
| 3434 | 2.0 mL, Graduated, Non-Sterile | Natural | 25,000 xg |
| 3453 | 2.0 mL, Graduated, Sterile | Natural | 25,000 xg |
| Packaging: 0.6 mL available 1000/unit, 10 units/case; 1.5 mL/2.0 mL available 500/unit, 10 units/case | | | |

Locking Lid Microcentrifuge Tubes



Locking lid tubes offer the extra sample security found with screw cap tubes, yet have the same ease of use as snap cap tubes. Autoclavable, boil proof and freezable to -80°C. Made from virgin polypropylene.

| Cat. No. | Description | Color | Rated at |
|---|---|---------|-----------|
| 3454 | 0.6 mL, Graduated, Locking Lid, Non-Sterile | Natural | 30,000 xg |
| 3455 | 0.6 mL, Graduated, Locking Lid, Sterile | Natural | 30,000 xg |
| 3456 | 1.5 mL, Graduated, Locking Lid, Non-Sterile | Natural | 24,000 xg |
| 3457 | 1.5 mL, Graduated, Locking Lid, Sterile | Natural | 24,000 xg |
| 3458 | 2.0 mL, Graduated, Locking Lid, Non-Sterile | Natural | 24,000 xg |
| 3459 | 2.0 mL, Graduated, Locking Lid, Sterile | Natural | 24,000 xg |
| Packaging: 0.6 mL available 1000/unit, 10 units/case; 1.5 mL/2.0 mL available 500/unit, 10 units/case | | | |

Capless Microcentrifuge Tubes

These tubes are ideal for experiments where caps are not needed and storage is not required. Made from virgin polypropylene and certified RNase, DNase, DNA and pyrogen free.

Packaging: 0.6 mL available 1000/unit, 10 units/case; 1.5 mL/2.0 mL available 500/unit, 10 units/case

| Cat. No. | Description | Color | Rated at |
|----------|--------------------------------|---------|-----------|
| 3477 | 0.6 mL, Graduated, Non-Sterile | Natural | 30,000 xg |
| 3479 | 1.5 mL, Graduated, Non-Sterile | Natural | 26,000 xg |



Screw Cap Microcentrifuge Tubes – Conical

Conical screw cap tubes provide the same convenience as a snap cap with the security of a screw cap. Color caps are available for purchase separately. Autoclavable, boil proof and freezable to -80°C. Made from virgin polypropylene.

| Cat. No. | Description | Color | Rated at |
|----------|-----------------------------------|---------|-----------|
| 3462 | 0.5 mL, Tube Only, Non-Sterile | Natural | 25,000 xg |
| 3460 | 0.5 mL, Tube and Cap, Non-Sterile | Natural | 25,000 xg |
| 3431 | 0.5 mL, Tube and Cap, Sterile | Natural | 25,000 xg |
| 3460A | 0.5 mL, Tube and Cap, Non-Sterile | Amber | 25,000 xg |
| 3431A | 0.5 mL, Tube and Cap, Sterile | Amber | 25,000 xg |
| 3466 | 1.5 mL, Tube Only, Non-Sterile | Natural | 20,000 xg |
| 3464 | 1.5 mL, Tube and Cap, Non-Sterile | Natural | 20,000 xg |
| 3461 | 1.5 mL, Tube and Cap, Sterile | Natural | 20,000 xg |
| 3464A | 1.5 mL, Tube and Cap, Non-Sterile | Amber | 20,000 xg |
| 3461A | 1.5 mL, Tube and Cap, Sterile | Amber | 20,000 xg |
| 3470 | 2.0 mL, Tube Only, Non-Sterile | Natural | 24,000 xg |
| 3468 | 2.0 mL, Tube and Cap, Non-Sterile | Natural | 24,000 xg |
| 3463 | 2.0 mL, Tube and Cap, Sterile | Natural | 24,000 xg |
| 3468A | 2.0 mL, Tube and Cap, Non-Sterile | Amber | 24,000 xg |
| 3463A | 2.0 mL, Tube and Cap, Sterile | Amber | 24,000 xg |

Packaging: 500/unit, 10 units/case



Screw Cap Microcentrifuge Tubes – Free Standing



Free standing tubes are suitable for sample preparation, centrifugation, and storage. Color caps available for purchase separately. Autoclavable, boil proof and freezable to -80°C . Made from virgin polypropylene.

| Cat. No. | Description | Color | Rating |
|----------|-----------------------------------|---------|-----------|
| 3422 | 0.5 mL, Tube Only, Non-Sterile | Natural | 30,000 xg |
| 3472 | 0.5 mL, Tube and Cap, Non-Sterile | Natural | 30,000 xg |
| 3465 | 0.5 mL, Tube and Cap, Sterile | Natural | 30,000 xg |
| 3472A | 0.5 mL, Tube and Cap, Non-Sterile | Amber | 30,000 xg |
| 3465A | 0.5 mL, Tube and Cap, Sterile | Amber | 30,000 xg |
| 3478 | 1.5 mL, Tube Only, Non-Sterile | Natural | 18,000 xg |
| 3474 | 1.5 mL, Tube and Cap, Non-Sterile | Natural | 18,000 xg |
| 3467 | 1.5 mL, Tube and Cap, Sterile | Natural | 18,000 xg |
| 3474A | 1.5 mL, Tube and Cap, Non-Sterile | Amber | 18,000 xg |
| 3467A | 1.5 mL, Tube and Cap, Sterile | Amber | 18,000 xg |
| 3490 | 2.0 mL, Tube Only, Non-Sterile | Natural | 18,000 xg |
| 3488 | 2.0 mL, Tube and Cap, Non-Sterile | Natural | 18,000 xg |
| 3469 | 2.0 mL, Tube and Cap, Sterile | Natural | 18,000 xg |
| 3488A | 2.0 mL, Tube and Cap, Non-Sterile | Amber | 18,000 xg |
| 3469A | 2.0 mL, Tube and Cap, Sterile | Amber | 18,000 xg |

Packaging: 500/unit, 10 units/case

Screw Caps for Microcentrifuge Tubes



Includes a frosted writing area on cap. Caps with an O-ring form a positive seal against the rim of the tube to prevent leakage. Made from virgin polypropylene and certified RNase, DNase, DNA and pyrogen free.

| Cat. No. | Description | Color |
|----------|------------------------------------|---------|
| 3471 | Screw Cap with O-Ring, Non-Sterile | Natural |
| 3471R | Screw Cap with O-Ring, Non-Sterile | Red |
| 3471Y | Screw Cap with O-Ring, Non-Sterile | Yellow |
| 3471B | Screw Cap with O-Ring, Non-Sterile | Blue |
| 3471G | Screw Cap with O-Ring, Non-Sterile | Green |
| 3471A | Screw Cap with O-Ring, Non-Sterile | Amber |
| 3471O | Screw Cap with O-Ring, Non-Sterile | Orange |

Packaging: 500/unit, 10 units/case

Microtiter Tubes

Microtiter tubes meet the needs of many laboratory procedures and are available in a racked, 96-well footprint for robotic workstations. Autoclavable, boil proof and freezable to -80°C. Made from virgin polypropylene.

| Cat. No. | Description | Color |
|----------|---|---------|
| 3492 | 1.2 mL, Bulk, Non-Sterile | Natural |
| 3496 | 1.2 mL, 96 Tube Robotic Rack, Non-Sterile | Natural |
| 3487 | 1.2 mL, 96 Tube Robotic Rack, Sterile | Natural |
| 3426 | Strip Caps, Non-Sterile* | Natural |
| 3425 | Strip Caps, Sterile* | Natural |

**Made from non-autoclavable LDPE*



Specialty Microcentrifuge Tubes

Certified RNase, DNase, DNA, and pyrogen free. Specialty tubes are designed for small volumes and can be easily cut to remove pellets.

| Cat. No. | Description | Color |
|----------|--|---------|
| 3485 | 0.4 mL, Specialty Tube, Non-Sterile* | Natural |
| 3483 | 0.25 mL, Specialty Tube, Non-Sterile** | Natural |

**Made from non-autoclavable LDPE.
** Made from autoclavable Polypropylene.*



Sterility Certifications:

Non-sterile products are certified free of RNase, DNase, DNA and pyrogens. Sterile products are certified free of RNase, DNase, DNA, pyrogen, endotoxin, ATP, bioburden and PCR Inhibitors.

Molecular BioProducts 4-Way Flipper Racks




Molecular BioProducts 4-Way Flipper* Racks hold thirty-two 0.5 mL, thirty-two 1.5 mL, twelve 15 mL or four 50 mL tubes – each on a different side.

Details

- Holds 0.5 mL centrifuge/PCR tubes, 1.5 mL centrifuge tubes, 15 or 50 mL conical tubes on a different side
- Maximizes bench top space
- Can be used in freezers and waterbaths
- Autoclavable
- Available in fluorescent colors: Green, orange and pink

| Cat. No. | Description | Unit Packaging | Color |
|----------|--------------------|----------------|--------------------|
| 8850 | 4-Way Flipper Rack | 1 | Fluorescent Green |
| 8860 | 4-Way Flipper Rack | 1 | Fluorescent Orange |
| 8870 | 4-Way Flipper Rack | 1 | Fluorescent Pink |



Molecular BioProducts 81-Well Cryogenic Rack with Lid

Molecular BioProducts 81-Well Cryogenic Storage Racks secure and organize cryogenic storage vials. Easy-to-read alphanumeric labels simplify sample identification.

Details

- Holds 0.5, 1.5 and 2.0 mL cryovials
- Clear rack cover increases visibility
- Available in natural and fluorescent colors
- Alphanumeric grid for easy tube identification
- Stackable and more durable compared to cardboard freezer boxes

| Cat. No. | Color | Unit Packaging |
|----------|--------------------|----------------|
| 8800 | Natural | 1 |
| 8810 | Fluorescent Green | 1 |
| 8820 | Fluorescent Orange | 1 |
| 8830 | Fluorescent Pink | 1 |

Molecular BioProducts FlipStrip Microtube Racks with Lids

Molecular BioProducts FlipStrip* Microtube Racks are the perfect choice for PCR setup, holding 1.5 and 2.0 mL PCR reagent tubes and PCR tubes all in one rack.


Details

- Holds 0.2 mL PCR strips and 1.5 mL tubes on one side, 0.5 mL tubes on other
- Stable, low-profile design with an alphanumeric grid for easy tube identification
- Lid protects samples and allows racks to be stacked
- Several colors available
- Rack comes with clear lid
- Made from 100% recycled plastic

| Cat. No. | Description | Unit Packaging | Color |
|----------|-----------------------------------|----------------|---|
| 8608 | FlipStrip Microtube rack with lid | 12 | Assorted (3 each: red, blue, yellow and green) |
| 8618 | FlipStrip Microtube rack with lid | 12 | Assorted (3 each: fl. green, fl. pink, fl. orange and purple) |
| 8628 | FlipStrip Microtube rack with lid | 1 | Green |
| 8638 | FlipStrip Microtube rack with lid | 1 | Blue |
| 8648 | FlipStrip Microtube rack with lid | 1 | Purple |
| 8658 | FlipStrip Microtube rack with lid | 1 | Yellow |
| 8668 | FlipStrip Microtube rack with lid | 1 | Fluorescent Pink |
| 8678 | FlipStrip Microtube rack with lid | 1 | Fluorescent Orange |
| 8688 | FlipStrip Microtube rack with lid | 1 | Fluorescent Green |
| 8698 | FlipStrip Microtube rack with lid | 1 | Red |



Molecular BioProducts Reversible Microtube Racks with Lids



Molecular BioProducts Reversible Microtube Racks are a great choice for setting up PCR reactions in 0.2 or 0.5 mL PCR tubes. Use one side to hold 0.2 mL PCR tubes, or flip it over for 0.5 mL microcentrifuge or PCR tubes.

Details

- Holds 0.2 mL tubes on one side, 0.5 mL tubes on the other
- Alphanumeric grid for easy tube identification
- Lid protects samples and allows stacking
- Made from recycled plastic
- Several colors available

| Cat. No. | Description | Unit Packaging | Color |
|----------|-------------------------------------|----------------|---|
| 8600 | Reversible microtube rack with lids | 12 | Assorted (3 each: Red, Blue, Yellow and Green) |
| 8601 | Reversible microtube rack with lids | 12 | Assorted (3 each: Fl. Green, Fl. Pink, Fl. Orange and Purple) |
| 8620 | Reversible microtube rack with lid | 1 | Green |
| 8630 | Reversible microtube rack with lid | 1 | Blue |
| 8640 | Reversible microtube rack with lid | 1 | Purple |
| 8650 | Reversible microtube rack with lid | 1 | Yellow |
| 8660 | Reversible microtube rack with lid | 1 | Fluorescent Pink |
| 8670 | Reversible microtube rack with lid | 1 | Fluorescent Orange |
| 8680 | Reversible microtube rack with lid | 1 | Fluorescent Green |
| 8690 | Reversible microtube rack with lid | 1 | Red |

Molecular BioProducts 96-Well Flipper Microtube Racks with Lids

Molecular BioProducts 96-Well Flipper Microtube Racks are designed to hold either 96- 0.5 mL or 96- 1.5/2.0 mL style microcentrifuge tubes. Each rack comes with a clear lid that protects and secures samples and allows racks to be stacked.

Details

- Holds 0.5 mL tubes on one side, 1.5/2.0 mL tubes on the other
- Molded with an alphanumeric grid for easy tube identification
- Frosted writing surface for smudge-resistant labeling
- Several colors available

| Cat. No. | Description | Unit Packaging | Color |
|----------|----------------------------|----------------|--------------------|
| 8760 | 96-Well Flipper Rack w/lid | 1 | Fluorescent Pink |
| 8770 | 96-Well Flipper Rack w/lid | 1 | Fluorescent Orange |
| 8780 | 96-Well Flipper Rack w/lid | 1 | Fluorescent Green |
| 8790 | 96-Well Flipper Rack w/lid | 1 | Red |



Compliance and Calibration Services

Committed to maintaining the accuracy of your data.

The accuracy and reliability of your data mean everything. That's why our compliance and calibration services serve a critical role in maintaining the operational integrity of your instruments and equipment.

Our compliance services include an array of key offerings – from pipetter calibration and documentation to standards, maintenance and qualification/validation services – for all of your laboratory instruments, regardless of manufacturer.



Multi-vendor Compliance Services

Our comprehensive services cover GML/GLP facilities, as well as facilities regulated by other national (e.g. EPA) and international regulations, including EMEA and ICH.

These services include:

- **Calibration** — We will verify the accuracy of a measuring system on your instrument against a reference standard, which may involve the adjustment of this system to minimize deviations from the standard and thereby improve its measurement accuracy.
- **Qualification** — Whether Installation Qualification (IQ) for your new or relocated instruments or Operational Qualification (OQ) for current instruments, we will perform a quantitative performance assessment or qualitative performance assessment.
- **Relocation Services** — We will properly uninstall, relocate and reinstall your instruments or equipment, along with end-calibrations and any other requested Compliance Services, if requested.
- **Standard Operating Procedure Development** — We will work closely with your lab to provide a written method for controlling a practice under a routine set of conditions, in accordance with predetermined specifications, to obtain a desired outcome.
- **Validation Services** — We have the validation support products, services and technical expertise to help you meet your objectives, and can evaluate all of your equipment, software, information systems, and entire processes to ensure compliance.

Extended Warranty and Preventive Maintenance Contracts

Performing routine preventive maintenance on laboratory instruments is the best way to keep them operating at peak performance for accurate results. It's also the best way of extending the life of the instruments, maximizing your investment.

Most of the Thermo Scientific instruments showcased in this catalog are available with optional Extended Warranty and Preventive Maintenance contracts that eliminate the need for you to track routine maintenance for each instrument. Plus, you enjoy the peace of mind knowing you're protected beyond the warranty period.

IQ/OQ/PQ Services

As with Extended Warranty and Preventive Maintenance contracts, we offer optional Instrument Qualification (IQ)/Operational Qualification (OQ)/Performance Qualification (PQ) services for many of the instruments presented in this catalog.

Whether the need is for GMP/GLP labs seeking regulatory compliance, or for research labs simply looking to make certain that their new (or used) instruments are performing well and within published specifications, we offer both IQ/OQ/PQ testing protocols specifically developed to test the capabilities of each instrument, and IQ/OQ/PQ field packages, which include sending a field engineer to your lab to perform this testing.



Multi-vendor Pipetter Calibration, Maintenance and Repair Services

Pipettors are probably the most ubiquitous tools used in the life sciences industry, and liquid handling in general is central to virtually all life sciences research and testing. A liquid handling tool that is not performing correctly and is delivering inaccurate or inconsistent volumes can have drastic consequences. To ensure that pipettors are aspirating and dispensing accurate volumes of liquid, as well as performing within specifications, it is generally accepted that they should undergo preventive maintenance and recalibration procedures at least once each year.


As a leading manufacturer of pipettors for more than 25 years, we possess pipetting knowledge and expertise that is second to none. And knowing how critical regular maintenance and calibration is to pipetting performance, we offer expert level pipettor service for both Thermo Scientific and non-Thermo Scientific brand pipettors through an expanding network of PipetteLab* service centers worldwide. We also certify as *factory authorized* many independent distributors and pipettor service providers who receive formal pipettor service training through our training program, and whose facilities and processes meet our standards.

When your pipettors are serviced through a PipetteLab service center or Factory Authorized service provider, you are receiving the very best service available.

Coming Soon: Look for exciting developments in 2011, as we continue to greatly expand our PipetteLab service network and deliver valuable new tools to support your pipetting success.



To learn more about any of the services described above, please contact your local sales representative.



Leasing

At Thermo Fisher Scientific, we won't let equipment financing stand between you and your next great discovery.

Thermo Fisher Scientific is your partner for the best laboratory products and analytical technologies available, plus the unique financing options you need to accelerate success in science or industry.

Cost-effective financing designed for each individual customer is key to any successful equipment solution.

Regardless of your liquid handling equipment needs, our finance professionals have the expertise and industry insight that can help you meet your financing needs.

Benefits of Leasing

Flexible financing solutions allow us to create a program that best fits the needs of your facility. Leasing with Thermo Fisher Financial Services is easy. Here are some of the best reasons to finance your equipment:

Use of the equipment

Leasing allows you use of the equipment for an agreed-upon monthly payment — so you're able to pay as you use.

Tax and accounting benefits

Your lease payment may be tax-deductible depending on the type of lease you select. Your accountant or tax attorney can advise you.

Technology refresh

Make technology upgrades at any time during the term of the lease as new technology becomes available through Thermo Fisher Scientific.

Flexibility

Select a payment plan that best fits your budget.

Conservation of capital budget

Since your money is not tied up in equipment costs, you are free to spend it on other items, such as supplies, personnel and training.

Easier cash flow forecasting

Fixed monthly payments help you to better budget money for the term of the agreement.

Fixed payments

By locking in your payments now, you can avoid fluctuation in interest rates and inflation risk in the future.

Preserves credit


Leasing doesn't tie up your line of credit, providing more capital at your disposal when you need it.

Positive cash flow

Accelerate your ROI: Acquire income-producing equipment today with our deferred payment plans.

Easy to grow

With an optional Master Lease line of credit, you can easily add equipment throughout the year with a simple one-page lease agreement.



Standard Lease and Finance Options

Our goal is to provide flexible, competitive and comprehensive financial solutions for all of your instrument acquisitions, tailored to your specific requirements. This includes the financial flexibility you need to acquire the very best laboratory products and analytical technologies.

Thermo Fisher Financial Services offers Fair Market Value and One Dollar Buyout lease options with payment structures ranging from 12 to 60 months, plus payment options including step and skip payment structures, as well as quarterly or annual payment terms.

Fair Market Value Lease

The Fair Market Value Lease provides the lowest monthly payments and meets your concerns about equipment obsolescence and new technology by offering the following choices at the end of the lease term:

- You can choose to return the equipment to us and upgrade to the newest products available
- You can purchase the equipment for its Fair Market Value
- You can continue to lease the equipment

One Dollar Purchase Option

The One Dollar Purchase Option combines the benefits of leasing with those of owning your equipment. Your monthly payment will be slightly higher than that of a Fair Market Value lease, but you will have ownership of the equipment at the end of the term for a mere \$1.

Master Lease Line of Credit

A Master Lease Line of Credit allows you to acquire multiple instruments over an extended period of time without having to go through the credit approval process more than once. You only have to sign one lease agreement.

Reagent Agreement Program

This all-inclusive lease program allows you one convenient monthly payment that includes the equipment, service agreement and consumables or reagents.

Future Funds Program

The Future Funds program allows you to take your time in deciding if your payment plan is right for you. You will have a three-to-nine month window to choose whether to buy or lease your equipment. During this time, you pay only one percent of the equipment cost per month.

Emerging Credit

A structured lease program is often a good fit for development stage companies seeking to help conserve capital. This all-inclusive lease program allows you one convenient monthly payment that includes the equipment, service agreement, and consumables or reagents.




Index

- Appliskan Multimode Reader 116
- Arktik Thermal Cyclor 152
- Ascent Software 136
- Benefits of Leasing 177
- BindIt Software for KingFisher Instruments 145
- Capless Microcentrifuge Tubes 167
- Catalyst Express Microplate Handler 108
- Compliance and Calibration Services 174
- Consumables for KingFisher Flex Systems 141
- Consumables for KingFisher mL Systems 143
- Consumables for KingFisher Systems 144
- FILLit Software for Multidrop Reagent Dispensers 82
- Finnpipette Accessories 21
- Finnpipette Dispensers 20
- Finnpipette F1 GLP Kits 11
- Finnpipette F1 Multichannel Pipetters 10
- Finnpipette F1 Single Channel Pipetters 8
- Finnpipette F2 GLP Kits 15
- Finnpipette F2 Multichannel Pipetters 14
- Finnpipette F2 Single Channel Pipetters 12
- Finnpipette Multistepper Pipetter 17
- Finnpipette Novus Multichannel Electronic Pipetters 19
- Finnpipette Novus Single Channel Electronic Pipetters 18

- Finnpipette Stepper Pipetter 16
- Finntip BioCon Pipette Tips 29
- Finntip Compatibility Chart 24
- Finntip Extended Length Pipette Tips 27
- Finntip Filter Pipette Tips 30
- Finntip Flex Filter Pipette Tips 33
- Finntip Flex Pipette Tips 31
- Finntip Multistepper Pipette Tips 35
- Finntip Pipette Tips 25
- Finntip Stepper Tips 34
- Finntip Wide Orifice Pipette Tips 28
- Fluoroskan Ascent FL Microplate Fluorometer and Luminometer 120
- Fluoroskan Ascent Microplate Fluorometer 122
- iEMS Incubator/Shaker 114
- iEMS Incubator/Shaker HT 115
- KingFisher Flex Magnetic Particle Processors 140
- KingFisher Kits 146
- KingFisher Magnetic Particle Processors 144
- KingFisher mL Magnetic Particle Processors 142
- Locking Lid Microcentrifuge Tubes 166
- Luminoskan Ascent Microplate Luminometer 124
- Matrix D.A.R.T.s Tip Transfer Tool 96
- Matrix D.A.R.T.s Tips 93



page 21



page 57



page 61

Index

| | | | |
|---|-----|---|-----|
| Matrix Deepwell Storage Blocks | 111 | Microtiter Tubes | 169 |
| Matrix Disposable Automation Reservoirs | 97 | Molecular BioProducts 4-Way Flipper Racks | 170 |
| Matrix Equalizer Electronic Multichannel Pipettors | 40 | Molecular BioProducts 81-Well Cryogenic Rack with Lid | 170 |
| Matrix EXP Electronic Pipettors | 39 | Molecular BioProducts 96-Well Flipper Microtube Racks with Lids | 172 |
| Matrix Filter Pipette Tips | 51 | Molecular BioProducts and Pure Pipette Tips with MicroPoint Tip Design | 57 |
| Matrix Filtered D.A.R.T.s Tips | 94 | Molecular BioProducts and Pure Pipette Tips with SoftFit Tip Design | 58 |
| Matrix Hybrid ClipTip Pipette Tips (for Matrix Hybrid Pipettes) | 53 | Molecular BioProducts ART Barrier Ultra Micro Tips with MicroPoint Design | 61 |
| Matrix Hybrid Multichannel Pipettors | 43 | Molecular BioProducts ART* Barrier Tips with MicroPoint Design | 60 |
| Matrix Hybrid Single-Channel Pipettors | 42 | Molecular BioProducts BioRobotix ART Filter Tips for Automated Workstations | 98 |
| Matrix Hydra DT Pipetting Workstation | 88 | Molecular BioProducts BioRobotix Pipet Tips, Standard, Black | 100 |
| Matrix Hydra II Liquid Handling System | 87 | Molecular BioProducts BioRobotix Pipette Tips: Standard, Natural, Nonsterile | 101 |
| Matrix Memowell Pipetting Aid | 55 | Molecular BioProducts DNA AWAY Surface Decontaminant | 161 |
| Matrix Multichannel Electronic Pipettors | 37 | Molecular BioProducts EasyStart PCR Mix-in-a-Tube | 162 |
| Matrix Pipette Tip Compatibility Chart | 47 | Molecular BioProducts Electroporation Cuvettes | 149 |
| Matrix Pipette Tips | 48 | Molecular BioProducts FlipStrip Microtube Racks with Lids | 171 |
| Matrix Pipettor Accessories | 45 | Molecular BioProducts HotStart Storage Reaction Tubes | 163 |
| Matrix PlateMate Plus/WellMate Stacker Chimneys | 86 | Molecular BioProducts Low Retention Pipette Tips | 64 |
| Matrix Reagent Reservoirs | 54 | Molecular BioProducts Low Retention Pipette Tips with ART Self-Sealing Barrier | 62 |
| Matrix Serological Pipet Filler | 44 | Molecular BioProducts PCR Plates and Caps | 160 |
| Matrix Single-Channel Electronic Pipettors | 36 | Molecular BioProducts PCR Strip Tubes and Caps | 160 |
| Matrix TallTip Extended Length Pipette Tips | 50 | Molecular BioProducts PCR Tubes | 159 |
| Matrix TallTip Filter Pipette Tips | 52 | | |
| Matrix WellMate Disposable Tubing Cartridges | 84 | | |
| Matrix WellMate Microplate Dispenser | 83 | | |
| Matrix WellMate Stacker Base Unit | 85 | | |



page 90



page 114




page 133



| | | | |
|---|-----|--|-----|
| Molecular BioProducts Pipette Tips Compatibility Chart | 66 | RapidStak Automated Microplate Stacker with Polara RS Software . . . | 103 |
| Molecular BioProducts Pure Pipette Tips with SoftFit L Design | 59 | RapidStak DLL Programming Kits | 105 |
| Molecular BioProducts Reversible Microtube Racks with Lids | 172 | RapidStak Instrument Drivers | 104 |
| Molecular Bioproducts RNase AWAY Surface Decontaminant | 161 | Reagent Reservoirs | 22 |
| Multidrop 384 Reagent Dispenser | 78 | Screw Cap Microcentrifuge Tubes – Conical | 167 |
| Multidrop Combi nL FILLit Software | 82 | Screw Cap Microcentrifuge Tubes – Free Standing | 168 |
| Multidrop Combi nL Reagent Dispenser | 76 | Screw Caps for Microcentrifuge Tubes | 168 |
| Multidrop Combi Reagent Dispenser | 77 | Sealing Films for Piko | 157 |
| Multidrop Combi, 384 and DW Dispensing Cassettes | 80 | Skant Software | 137 |
| Multidrop DW Reagent Dispenser | 79 | Small Tube Dispensing Cassette | 80 |
| Multiskan EX Microplate Photometer | 126 | SMART dispensing cassettes | 80 |
| Multiskan FC Microplate Photometer | 127 | Snap-Cap Centrifuge Tubes | 166 |
| Multiskan GO Microplate Spectrophotometer | 128 | Specialty Microcentrifuge Tubes | 169 |
| Multiskan Spectrum Microplate Spectrophotometer | 130 | Standard Lease and Finance Options | 178 |
| Nalgene Disposable Robotic Reservoirs | 109 | Standard tube dispensing cassette | 80 |
| Nunc Disposable Plastic Reservoirs | 110 | Varioskan Flash Multimode Reader | 118 |
| Nunc Serological Pipets | 46 | Versette Accessories | 91 |
| Orbitor RS Microplate Mover | 107 | Versette ClipTip Automation Tips | 92 |
| Piko PCR Plates | 156 | Versette Pipetting Heads | 90 |
| Piko Plate Illuminator | 157 | Versette Pipetting Workstation | 89 |
| Piko Thermal Cycler | 153 | Wellwash 4 Mk2 Microplate Washer | 132 |
| PikoReal Real-Time PCR System | 154 | Wellwash Microplate Washer | 133 |
| PocketTip D.A.R.T.s Automation Tips | 95 | Wellwash Versa Microplate Washer | 134 |
| Polara RS | 105 | | |
| RapidStak Accessories: Microplate Stacks | 104 | | |



page 137



page 141



page 171

About Thermo Fisher Scientific

Thermo Fisher Scientific Inc. (NYSE: TMO) is the world leader in serving science. Our mission is to enable our customers to make the world healthier, cleaner and safer. With revenues of more than \$10 billion, we have approximately 35,000 employees and serve customers within pharmaceutical and biotech companies, hospitals and clinical diagnostic labs, universities, research institutions and government agencies, as well as in environmental and process control industries. We create value for our key stakeholders through two premier brands, Thermo Scientific and Fisher Scientific, which offer a unique combination of continuous technology development and the most convenient purchasing options. Our products and services help accelerate the pace of scientific discovery, and solve analytical challenges ranging from complex research to routine testing to field applications. Visit www.thermofisher.com.

Trademark Information

The Thermo Scientific brand name, the Thermo Scientific logo and the following trademarks are the property of Thermo Fisher Scientific, Inc. and/or its subsidiaries.

| | | | | |
|------------------|-------------|-------------|------------|-----------|
| ALPS | Dimension | HotStart | Multidrop | SmartMove |
| Appliskan | DNA AWAY | Hydra | Multiskan | SoftFit L |
| Arktik | Duraflex | iEMS | Nalgene | Spectrum |
| ART | EasyStart | Impact | Nunc | TallTip |
| Ascent Software | Equalizer | KingFisher | Orbitor | Variomag |
| Auto-Lock | Ergo | Luminoskan | Piko | Varioskan |
| BindIt | Fiberlite | Matrix | PikoReal | Versette |
| BioRobotix | FillIt | Memowell | Platemate | WellMate |
| Catalyst Express | Finnpipette | MicroPoint | Polaris | Wellwash |
| ClickSeal | Finntip | Microscan | RapidStak | |
| ClipTip | Flipper | Molecular | REACH | |
| Controlmate | FlipStrip | BioProducts | RNase AWAY | |
| D.A.R.T.s | Fluoroskan | Momentum | Skant | |

The following brands, trademarks or service marks are the property of the listed company and/or its subsidiaries. Every effort has been taken to ensure this list is accurate at the time of printing.

Agilent Technologies, Inc.

Velocity11

American National Standards Institute

ANSI

Beckman Coulter, Inc.

Biomek
FX
Span
Multimek

Biohit Oyj.

Proline

Biotek Instruments, Inc.

Precision

Brandtech Scientific, Inc.

Transferpette

Caliper Corporation

Zymark

CAPP aPs

Aero

Corning Incorporated

Eppendorf AG

Research
Response
Titerman

Gilson, Inc.

Pipetman

Hamilton Robotics, Inc.

Precision

Microsoft Corporation

Excel
Windows
Vista

Vista Lab Technologies

mLA

Molecular Devices, Inc.

FLIPR

Nichiryu Corp.

Benchmark
Oxford

PerkinElmer Inc.

CCS
Plate Trak

Promega Corp.

DLReady

Qiagen Group

Mettler-Toledo International Inc.

LTS
LiteTouch
Pipet Lite
Rainin

Simport Scientific

CryoVial
EDP


Socorex isba S.A.

Calibra

Tecan Trading AG

Te-MO
Tecan





© 2011 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

North America: USA/Canada 800 345 0206 or 800 522 7763

Europe: Austria +43 1 801 40 0, Belgium +32 53 73 42 41, France +33 2 2803 2180, Germany national toll free 08001-536 376, Germany international +49 6184 90 6940, Italy +39 02 02 95059 448, Netherlands +31 76 571 4440, Nordic/Baltic countries +358 9 329 100, Russia/CIS +7 (495) 739 76 41, Spain/Portugal +34 93 223 09 18, Switzerland +41 44 454 12 12, UK/Ireland +44 870 609 9203

Asia: China +86 21 6865 4588 or +86 10 8419 3588, India toll free 1800 22 8374, India +91 22 6716 2200, Japan +81 45 453 9220, Other Asian countries +852 2885 4613

Countries not listed: +49 6184 90 6940 or +33 2 2803 2180

www.thermoscientific.com